Dentsply Sirona

ProTaper Ultimate™

EN Instructions For Use 2
FR Instructions d'utilisation 8
IT Istruzioni per l'uso 15
DE Gebrauchsanweisung 21
ES Instrucciones de uso 28
PT Instruções de utilização 34
NL Gebruiksinstructies 40
SV Användarinstruktioner 46
DA Brugsanvisning 52
NO Brukerinstruksjoner 58
LT Naudojimo instrukcija 64
LV Lietošanas norādījumi 70
ET Kasutusjuhised 76
CS Návod k použití 82
SK Návod na použitie 88
HU Használati utasítás 94
RO Instrucţiuni de utilizare 100
SR Uputstva za upotrebu 106
BS Uputstva za upotrebu 112
CNR Uputstva za upotrebu 118
SQ Udhëzimet për përdorim 124
EL Οδηγίες χρήσης 130
TR Kullanım Kılavuzu 137
PT-BR Instruções de uso 143
KO 사용 지침 149

Dentsply Sirona

EN

ProTaper Ultimate™

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
STERILE – STERILIZED BY RADIATION
LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

REUSABLE STERILE ENDODONTIC INSTRUMENTS FOR HAND USE

1 PRODUCT DESCRIPTION

- Product description:
 ProTaper Ultimate™ are manual endodontic instruments, multiple use, made of nickel titanium alloy for the preparation of the root canal (glide path, shaping and finishing) during endodontic therapy. They are designed with an operative part of Nickel Titanium alloy heat treated, with cutting edges and a silicone handle. ProTaper Ultimate™ manual reusable files are available in sterile version only.
- Clinical benefits:
 The clinical benefits for the patient being treated with ProTaper Ultimate™ files is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Manual instruments intended for root canal preparation (glide path, shaping and debridement of the root canal).
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

ProTaper Ultimate™ instruments are delivered in sterile blisters.

The ProTaper Ultimate™ range comprises the following endodontic instruments:

Code **xx-length** **y-file size**	**Product**	**ISO XXX YYYv** **XXX-ISO size** **YYY-taper size** **v-variable taper**	**Available lengths**	**Packaging and content**
B ST PULH 619 0SX	ProTaper Ultimate™ Hand Use SX	ISO 020.003v	19 mm	Blister with 6 instruments (sterile)
B ST PULH 6 xx 0SL	ProTaper Ultimate™ Hand Use SLIDER	ISO 016.002v	21/25/31 mm	Blister with 6 instruments (sterile)
B ST PULH 6 xx 0SH	ProTaper Ultimate™ Hand Use SHAPER	ISO 020.004v	21/25/31 mm	Blister with 6 instruments (sterile)
B ST PULH 6 xx 0 yy	ProTaper Ultimate™ Hand Use FINISHER	F1: ISO 020.007v F2: ISO 025.008v F3: ISO 030.009v FX: ISO 035.012v	21/25/31 mm	Blister with 6 instruments (sterile)
B ST PULH 6 xx FXL	ProTaper Ultimate™ Hand Use FINISHER	FXL: ISO 050.010v	21/25/31 mm	Blister with 6 instruments (sterile)
B ST PULH 5 xx SQE	ProTaper Ultimate™ Hand Use Sequence	SLIDER: ISO 016.002v SHAPER: ISO 020.004v F1: ISO 020.007v F2: ISO 025.008v F3: ISO 030.009v	21/25/31 mm	Blister with 1 instrument of each instrument (sterile)
B DD PULH 525 SQE	ProTaper Ultimate™ Hand Use Demo Pack	SLIDER: ISO 016.002v SHAPER: ISO 020.004v F1: ISO 020.007v F2: ISO 025.008v F3: ISO 030.009v	25 mm	Blister with 1 instrument of each instrument (sterile)

5 COMPOSITION

- ProTaper Ultimate™ instruments are made of three main components: the working part made of a nickel-titanium alloy, a colored silicone stopper (except for instrument SX) and a colored silicone handle.
- The files are manufactured with a process that results in a file that has a colored appearance. Due to this proprietary processing, the files may appear slightly curved. This is not a manufacturing defect. While the file can be easily straightened using only your fingers, it is not necessary to straighten the file prior to use. Once inside the canal, the file will follow the root canal anatomy.

6 COMPATIBLE DEVICES

- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicone stopper to the correct working length is the appropriate method of working length determination.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use and the General Processing Instructions For Endodontic Products (see section [9]) to minimize the following risks to the device, the patient and/or the user:
 - Breakage of instrument.
 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Toxic or allergic reactions caused by processing residues.
- Sterility cannot be guaranteed if the packaging is open, damaged or wet.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).

- Use small size, flexible and/or NiTi files (this will help avoid canal transportation).
- Visually inspect the working part for all the defects listed in the former paragraph during use (i.e. after each wave).
- Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

All ProTaper Ultimate™ MANUAL instruments can be used with a watch-winding motion to progress into the canal. When light resistance is encountered, gently turn the file handle clockwise (CW) to engage and cut dentin. Continue this light inward engaging and cutting cycle until the flutes are loaded and the file simply spins freely. If the file is over-engaged, turn the handle counterclockwise to loosen, then use a CW rotation in order to load debris and haul it out of the canal.

1) Before use, check if the instrument is sterile or if processing is needed (refer to Processing Instructions).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root and its canal(s).
3) Prepare an access that enables the easy identification of each canal orifice.
4) The auxiliary shaping file, SX, may be used when there is restrictive space, or to pre-enlarge the body of a canal, or to relocate the coronal-most aspect of a canal away from an external root concavity.
5) In the presence of a lubricant, select the Slider and PASSIVELY follow the canal, in one or more passes, to its terminus. Determine WL using an electronic apex locator in combination with a radiographic image, then confirm patency.
6) If the Slider will not easily reach the canal terminus, select a small-sized manual Stainless Steel (SS) hand file. In the presence of a lubricant catheterize this canal, establish working length, confirm patency, and verify the Slidepath. Now, repeat step #4 above.
7) In the presence of NaOCl, select the Shaper and advance along the Slidepath, in one or more passes, until the WL is reached.
8) Upon removing the Shaper, irrigate, EndoActivate to break up debris and move it into solution, then re-irrigate to liberate this debris.
9) Reconfirm WL, especially in more curved canals.
10) Select the FINISHER F1 (020.07) and passively follow the canal to the WL, in one or more passes. Remove and inspect its apical flutes. When loaded with dentinal debris, the preparation is finished.

11) If the FINISHER F1 is loose at length and its apical flutes are not loaded with debris, select the FINISHER F2 (025.08) and use this file in the same manner described above for FINISHER F1.
12) If the FINISHER F2 is loose at length and its apical flutes are not loaded with debris, select the FINISHER F3 (030.09) and use this file in the same manner described above for FINISHER F1 and F2.
13) Upon removing any given file, clean and inspect its cutting flutes, irrigate, recapitulate with either a size 10 file or Endo Activator to break up debris, then re-irrigate.
14) Inspect the file's cutting flutes upon removal for the presence of unwinding, straightening or stretching. If deformation is noted, discard and use a new ProTaper Ultimate™ file.
15) The preparation is finished when the apical extent of any Finisher is loaded with debris, and the dedicated gutta-percha point or Size Verifier fits at the WL.

Use the auxiliary Finishers in larger and straighter canals only, such as maxillary central incisors, some palatal or distal canals of molars, or when there is a pathologic or iatrogenic defect:

16) Select manual auxiliary file, FINISHER FX (035.12), when working length is established and patency is confirmed. Use this stand-alone file directly and passively to follow the canal to the WL using the same technique as described above. Remove and inspect its apical flutes. When loaded with dentinal debris, the preparation is finished.
17) If the FINISHER FX is loose at length and its apical flutes are not loaded with debris, select manual auxiliary FINISHER FXL (050.10) and use this file in the same manner described above for FINISHER FX.
18) The preparation is finished when the apical extent of any auxiliary Finisher is loaded with debris, and the dedicated gutta-percha point or Size Verifier fits at the WL.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.
- Do not re-use silicone stops. Remove and discard silicone stops after each use.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE 2797

Authorised representatives

EU REP
DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP
DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH