Dentsply Sirona

ProTaper Next®

EN Instructions For Use 2
FR Instructions d'utilisation 8
IT Istruzioni per l'uso 14
DE Gebrauchsanweisung 20
ES Instrucciones de uso 26
PT Instruções de utilização 32
NL Gebruiksinstructies 38
SV Användarinstruktioner 44
DA Brugsanvisning 50
NO Brukerinstruksjoner 56
LT Naudojimo instrukcija 62
LV Lietošanas norādījumi 68
ET Kasutusjuhised 74
CS Návod k použití 80
SK Návod na použitie 86
HU Használati utasítás 92
RO Instrucţiuni de utilizare 98
SR Uputstva za upotrebu 104
BS Uputstva za upotrebu 110
CNR Uputstva za upotrebu 116
SQ Udhëzimet për përdorim 122
EL Οδηγίες χρήσης 128
TR Kullanım Kılavuzu 134
PT-BR Instruções de uso 140
KO 사용 지침 146

Dentsply Sirona

EN

ProTaper Next®

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
SINGLE USE – REPROCESSING PROHIBITED
STERILE – STERILIZED BY RADIATION

INSTRUCTIONS FOR USE

SINGLE USE STERILE ENDODONTIC INSTRUMENTS

1 PRODUCT DESCRIPTION

- Product description:
 ProTaper Next® are engine-driven tapered endodontic instruments made of nickel titanium alloy with cutting edges. These instruments are used to mechanically shape and prepare the root canals (shaping and debridement of the root canal) during endodontic therapy. ProTaper Next® files are single use instruments and are available in sterile version. Basic sequence consists of files (X1-X5) whereas an accessory file is available (XA), if necessary, for preflaring the orifice and improve radicular access when needed.
- Clinical benefits:
 The clinical benefit for the patient being treated with ProTaper Next® is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation (shaping and debridement of the root canal).
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

ProTaper Next® instruments are delivered in sterile blisters.

The ProTaper Next® range comprises the following endodontic instruments:

Code	Product	**ISO XXX YYYv** **XXX-ISO size** **YYY-taper size** **v-variable taper**	Available lengths	Packaging and content
A 0803	ProTaper Next®	X1: ISO 017 004v X2: ISO 025 006v X3: ISO 030 007v X4: ISO 040 006v X5: ISO 050 006v	21/25/31 mm	Blister with 3 instruments (sterile)
	ProTaper Next® Accessory Opener File	XA: ISO 019 035	19 mm	Blister with 3 instruments (sterile)
	ProTaper Next® Assortment	X1 X2 X3	21/25/31 mm	Blister with 1 instrument of each ISO (sterile)
A 0800	ProTaper Next®	XA: ISO 019 035 X1: ISO 017 004v X2: ISO 025 006v X3: ISO 030 007v X4: ISO 040 006v X5: ISO 050 006v	21/25/31 mm	Blister with 6 instruments (sterile)
	ProTaper Next® Accessory Opener File	XA: ISO 019 035	19 mm	Blister with 3 instruments (sterile)
	ProTaper Next® Assortment	X1 X2 X3	21/25/31 mm	Blister with 2 instruments of each ISO (sterile)
A 0802	ProTaper Next® DEMO	X1 X2 X3	25 mm	Blister with 1 instrument of each ISO (sterile)

5 COMPOSITION

ProTaper Next® instruments are made of three main components: the working part made of a nickel-titanium alloy (M-Wire®) (except for instrument XA which is made of standard NiTi), a silicone stopper and a gold-plated brass shank with colored rings (except for instrument XA).

6 COMPATIBLE DEVICES

- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).
- Torque control devices ensure optimal usage.
- Use with endodontic motors in a constant rotation at a speed of 300 rpm.
- Set torque at:

File size	Torque [Ncm]	Speed [rpm]
ProTaper Next® X1 ProTaper Next® X2 ProTaper Next® X3 ProTaper Next® X4 ProTaper Next® X5	2 - 3	300
ProTaper Next® XA	3 - 4	

- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicone stopper to the correct working length is the appropriate method of working length determination.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use to minimize the following risks to the device, the patient and/or the user:
 - Breakage of instrument.
 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Toxic or allergic reactions caused by processing residues.
- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Similar to all machine-driven root canal instruments, ProTaper Next® instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.

- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible and/or NiTi files (this will help avoid canal transportation).
 - Visually inspect the working part for all the defects listed in the former paragraph during use (i.e. after each wave).
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

1) Before use, make sure the instrument is sterile (sterility cannot be guaranteed if packaging is open, damaged or wet).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
3) Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.
4) Irrigate abundantly and frequently after each instrument used (according to good dental practices).
5) Scout the whole root canal with small-sized hand files, determine working length by using an apex locator in combination with radiographs, irrigate and confirm patency.

6) If preflaring of the orifice is necessary, use the ProTaper Next® XA (Accessory) file to improve radicular access by removing triangles of dentin, relocating the coronal aspect of a canal away from an external root concavity and creating more shape as desired.
7) Always irrigate and if necessary, expand the glide path using small-sized hand files or dedicated mechanical glide path files.
8) In the presence of NaOCl, brush and follow along the glide path with the ProTaper Next® X1 (ISO 017 004), in one or more passes, alternatively with small-sized hand files if necessary, until full working length is passively reached.
9) During protocol of use, irrigate, recapitulate with a small-sized hand file after each sequential ProTaper Next® instrument, then re-irrigate.
10) Use ProTaper Next® X2 (ISO 025 006), exactly as described for ProTaper Next® X1, until full working length is passively reached.
11) Inspect the apical flutes of the ProTaper Next® X2:
 ⇨ If they are loaded with dentin, then the shape is finished, the correspondingly sized gutta-percha master cone or size verifier may be fitted, and the canal is ready for disinfection.
12) Alternatively, gauge the foramen with an ISO 025 hand file:
 ⇨ If the instrument is snug at length, preparation is finished and the canal is ready to be obturated.
 ⇨ If the instrument is loose at length, continue shaping with ProTaper Next® X3 (ISO 030 007) and, when necessary, proceed with ProTaper Next® X4 (ISO 040 006) or ProTaper Next® X5 (ISO 050 006), gauging after each instrument respectively with a hand file of the same size. Use the same shaping technique up to working length.
13) Use dedicated paper points to dry and dedicated gutta-percha points to obturate the root canals.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- ProTaper Next® instruments are single use devices and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
- Disinfection, cleaning and sterilization are not applicable.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE 2797

Authorised representatives

EU REP
DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP
DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH