Dentsply Sirona

ProTaper Gold®

EN Instructions For Use 2
FR Instructions d'utilisation 8
IT Istruzioni per l'uso 15
DE Gebrauchsanweisung 22
ES Instrucciones de uso 29
PT Instruções de utilização 36
NL Gebruiksinstructies 43
SV Användarinstruktioner 50
DA Brugsanvisning 56
NO Brukerinstruksjoner 62
LT Naudojimo instrukcija 68
LV Lietošanas norādījumi 74
ET Kasutusjuhised 80
CS Návod k použití 86
SK Návod na použitie 92
HU Használati utasítás 98
RO Instrucţiuni de utilizare 105
SR Uputstva za upotrebu 112
BS Uputstva za upotrebu 118
CNR Uputstva za upotrebu 124
SQ Udhëzimet për përdorim 130
EL Οδηγίες χρήσης 136
TR Kullanım Kılavuzu 143
PT-BR Instruções de uso 149
KO 사용 지침 156

Dentsply Sirona

EN

ProTaper Gold®

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
STERILE – STERILIZED BY RADIATION
LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

REUSABLE STERILE ENDODONTIC INSTRUMENTS

1 PRODUCT DESCRIPTION

- Product description:
 ProTaper Gold® are engine-driven tapered endodontic instruments made of nickel titanium alloy followed by a heat treatment. These instruments are used to mechanically shape and prepare the root canals during endodontic therapy. These devices are multiple use root canal instruments available in sterile versions. ProTaper Gold® Shaping files pre-enlarge and optimally shape the coronal 2/3 of the root canal. Finishing files finish the apical 1/3.
- Clinical benefits:
 The clinical benefit for the patient being treated with ProTaper Gold® is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation (shaping and debridement of the root canal).
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

ProTaper Gold® instruments are delivered in sterile blisters.

The ProTaper Gold® range comprises the following endodontic instruments:

Code **xx-length** **y-file size**	**Product**	**ISO XXX YYYv** **XXX-ISO size** **YYY-taper size** **v-variable taper**	**Available lengths**	**Packaging and content**
A 0410 2xx G0y 03	ProTaper Gold® Shaping files	SX: ISO 019 004v	19 mm	Blister with 6 instruments (sterile)
		S1: ISO 018 002v S2: ISO 020 004v	21/25/31 mm	
A 0411 2xx G0y 03	ProTaper Gold® Finishing files	F1: ISO 020 007v F2: ISO 025 008v F3: ISO 030 009v F4: ISO 040 006v F5: ISO 050 005v	21/25/31 mm	Blister with 6 instruments (sterile)
A 0409 2xx G01 03	ProTaper Gold® Assortment	SX S1 S2 F1 F2 F3	19 mm 21/25/31 mm	Blister with 6 instruments (sterile)
A 0409 2xx GD1 03	ProTaper Gold® DEMO PACK / Assortment	SX S1 S2 F1 F2 F3	19 mm 21/25/31 mm	Blister with 6 instruments (sterile)

5 COMPOSITION

- ProTaper Gold® instruments are made of three main components: the working part made of a nickel-titanium alloy, a colored silicone stopper and a plated brass shank with a colored ring (except for instrument SX).
- The files are manufactured with a process that results in a file that has a colored appearance. Due to this proprietary processing, the files may appear slightly curved. This is not a manufacturing defect. While the file can be easily straightened using only your fingers, it is not necessary to straighten the file prior to use. Once inside the canal, the file will follow the root canal anatomy.

6 COMPATIBLE DEVICES

- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).
- Torque control devices ensure optimal usage.
- Use with endodontic motors in a constant rotation at a speed of 300 rpm.
- Set torque at:

File size	Torque [Ncm]	Speed [rpm]
ProTaper Gold® S1 ProTaper Gold® SX	5 - 6	300
ProTaper Gold® S2 ProTaper Gold® F1	1.5 - 2	
ProTaper Gold® F2 ProTaper Gold® F3 ProTaper Gold® F4 ProTaper Gold® F5	3 - 4	

- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicone stopper to the correct working length is the appropriate method of working length determination.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use and the General Processing Instructions For Endodontic Products (see section [9]) to minimize the following risks to the device, the patient and/or the user:
 - Breakage of instrument.
 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Toxic or allergic reactions caused by processing residues.
- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Similar to all machine-driven root canal instruments, ProTaper Gold® instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.

- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible and/or NiTi files (this will help avoid canal transportation).
 - Visually inspect the working part for all the defects listed in the former paragraph during use (i.e. after each wave).
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

1) Before use, check if the instrument is sterile or if processing is needed (refer to Processing Instructions).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
3) Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.
4) Irrigate abundantly and frequently after each instrument used (according to good dental practices).
5) The ProTaper Gold® sequence is the same regardless of the length, diameter or curvature of the canal.

6) Passively scout the coronal 2/3 with hand files ISO 010 and ISO 015, in the presence of a viscous chelator (such as Glyde File Prep™ root canal conditioner).
7) Gently work these instruments until a smooth, reproducible glide path is confirmed. Alternatively, mechanized glide path files (such as ProGlider® or PathFile®) may be used after an ISO 010 hand file.
8) In the presence of NaOCl, brush and follow along the glide path with Shaping File S1. Before light resistance is encountered, laterally brush and cut dentin on the outstroke to improve straight-line access and apical progression. Always brush away from the furcation.
9) Use passively an ISO 015 hand file to resistance.
10) Use Shaping File S1 with brushing action until the same depth as the ISO 015 hand file is reached.
11) Repeat this sequence with Shaping File S2, exactly as described for the Shaping File S1, until the depth of the ISO 015 hand file is reached.
12) In the presence of a viscous chelator or NaOCl, scout the apical 1/3 with ISO 010 and 015 hand files and gently work them until they are loose at length.
13) Establish working length, confirm patency and verify the presence of a smooth reproducible glide path in the apical 1/3.
14) Use first Shaping File S1 and then Shaping File S2, with a brushing action, until working length is reached.
15) Reconfirm working length, irrigate, recapitulate and re-irrigate, especially in more curved canals.
16) Use Finishing File F1, in a non-brushing action, with each insertion deeper than the previous insertion until working length is reached. Do not leave the file at working length for longer than one second.
17) Gauge the foramen with an ISO 020 hand file:
 - If the instrument is snug at length, preparation is finished and the canal is ready to be obturated.
 - If the instrument is loose at length, continue with Finishing File F2 and, when necessary proceed with Finishing Files F3, F4 and F5, with the same non-brushing motion to working length, gauging after each Finishing File respectively with a hand files of the same size (ISO 025, ISO 030, ISO 040 or ISO 050).
18) If necessary, use Shaping File SX with a brushing action to move the coronal aspect of the canal away from furcation and/or to create more coronal shape. SX can also be used to optimally shape canals in shorter roots.
19) Use dedicated paper points to dry and dedicated gutta-percha points to obturate the root canals.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.
- Do not re-use silicone stops. Remove and discard silicone stops after each use.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

Authorised representatives

CE 2797

EU REP
DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP
DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB
United Kingdom

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Dentsply Sirona Europe GmbH