Dentsply Sirona

ProTaper® Universal-Retreatment

EN	Instructions For Use	2
FR	Instructions d'utilisation	8
IT	Istruzioni per l'uso	14
DE	Gebrauchsanweisung	20
ES	Instrucciones de uso	26
PT	Instruções de utilização	32
NL	Gebruiksinstructies	38
SV	Användarinstruktioner	44
DA	Brugsanvisning	50
NO	Brukerinstruksjoner	56
LT	Naudojimo instrukcija	62
LV	Lietošanas norādījumi	68
ET	Kasutusjuhised	74
CS	Návod k použití	80
SK	Návod na použitie	86
HU	Használati utasítás	92
RO	Instrucţiuni de utilizare	98
SR	Uputstva za upotrebu	104
BS	Uputstva za upotrebu	110
CNR	Uputstva za upotrebu	116
SQ	Udhëzimet për përdorim	122
EL	Οδηγίες χρήσης	128
TR	Kullanım Kılavuzu	134
PT-BR	Instruções de uso	140
KO	사용 지침	146

Dentsply Sirona

EN

ProTaper® Universal-Retreatment

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
STERILE – STERILIZED BY RADIATION
NON-STERILE – STERILIZE BEFORE FIRST USE
LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

REUSABLE STERILE AND NON-STERILE ENDODONTIC INSTRUMENTS

1 PRODUCT DESCRIPTION

- Product description:
 ProTaper® Universal-Retreatment files are engine-driven endodontic instrument made of nickel-titanium alloy used to mechanically remove the root canal obturating material while performing retreatment. These devices are multiple use root canal instruments available in sterile and non-sterile versions.
- Clinical benefits:
 The clinical benefits for the patient being treated with ProTaper® Universal-Retreatment files is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal retreatment (shaping and debridement of the root canal).
- Intended patient population:
 Patients who need root canal retreatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

ProTaper® Universal-Retreatment instruments are delivered in both versions: sterile blisters or non-sterile plastic boxes.

The ProTaper® Universal-Retreatment range comprises the following endodontic instruments:

Code	Product	ISO XXX YYYv XXX-ISO size YYY-taper size v-variable taper	Available lengths	Packaging and content
A 1410	ProTaper® Universal-Retreatment	D1: ISO 030 009v D2: ISO 025 008v D3: ISO 020 007v	16/18/22 mm	Blister with 6 instruments (sterile)
				Plastic box with 6 instruments (non-sterile)
	ProTaper® Universal-Retreatment Assortment	D1 D2 D3	16/18/22 mm	Blister with 2 instruments of each ISO (sterile)
				Plastic box with 2 instruments of each ISO (non-sterile)

5 COMPOSITION

ProTaper® Universal-Retreatment instruments are made of three main components: the working part made of a nickel-titanium alloy, a silicone stopper and a plated brass shank with identifier rings.

6 COMPATIBLE DEVICES

- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).
- Torque control devices ensure optimal usage.
- Use with endodontic motors in a constant rotation at a speed of 250 - 300 rpm for removing zinc oxide – eugenol based soluble paste and of 500 rpm for removing gutta-percha or Thermafil obturators.
- Set torque at:

File size	Torque [Ncm]	Speed [rpm]
ProTaper® Universal-Retreatment D1 ProTaper® Universal-Retreatment D2	2.0	500 (for removing gutta percha or Thermafil obturators) 250 - 300 (for removing zinc oxide – eugenol based soluble paste)
ProTaper® Universal-Retreatment D3	1.5	

- The use of radiographs in combination with an apex locator and a tool for adjusting the silicone stopper to the correct working length is the appropriate method of working length determination.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use and the General Processing Instructions For Endodontic Products (see section [9]) to minimize the following risks to the device, the patient and/or the user:
 - Breakage of instrument.
 - Cross-contamination.

- Heat generation due to insufficient lubrication and irrigation.
- Swallowing of working part of the instrument.
- Toxic or allergic reactions caused by processing residues.

- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Similar to all machine-driven root canal instruments, ProTaper® Universal-Retreatment instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.
- Do not use instrument D1 file around a canal curvature.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.

- Use small size, flexible and/or NiTi files (this will help avoid canal transportation).
- Visually inspect the working part for all the defects listed in the former paragraph during use (i.e. after each wave).
- Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

Before removing gutta-percha, carrier-based obturators or paste from a root canal:

1) Before use, check if the instrument is sterile or if processing is needed (refer to Processing Instructions).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
3) Visualize the density of obturation material relative to the width, length, and curvature of the canal.
4) Access the pulp chamber and note the circumferential dimensions of the obturation material at the orifice(s).
5) Select the best removal technique after radiographic and clinical assessment.
6) Without cutting dentin, remove obturation material in a progressive crown-down manner.
7) When the rotary removal method is utilized, select the speed dependent on the obturation material for effectively engaging and removal of the obturation material from the canal.
8) Repeat preparation depending on the obturation material as follow:

For gutta-percha / carrier-based Obturator removal	For soluble Zinc oxide – eugenol based paste removal
Without engaging dentin, gently press the spinning ProTaper® Universal-Retreatment D1 file into the gutta-percha to create friction, generate a heat wave, and extract material out of the canal. Never engage ProTaper® Universal-Retreatment D1 around a canal curvature.	Flood the pulp chamber with the appropriate solvent (DMS IV from Dentsply Sirona) and probe the canal orifice with an explorer to check if the paste has been effectively softened.
In the instance of carrier removal, select the appropriately tapered ProTaper® Universal-Retreatment file that can be carried sufficiently deep into the canal and lateral to the carrier. A long engagement zone will more effectively auger the entire length of the carrier out of the canal.	Without engaging dentin, gently press the spinning ProTaper® Universal-Retreatment D1 file into the material and use a short pecking motion to extract material out of the canal. Never engage ProTaper® Universal-Retreatment D1 around a canal curvature.

9) Continue with the D1 file, or the ProTaper® Universal-Retreatment file that fits passively between the dentinal walls, until material is removed from the coronal one-third of the canal.
10) Select the ProTaper® Universal-Retreatment D2 file and, using one or more passes, extract obturation material from the middle one-third of the canal. Use a brushing outstroke motion to remove material from the canal walls.

11) When appropriate, choose the ProTaper® Universal-Retreatment D3 file and lightly press into the more deeply positioned material and auger obturation material out of the apical one-third of the canal.
12) Continue with the D3 file as long as the flutes of the instrument, upon removal, are loaded with obturation material.
13) When the obturation material is short of the canal terminus, use small sized hand files in the presence of a viscous chelator to negotiate and secure the rest of the canal.
14) After assessing the glide path, select either manual or mechanized shaping files to shape and finish the canal.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.
- Do not re-use silicone stops. Remove and discard silicone stops after each use.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE 2797

Authorised representatives

EU REP

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP

DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH