Dentsply Sirona

WaveOne® Gold

EN	Instructions For Use	2
FR	Instructions d'utilisation	9
IT	Istruzioni per l'uso	17
DE	Gebrauchsanweisung	25
ES	Instrucciones de uso	33
PT	Instruções de utilização	41
NL	Gebruiksinstructies	49
SV	Användarinstruktioner	57
DA	Brugsanvisning	64
NO	Brukerinstruksjoner	71
LT	Naudojimo instrukcija	78
LV	Lietošanas norādījumi	85
ET	Kasutusjuhised	92
CS	Návod k použití	99
SK	Návod na použitie	106
HU	Használati utasítás	114
RO	Instrucţiuni de utilizare	122
SR	Uputstva za upotrebu	130
BS	Uputstva za upotrebu	137
CNR	Uputstva za upotrebu	144
SQ	Udhëzimet për përdorim	151
EL	Οδηγίες χρήσης	158
TR	Kullanım Kılavuzu	166
PT-BR	Instruções de uso	173
KO	사용 지침	180

Dentsply Sirona

EN

WaveOne® Gold

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
SINGLE USE – REPROCESSING PROHIBITED
STERILE – STERILIZED BY RADIATION

INSTRUCTIONS FOR USE

SINGLE USE STERILE ENDODONTIC INSTRUMENTS

1 PRODUCT DESCRIPTION

- Product description:
 WaveOne® Gold are engine-driven tapered endodontic instruments made of nickel titanium alloy.
 WaveOne® Gold Glider files are recommended to be used in endodontic treatment to expand the secured glide path in preparation for the reciprocating mechanized WaveOne® Gold files system.
 WaveOne® Gold instruments are used to mechanically shape and prepare the root canals during endodontic therapy.
 These devices are single use root canal instruments available in sterile version.
- Clinical benefits:
 The clinical benefit for the patient being treated with WaveOne® Gold is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation.
 - Glide path for the WaveOne® Gold Glider files.
 - Shaping and debridement of the root canal for the WaveOne® Gold instruments.
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

WaveOne® Gold instruments are delivered in sterile blisters.

The WaveOne® Gold range comprises the following endodontic instruments:

Code z-quantity xx-length y-file size	Product	ISO XXX YYYv XXX-ISO size YYY-taper size v-variable taper	Available lengths	Packaging and content
BST W1G1 316 OSH	Orifice Shaper for WaveOne® Gold files	ISO 018 004v	16 mm	Blister with 3 instruments (sterile)
BST W1GG zxx 015	WaveOne® Gold Glider files	GLIDER: ISO 015 002v	21/25/31 mm	Blister with 3 instruments (sterile) Blister with 6 instruments (sterile)
A 075z 2xx G0y	WaveOne® Gold shaping files	SMALL: ISO 020 007v PRIMARY: ISO 025 007v MEDIUM: ISO 035 006v LARGE: ISO 045 005v	21/25/31 mm	Blister with 3 instruments (sterile) Blister with 6 instruments (sterile)
A 0754 2xx G9A 03	WaveOne® Gold Assortment	SMALL: ISO 020 007v PRIMARY: ISO025007v MEDIUM: ISO 035 006v LARGE: ISO 045 005v	21/25/31 mm	Blister with 1 instrument of each ISO (sterile)
BST W1GS 325 AST	WaveOne® Gold Sequence	GLIDER: ISO 015 002v PRIMARY: ISO025007v	25 mm	Blister with 1 instrument of each ISO (sterile): K-File 010 GLIDER ISO 015 002 PRIMARY ISO 020 007v
BST W1G0 5xx SQE	WaveOne® Gold Sequence	GLIDER: ISO 015 002v SMALL: ISO 020 007v PRIMARY: ISO025007v MEDIUM: ISO 035 006v LARGE: ISO 045 005v	21/25/31 mm	Blister of 5 instruments with 1 instrument of each ISO (sterile)

5 COMPOSITION

- WaveOne® Gold instruments are made of three main components: the working part made of a nickel-titanium alloy (M-Wire®), a colored silicone stopper and a Nickel-Gold plated brass shank with an Acrylonitrile-Butadiene-Styrene (ABS) colored ring.
- The files are manufactured with a process that results in a file that has a colored appearance. Due to this proprietary processing, the files may appear slightly curved. This is not a manufacturing defect. While the file can be easily straightened using only your fingers, it is not necessary to straighten the file prior to use. Once inside the canal, the file will follow the root canal anatomy.

6 COMPATIBLE DEVICES

- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).
- Use in a reciprocating motion with a designated drive system having the approved WaveOne® / WaveOne® Gold settings.
- The Orifice Shaper for WaveOne® Gold can also be used with RECIPROC® / RECIPROC® Blue motor settings.
- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicone stopper to the correct working length is the appropriate method of working length determination.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use to minimize the following risks to the device, the patient and/or the user:
 - Breakage of instrument.
 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Toxic or allergic reactions caused by processing residues.
- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Similar to all machine-driven root canal instruments, WaveOne® Gold instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.

- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible and/or NiTi files (this will help avoid canal transportation).
 - Visually inspect the working part for all the defects listed in the former paragraph during use (i.e. after each wave).
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

Orifice preparation with the Orifice Shaper for WaveOne® Gold

The Orifice Shaper for WaveOne® Gold may be used when there is restrictive space, or to pre-enlarge the body of a canal, or to relocate the coronal-most aspect of a canal away from an external root concavity.

1) Before use, make sure the instrument is sterile (sterility cannot be guaranteed if packaging is open, damaged or wet).
2) Review different horizontally angulated radiographs, or CBCT images (if required), to diagnostically determine the anatomy, width, estimate length, and curvature of any given root canals.

3) Prepare an adequate access cavity to allow easy identification of the canal orifices.
4) In the presence of a viscous chelator, explore the canal up to a size 010 hand file.
5) Confirm patency and verify a smooth reproducible glide path up to a size 010 hand file.
6) Select WaveOne® Gold motor settings (*) on the endodontic motor.
7) Introduce the Orifice Shaper for WaveOne® Gold instrument into the canal in the presence of irrigant.
8) Advance the Orifice Shaper for WaveOne® Gold in 2-3 gentle amplitudes approximately 2-3 mm in-and-out of the canal. If necessary, use a brushing action to create more coronal space. Irrigate the canal and clean the flutes routinely. Repeat until the access to the root canal is smoother for the following instruments.

(*) the Orifice Shaper for WaveOne® Gold is also compatible with RECIPROC® / RECIPROC® Blue motor settings.

Glide Path creation with WaveOne® Gold Glider

1) Before use, make sure the instrument is sterile (sterility cannot be guaranteed if packaging is open, damaged or wet).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
3) If needed, create a straight-line access to the canal orifice(s) using the Orifice Shaper for WaveOne® Gold, following the instructions above (see *Orifice preparation with the Orifice Shaper for WaveOne® Gold*).
4) Irrigate abundantly and frequently after each instrument used (according to good dental practices).
5) In the presence of a viscous chelator, scout the root canal with hand files, at least an ISO 010, determine the working length by using an electronic apex locator in combination with radiographs, irrigate and confirm patency.
6) With gentle inward pressure, let the WaveOne® Gold Glider file passively progress through any secured region of the canal. Use the WaveOne® Gold Glider file in one or more passes until the full working length is reached.
7) Irrigate the expanded glide path, recapitulate and irrigate and again.
8) Reconfirm the working length, before shaping the canal with the WaveOne® Gold shaping system.

Shaping procedure with WaveOne® Gold

1) Before use, make sure the instrument is sterile (sterility cannot be guaranteed if packaging is open, damaged or wet).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
3) If needed, create a straight-line access to the canal orifice(s) using the Orifice Shaper for WaveOne® Gold, following the instructions above (see *Orifice preparation with the Orifice Shaper for WaveOne® Gold*).
4) Irrigate abundantly and frequently after each instrument used (according to good dental practices).
5) With irrigant in the canal, create and confirm a reproducible glide path using an ISO 010 hand file until full working length is reached.

6) Expand and secure this glide path with at least an ISO 015 file by using either a manual or a dedicated mechanical file, such as PathFile®, ProGlider® or the dedicated WaveOne® Gold Glider file.
7) ALWAYS initiate the shaping procedure with PRIMARY file (025.07 red) in the presence of sodium hypochlorite.
8) Use gentle inward pressure and let PRIMARY file passively progress through any region of the canal that has a confirmed glide path. After shaping 2-3 mm of any given canal, remove and clean the PRIMARY file.
9) Irrigate, recapitulate with an ISO 010 hand file and re-irrigate.
10) Continue with PRIMARY file, in 2-3 passes, to pre-enlarge the coronal 2/3 of the canal.
11) Use a brushing motion on the outstroke to eliminate coronal interferences or to enhance shaping results in canals that exhibit irregular cross-sections.
12) In more restrictive canals, use an ISO 010 hand file, in the presence of viscous chelator, negotiate to the terminus of the canal. Gently work this file until it is completely loose at length.
13) Establish working length, confirm patency and verify the glide path.
14) Expand and secure this glide path with at least an ISO 015 file by using either a manual or a dedicated mechanical glide path file.
15) Carry PRIMARY file to full working length in one or more passes. Upon reaching length, remove the file to avoid over-enlarging the foramen.
16) Inspect the apical flutes:
 - If they are loaded with dentinal debris, then the shape is finished.
 - If PRIMARY file is loose at length with no dentinal debris in the apical flutes, continue shaping with MEDIUM file (035.06 green) and/or LARGE file (045.05 white) until the apical flutes are loaded.
17) If the PRIMARY does not progress, then use SMALL file (020.07 yellow) in one or more passes to working length and then use PRIMARY file to working length to optimize the shape.
18) Use dedicated paper points to dry and dedicated gutta-percha points to obturate the root canals.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- WaveOne® Gold instruments are single use devices and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
- Disinfection, cleaning and sterilization are not applicable.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE 2797

Authorised representatives

EU REP
DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP
DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB
United Kingdom

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Dentsply Sirona Europe GmbH