Dentsply Sirona

Aquasil® Ultra+

Smart Wetting® Impression Material – Putty

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© 2023 Dentsply Sirona Inc. All Rights Reserved. **578636WEB** (R 4/20/23)

Dentsply Sirona

Aquasil® Ultra+

Smart Wetting® Impression Material Putty

Regular Set (Hand Mix) • Soft Fast Set (Hand Mix) • Soft Regular Set (DECA™ delivery)

INSTRUCTIONS FOR USE — ENGLISH

Caution: This is a medical device. For dental use only by dental professionals.
U.S.A.: Rx only.

1. PRODUCT DESCRIPTION

Aquasil® Ultra+ Putty Smart Wetting® Impression Materials are vinyl polysiloxane addition reaction silicone elastomeric dental impression materials with excellent dimensional accuracy, high tear strength, and resistance to permanent deformation.

Aquasil® Ultra+ Material complies with ISO 4823:2021 Elastomeric impression materials.

1.1 Indications for use
Aquasil® Ultra+ Putty materials are thick, vinyl polysiloxane dental impression materials intended to be used in conjunction with more fluid vinyl polysiloxane impression materials in order to make impressions. The resulting impressions are used to make plaster models of the teeth.

1.2 Contraindications
None known.

1.3 Delivery Forms (some delivery forms may not be available in all countries)
Aquasil® Ultra+ Putty is available in:
- Soft or Firm consistency
- Base and Catalyst Tubs for 1:1 hand mixing or 5:1 DECA™ delivery for Dynamic Mixing Machines
- Fast Set or Regular Set

1.4 Composition
Aquasil® Ultra+ Putty is an addition-reaction silicone (vinyl polysiloxane) composed of:
Base: Polyvinyl siloxane, Silicium dioxide, Aluminum hydroxide, Methylhydrogen siloxane, Pigments
Catalyst: Polyvinyl siloxane, Silicium dioxide, Aluminum hydroxide, Organoplatinum complex

1.5 Technical specifications
- Classification: ISO 4823
- Type 0: Putty consistency: Aquasil® Ultra+ Putty (all versions)
- Working Time (22 °C/72 °F).....**Hand Mix Fast Set:** 1'30" minimum
 Hand Mix Regular Set: 2'00" minimum
 Machine Mix Regular Set: 2'00" minimum
 (room temperature mixing and preparation, maximum insertion time)

- Setting Time (Mouth Removal Time)**Hand Mix Fast Set:** 3'00" from start of mix
 Hand Mix Regular Set: 5'00" from start of mix
 Machine Mix Regular Set: 5'00" from start of mix

- % Linear Dimensional Change≤ 1.5
- % Recovery from deformation.....≥ 96.5
- Detail Reproduction.....75 microns

1.6 Compatible tray adhesive
All standard and/or special trays can be used, provided they are rigid and non-deformable. Proper adhesion of the material to the impression tray is essential to avoid errors, especially during removal from the mouth. To ensure a strong bond, make sure that the adhesive is appropriate for the type and material of the tray. (See complete Instructions for Use).

1.7 Compatible Die Stones
- The material is compatible with a range of dental plasters and die stones including Glastone® Type III dental stone and Glastone® 3000 Type V high strength die stone.
- Following impression casting, reprocess reusable impression trays according to tray manufacturer's instructions.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.
This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings
1. Aquasil® Ultra+ Material contains polymerizable monomers and surfactants which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
 - **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - **Avoid ingestion/swallowing of material** to prevent irritation or obstruction. If accidental swallowing occurs, drink lots of water. This material is not hazardous when small quantities are ingested. Larger quantities may cause bowel obstruction. Seek medical attention in the event of digestive irregularities.
2. The product is intended for transient contact. Do not use in applications other than the intended use (e.g., do not use as a temporary reliner) because it may cause irritation.
3. Avoid multiple impressions at one sitting (more than three of the same arch/area) to prevent irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.
4. Aquasil® Ultra+ Material should not be used with patients who have a known hypersensitivity to any of the components.
5. Due to the stiffness of Aquasil® Ultra+ Putty Material, severe undercuts should be blocked out prior to tray insertion. For full arch dentate impressions, select other Aquasil® Ultra+ Material viscosities for tray material.
6. Avoid use on periodontally involved teeth or teeth with compromised attachments, severe undercuts or open furcations to prevent injury to periodontal tissues or tooth extraction. Do not use excessive force to remove impression. If necessary, section tray and material for removal.
7. To avoid periodontal tissue harm, check the integrity of the impression and check that there is no residue of material to be removed from the patient's mouth. Verify that there are no residues of material in the gingival area.

2.2 Precautions
- This product is intended to be used only as specifically outlined in the Instructions For Use. Any use of this product inconsistent with the Instructions For Use is at the discretion and sole responsibility of the practitioner.
- Insufficient data exist to support the use of Aquasil® Ultra+ Material for impressions on exposed bone. Aquasil® Ultra+ Material should only be used on or around intact mucosa.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- The jars cannot be reprocessed. To prevent the jars from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the jars are handled with clean/disinfected gloves.
- As additional precautionary measure, jars may be protected from gross debris but not from all contamination by applying a protective barrier.
- Never interchange different lots of Aquasil® Ultra+ Putty base and catalyst. Unpredictable work and/or set times may result.
- Never interchange measuring spoons or jar lids of base and catalyst. Contamination and partial polymerization of the product may occur. Unpredictable work and/or set times may result.
- Use only clean and properly disinfected and sterilized, or ready to use accessories.
- All products should be used at room temperature. Allow the product to reach room temperature before use. Variations in temperature and mixing speed will affect work time and set time. Higher temperatures and rapid, vigorous mixing will reduce work times and laboratory bench set times (faster) lower temperatures increase them (slower). Clinically, intraoral set time is unaffected within the range of allowable storage and initial mixing temperatures. The material sets faster in the mouth than on the bench. Check impression in the mouth (not on the bench) to be sure it is resilient before removal.
- Do not use if the mix is not homogeneous. To obtain optimum physical properties, the material should be mixed in equal volumes (1:1). The mix should be completely homogeneous (streak free). Base and Catalyst are contrasting colors. When dispensing, observe both components entering the mix tip, and exiting as a homogeneous mix. If not, discard the mixing tip and ensure flow from both ports. Replace mix tip or cartridge and tip and return to dispensing.
- Aquasil® Ultra+ Material should be used with properly fitted, rigid trays with sides for most applications. Trays with retentive features such as Caulk Rim-Lock® Trays or perforations, pretreated with tray adhesive, are recommended. To reduce gagging, choking or swallowing, do not overfill tray.
- Check the integrity of the impression and check that there is no residue of material to be removed from the patient's mouth.
- Do not use the product if package is broken.
- Do not reuse the original cartridge cap. Store the cartridge with the disinfected used mixing tip.
- Two-step technique is NOT recommended for double arch trays.
- **Interactions:**
 - Use of latex or "sulfur-containing" polymeric gloves, some hand soaps and lotions can lead to interference with setting reaction. Contamination may accumulate on the surface of the skin and may not easily wash off. Non-latex gloves worn during use will avoid this possibility. The use of non-latex dental dams is recommended on patients.
 - Astringents based on aluminum or iron salts, or Hydrogen peroxide-containing dental materials may interfere with the setting reaction of vinyl polysiloxanes. Rinse carefully with water before applying the material. Please consult Instructions for Use of astringent to be used.
 - Handle cords with college pliers. Do not touch prepared tooth with gloves. Rinse preparation with copious amounts of water before making the impression.
 - Aquasil® Ultra+ Material should not be intermixed with, or used in conjunction with, other manufacturers' vinyl polysiloxanes or any polyether, polysulfide or conventional (condensation cured) silicones.
 - When used as a primary impression, material should be covered with a polyethylene spacer sheet, and protected from direct contact with oral tissue.

2.3 Adverse Reactions
1. Product may irritate the eyes and skin.
 - **Eye contact:** irritation and possible corneal damage.
 - **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin.
 - **Mucous Membranes:** inflammation
2. Bowel obstruction or other digestive distress may result from ingestion of mixed impression material (see Warnings).
3. Residual material may lead to soft tissue injury or irritation (See Warnings and Step-by-Step instructions).
4. Long term or repeated contact may result in tissue irritation (see Warnings).

2.4 Storage conditions
Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Containers should be tightly closed immediately after use. Jars should be tightly closed by replacing the original lids immediately after use. Jars cannot be processed and must be properly disposed of if contaminated.
- Store away from direct sunlight at temperatures between 5 °C-27 °C/41 °F-81 °F.
- Use the product at room temperature.
- Do not store used cartridges with new (unused) mix tips attached.
- Keep out of sunlight.
- Do not freeze.
- Protect from moisture.
- Do not use after expiration date.
- Do not use if the lot number and/or expiration date are missing or illegible.
- Work time may be affected by temperature. See Step-by-Step Instructions For Use.

3. STEP-BY-STEP INSTRUCTIONS

3.1 DECA™ cartridge loading
NOTE: Heavy Cartridge – Dropping may result in damage to the cartridge rendering it unusable. Never use a damaged cartridge.

The DECA™ cartridge fits duomix™, duomix II, Pentamix®, Dynamix and Volume Mix dynamic mixers*. Use with MixStar dynamic mixer is not recommended. Cartridges are not compatible with 5:1 Pentamix Lite. Check compatibility with the manufacturer instructions.

Carefully remove dispensing cartridge from the package. Also refer to separate dynamic mixer manufacturer's instructions.

3.1.1 Remove the cap from the cartridge by gently turning it counterclockwise and lifting the cap from the cartridge.
3.1.2 Affix the dynamic mixing tip by turning it clockwise to the cartridge. Ensure that the alignment key of the dynamic mixing tip slides into the port of the cartridge. Turn the mix tip until tip is locked in place. No locking ring (or bayonet ring) is required to retain the dynamic mixing tip. Use only Dentsply Sirona mix tips. Other tips may not seal properly.
3.1.3 Insert the cartridge into the dynamic mixer.
3.1.4 Engage machine plunger until the machine's drive shaft is engaged in the mixer drive socket (mixing element is turning). Do not extrude material until ready for use.
3.1.5 When ready to dispense, activate the machine, dispense and discard initial mixed material to ensure homogeneous mix.
NOTE: For first use, always dispense and discard approximately 75 mm of material before filling tray. For all subsequent uses, always dispense and discard approximately 30 mm of material before filling tray.

CAUTION
To reduce risk of incomplete mix – work time or set time discrepancies.
- For first use, always dispense approximately 75 mm of material through the mixing tip and discard.
- For subsequent uses, always dispense approximately 30 mm of material through the mixing tip and discard.
- Material must be bled through the mixing tip (not the cartridge before tip attachment) before loading tray.
- Once dispensing and bleeding through the mixing tip is completed, immediately load tray. Observe material work time starting from cartridge bleeding.

3.1.6 Immediately after activation and bleeding material through the tip, proceed to dispensing material to tray as desired.
3.1.7 Disinfect cartridge with attached used mixing tip. It will work as cap for storage.
3.1.8 Immediately prior to subsequent usage, remove the used mixing tip, wipe end of cartridge, install a new mixing tip, bleed and discard material through the tip and proceed to dispensing.

3.2 Hand Mix (Jars)
3.2.1 Take the 2 components out of the jars using the color-coded dosing spoons.

CAUTION
To reduce risk of incomplete mix – work time or set time discrepancies.
- Never interchange measuring spoons or jar lids of base and catalyst.
- Contamination and partial polymerization of the product may occur.
- Unpredictable work and/or set times may result.

3.2.2 With clean hands, or wearing non-latex gloves, knead identical proportions of Base and Catalyst material together until homogeneous, with no streaks (about 30 sec.).
3.2.3 Place an even layer into the impression tray. Avoid trapping air.

3.3 Impression techniques
Handling: Use only clean and properly disinfected and sterilized, or ready to use accessories.

3.3.1 Putty/Wash Dual Phase one-step ("double-mix") technique
Overview: In this technique both the putty and wash material are simultaneously placed into the mouth. It is recommended that two operators be used for this technique. The timing of events in this technique is critical. A delay inserting the tray loaded with putty material beyond the mixed putty tray **OR** injected wash work time will lead to distortion and ill fitting crowns and bridges.

Tray and wash material mix timing: It is highly recommended to utilize an assistant to assure simultaneous syringing and tray loading.

1. Assure hemostasis and clean the field with air water spray.
2. Using supplied spoons, dispense equal amounts of Aquasil® Ultra+ Putty base and catalyst onto a mixing pad. Close tubs immediately. Do not interchange base and catalyst screw caps or spoons.
 OR
 Prepare DECA™ cartridge and dynamic mixer.
3. Syringe material: Backfill the impression syringe or attach an intraoral tip to the end of the mix tip or prepare digit® cartridge syringe/Dispenser.
4.(a) With clean hands, or wearing non-latex gloves, knead material together (assistant) approximately 30 seconds until a uniform, streak-free color is achieved. Alternatively, material may be mixed on the pad with a stiff spatula. Apply an even layer into tray. Avoid incorporating air. Create a smooth dimple in the material in the area if the preparation(s). **NOTE:** Do NOT apply polyethylene spacer.
 OR
 Dispense material from dynamic mixer into impression tray. Do not overfill tray.
4.(b) Simultaneous with tray loading, syringe (dentist) completely around the preparation/abutment.

Do not use if the mix is not homogeneous. To obtain optimum physical properties, the material should be mixed in equal volumes (1:1). The mix should be completely homogeneous (streak free). Base and Catalyst are contrasting colors. When dispensing, observe both components entering the mix tip, and exiting as a homogeneous mix. If not, discard the mixing tip and ensure flow from both ports. Replace mix tip or cartridge and tip and return to dispensing.

CAUTION
To reduce risk of incomplete mix – tissue irritation or compromised properties.
- Observe both base and catalyst entering mix tip.
- If not, remove and replace tip, or replace cartridge and tip.

Tray and wash material insertion timing – Loaded tray must be inserted within the intraoral work time of first syringing of wash material intraorally.

5. Insert filled tray into the mouth.
6. Retain impression in position until firmly set. The material sets faster in the mouth than on the bench. Check impression in the mouth (not on the bench) to be sure it is firm, resilient, and non-tacky before removal.
 - Minimum Removal Time (from the start of mix)
 Fast Set: 3 minutes
 Regular Set: 5 minutes

Observe minimum removal time of selected wash material if greater than Aquasil® Ultra+ Putty minimum removal times.

7. Carefully release the impression and remove.
8. Check for residual material in undercut areas.

CAUTION
To reduce risk of Residual material – Periodontal irritation.
- Undercuts should be blocked out prior to use.
- Check the integrity of the impression.
- Verify that there is no material residue.

9. Rinse impression thoroughly under running water. Follow Hygiene processing procedures below before shipping and/or use.

3.3.2 Dual Phase two-step ("putty-wash") technique with a primary tray impression
Overview: The two step putty wash technique separates the impression technique into two steps.
- Formation of the custom putty wash with relief and;
- The injection of wash material around the prepared teeth and placement of wash material into the custom putty tray for seating into the mouth.

1. Using supplied spoons, dispense equal amounts of Aquasil® Ultra+ Putty base and catalyst onto a mixing pad. Close tubs immediately. Do not interchange base and catalyst screw caps or spoons.
 OR
 Prepare DECA™ cartridge and dynamic mixer
2. With clean hands, or wearing non-latex gloves, knead material together (assistant) approximately 30 seconds until a uniform, streak-free color is achieved. Alternatively, material may be mixed on the pad with a stiff spatula. Apply an even layer into tray. Avoid incorporating air.
 OR
 Dispense material from dynamic mixer into impression tray. Do not overfill tray.
3. Place one or more polyethylene plastic sheets over the impression material surface. This sheet creates the uniform spacing required for this technique. It also protects the surface from contamination.
4. Insert filled tray into the mouth. Move tray slightly from side to side and back and forth. This creates the spacing for the wash material.
5. Remove tray from the mouth when material has reached a rubbery consistency and set aside. Do not disturb impression until fully set.
6. Remove polyethylene sheet, rinse impression with water and air dry.
7. It is recommended to utilize an assistant and two cartridge guns, or one cartridge gun and digit® Targeted Delivery System Unit Dose syringe or digit Power® Dispenser to assure simultaneous syringing and tray loading.
8. Follow impression steps above, using wash material for both intraoral syringing and loading into primary tray impression.
9. Carefully release the impression and remove. Observe the minimum removal time of the light-bodied material selected.
10. Check for residual material in undercut areas.

CAUTION
To reduce risk of Residual material – Periodontal irritation.
- Undercuts should be blocked out prior to use.
- Check the integrity of the impression.
- Verify that there is no material residue.

11. Rinse impression thoroughly under running water. Follow Hygiene processing procedures below before shipping and/or use.

3.3.3 Border Molding Technique
1. Dispense material on to the borders of the tray.
2. Insert the tray and perform conventional border molding techniques.
3. When material has reached a firm consistency, remove tray and examine borders. The addition of a small mix can be used if voids or deficient areas are present. Over-extensions can be easily removed with a scalpel blade.
4. The final impression can be made with Aquasil® Ultra+ XLV, LV, Medium, or Heavy material.
5. Insert filled tray into the mouth.
6. Observe the Minimum Removal Time of the selected Aquasil® Ultra+ final impression material.
7. Rinse impression thoroughly under running water. Follow Hygiene processing procedures below before shipping and/or use.

3.3.4 Impression Technique: duplication of models
1. Load the impression tray with Aquasil® Ultra+ Putty material.
2. It is suggested to simultaneously syringe, using Aquasil® Ultra+ wash material, completely around the model teeth.
3. Seat model into the loaded tray. Retain impression in position until firmly set.
4. Remove model from impression by pulling slowly to break seal.

* Pentamix®, Dynamix and Volume Mix dynamic mixers are not registered trademarks of Dentsply Sirona.

4. HYGIENE, PROCESSING AND DISPOSAL

Cross-contamination.
Infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Protective barriers are designed for single use and must be disposed of after each use in accordance with local regulations. The barrier is not a replacement for applicable cleaning, disinfection, and sterilization.
- Do not attempt to refill or consolidate used cartridges. Dispose of empty cartridges in accordance with local regulations.
- Reprocess reusable products as described below.
- Refer to Dynamic Mixer manufacturer's instructions for dispenser cleaning and disinfection instructions.
- Used cartridge and mix tips may be contaminated and should be processed as described below before storage. Keep used, processed mix tip on cartridge as a cap.
- Process impression as described below before shipping and/or use.
- The jars and dispensing spoons cannot be reprocessed.
- To prevent the jars and spoons from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the cartridges, tubs and spoons are handled with clean/disinfected gloves.
- Dispose of contaminated jars, spoons in accordance with local regulations.
- Do not reuse jars or spoons if contaminated.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- Dispensing materials with clean/disinfected gloves or hands or suitable transfer forceps in a separate room, bringing only what will be used into the operatory is strongly recommended.
- It is advisable to wrap disposable cartridges and automatic mixers in cellophane after each disinfection stage.

Incidental contact of the jars and dispensing spoons with water, soap or a water-based hospital-level disinfection solution will not damage jars or spoons. Do not allow any solution contact with contained material. Discard material that has been in contact with any fluid or non-sterile instrument.
Repeated liquid contact may damage label. Dry the jars and dispensing spoons with a lint-free single-use cloth.

NOTE: Vigorous wiping can destroy the label. Wipe cartridge gently.

4.1 Cartridge/Mix Tip

Process		**✔ Acceptable**
Cleaning	Mechanical (Ultrasonic)	
	Manual	✔
	Automated (Washer/Disinfector)	
Disinfection	Automated (Washer/Disinfector)	
	Manual	✔
	Immersion	
Sterilization	Steam Autoclave	
	Immersion	
Single Use	Dispose of after Use	✔
Reusable/Not Reprocessable	Dispose of if contaminated	

4.2 VPS Bite Registration Impression

Process		**✔ Acceptable**
Cleaning	Mechanical (Ultrasonic)	
	Manual	✔
	Automated (Washer/Disinfector)	
Disinfection	Automated (Washer/Disinfector)	
	Manual	
	Immersion	✔
Sterilization	Steam Autoclave	
	Immersion	
Single Use	Dispose of after Use	
Reusable/Not Reprocessable	Dispose of if contaminated	

For complete processing instructions and validated parameters, please refer to the "50 ml Dispenser-Cartridge-VPS Bite Registration/Impression Hygiene" Instructions, which is available on our webpage at www.dentsplysirona.com/en/service-contact/ifu.html. If requested, we will send you a free printed copy of the Instructions for Use in the language you require within 7 days. In the USA, call 1-800-532-2855. Outside the USA, use the order form provided on www.dentsplysirona.com/en/service-contact/ifu.html#orderform-ifu for this purpose.

Unset Aquasil® Ultra+ Material can be removed from clothing with dry cleaning solvent.

4.3 Casting
- The impression should be removed from the disinfectant, rinsed with water, and lightly air dried.
- The impression may be poured in plaster or stone 60 minutes after disinfection or delayed up to fourteen days. For epoxy dies, delay pour 60 minutes.
- The material is compatible with a range of dental plasters and die stones including Glastone® Type III dental stone and Glastone® 3000 Type V high strength die stone. Impressions may be silver or copper-plated.
- Store the impressions at ambient temperature, in a dry place away from direct sunlight.
- If the impression is to be shipped, use suitable packaging to prevent distortion.
- Impression scanning can be performed without the use of opaque sprays using specific scanners. Follow the instructions of the scanner manufacturer.

4.4 Disposal
- Based on the information in the SDS, the waste generated by the disposal of the product does not present characteristics of physical danger, for the person or the environment; it can therefore be managed, according to local regulations, as a non-hazardous waste.
- Contaminated containers, accessories and finished impressions that have been in contact with the patient must be cleaned and disinfected as outlined in Hygiene above prior to disposal or disposed of as hazardous waste at risk of biological contamination.
- Always wear gloves when handling packaging, accessories and finished impressions.

5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

1. Do not use after expiration date. Do not use if the lot number and/or expiration date are missing or illegible. ISO standard uses: "YYYY-MM-DD."
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date
3. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
4. Device Identification (Basic UDI-DI): ++D002IMPRESSIONMATSGC

© 2023 Dentsply Sirona Inc. All Rights Reserved. **578636WEB** (R 4/20/23)

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