Dentsply Sirona

TruNatomy®

EN	Instructions For Use	2
FR	Instructions d'utilisation	8
IT	Istruzioni per l'uso	15
DE	Gebrauchsanweisung	22
ES	Instrucciones de uso	29
PT	Instruções de utilização	36
NL	Gebruiksinstructies	43
SV	Användarinstruktioner	50
DA	Brugsanvisning	56
NO	Brukerinstruksjoner	62
LT	Naudojimo instrukcija	68
LV	Lietošanas norādījumi	74
ET	Kasutusjuhised	80
CS	Návod k použití	86
SK	Návod na použitie	92
HU	Használati utasítás	98
RO	Instrucţiuni de utilizare	105
SR	Uputstva za upotrebu	111
BS	Uputstva za upotrebu	117
CNR	Uputstva za upotrebu	123
SQ	Udhëzimet për përdorim	129
EL	Οδηγίες χρήσης	135
TR	Kullanım Kılavuzu	142
PT-BR	Instruções de uso	148
KO	사용 지침	155

Dentsply Sirona

EN

TruNatomy®

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
SINGLE USE – REPROCESSING PROHIBITED
STERILE – STERILIZED BY RADIATION

INSTRUCTIONS FOR USE

SINGLE USE STERILE ENDODONTIC INSTRUMENTS

1 PRODUCT DESCRIPTION

- Product description:
 TruNatomy® instruments are engine-driven tapered endodontic instruments made of nickel titanium alloy. These instruments are used to mechanically shape and prepare the root canals during endodontic therapy. These devices are single-use root canal instruments available in sterile version.
- Clinical benefits:
 The clinical benefit for the patient being treated with TruNatomy® is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation (glide path, shaping and debridement of the root canal).
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

TruNatomy® instruments are delivered in sterile blisters.

The TruNatomy® range comprises the following endodontic instruments:

Code **xx-length** **yyy-file size**	**Product**	**ISO XXX YYYv** **XXX-ISO size** **YYY-taper size** **v-variable taper**	**Available lengths**	**Packaging and content**
BST TNMY zxx GLR	TruNatomy® Glider	ISO 017 002v	21/25/31 mm	Blister with 3 instruments (sterile) Blister with 6 instruments (sterile)
BST TNMY zxx OMD	TruNatomy® Orifice Modifier	ISO 020 008	16 mm	Blister with 3 instruments (sterile) Blister with 6 instruments (sterile)
BST TNMY zxx yyy	TruNatomy® files	Small: ISO 020 004v Prime: ISO 026 004v Medium: ISO 036 003v Large: ISO 046 002v	21/25/31 mm	Blister with 3 instruments (sterile) Blister with 6 instruments (sterile)
BST TNMY 5xx AST	TruNatomy® Assortment	ISO 017 002v ISO 020 008 ISO 020 004v ISO 026 004v ISO 036 003v	21/25/31 mm	Blister with 1 instrument of each ISO (sterile)
BST TNMY 6xx AST	TruNatomy® Assortment	ISO 017 002v ISO 020 008 ISO 020 004v ISO 026 004v ISO 036 003v ISO 046 002v	21/25/31 mm	Blister with 1 instrument of each instrument (sterile)
BST TNMY 3xx SQC	TruNatomy® Sequence	ISO 017 002v ISO 020 008 ISO 026 004v	21/25/31 mm	Blister with 1 instrument of each ISO (sterile)

5 COMPOSITION

- TruNatomy® instruments are made of three main components: the working part made of a nickel-titanium alloy, a colored silicone stopper and a plated brass shank with a colored ring (except for Orifice Modifier).
- The files are manufactured with a process that results in a file that has a colored appearance. Due to this proprietary processing, the files may appear slightly curved. This is not a manufacturing defect. While the file can be easily straightened using only your fingers, it is not necessary to straighten the file prior to use. Once inside the canal, the file will follow the root canal anatomy.

6 COMPATIBLE DEVICES

- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).
- Torque control devices ensure optimal usage.

- Use with endodontic motors in a constant rotation at a speed of 500 rpm.
- Set torque at: 1.5 Ncm.
- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicone stopper to the correct working length is the appropriate method of working length determination.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use to minimize the following risks to the TruNatomy® devices, the patient and/or the user:
 - Breakage of instruments.
 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Toxic or allergic reactions caused by processing residues.
- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Similar to all machine-driven root canal instruments, TruNatomy® instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.

- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible and/or NiTi files (this will help avoid canal transportation).
 - Visually inspect the working part for all the defects listed in the former paragraph during use (i.e. after each wave).
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

1) Before use, make sure the instrument is sterile (sterility cannot be guaranteed if packaging is open, damaged or wet).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
3) Estimate the working length using well-angulated preoperative radiographs.
4) Prepare a conservative access cavity sufficient enough to reveal all root canal orifices.
5) Scout the coronal 2/3 of canals with an ISO 010 K-File in the presence of lubricant such as ProLube® and Glyde File Prep™ and irrigate.
6) Followed by a TruNatomy® Orifice Modifier at 500 rpm and 1.50 Ncm. With irrigant in the canal advance TruNatomy® Orifice Modifier in 2-3 gentle amplitudes approximately 2-5 mm in-and-out of the canal. Repeat until the coronal third is shaped. The instrument has 7 mm of cutting flutes, insertion into the canal should not exceed the length of the cutting portion. Irrigate the canal and clean cutting flutes routinely.
7) Scout the whole root canal with an ISO 010 K-File, determine working length by using an electronic apex locator (EAL) in combination with radiographs, irrigate and confirm patency.

8) With irrigant in the canal, create and confirm a reproducible glide path using a TruNatomy® Glider in 2-3 gentle amplitudes approximately 2-5 mm. Irrigate and repeat until previously confirmed working length with an EAL has been reached.
9) ALWAYS begin shaping with TruNatomy® PRIME (500 rpm / 1.5 Ncm) file passively in the presence of sodium hypochlorite with no more than 2-3 gentle amplitudes approximately 2-5 mm in-and-out of the canal. Irrigate and repeat as necessary until full working length is reached. Upon reaching length, remove the file to avoid over-enlarging the apical canal portion.
10) Routinely irrigate the canal and clean the instruments cutting flutes of debris upon removal.
11) If TruNatomy® PRIME does not progress easily, remove, irrigate, and recapitulate with an ISO 010 K-File to confirm canal patency and move to the TruNatomy® SMALL.
12) Inspect cutting flutes routinely upon removal for presence of unwinding and straightening. If deformation is noted, discard and use a new TruNatomy® instrument.
13) Advance TruNatomy® SMALL passively in the presence of sodium hypochlorite with no more than 2-3 gentle amplitudes approximately 2-5 mm in-and-out and remove instrument. Irrigate and repeat as necessary to working length in a gentle / passive in-and-out motion and then use TruNatomy® PRIME to working length to optimize the shape.
14) If TruNatomy® PRIME is loose at length with no dentinal debris in the apical flutes, continue shaping with TruNatomy® MEDIUM and subsequently with LARGE if needed. Care should be taken to guide each instrument gently to full working length to ensure completeness of shape. Use TruNatomy® LARGE in larger and straighter canals only, such as maxillary central incisors, some palatal or distal canals of molars.
15) Following proper cleaning and shaping, confirm the final working length measurement using the last TruNatomy® instrument by hand. Verifying the shape at working length and passive fit ensures that TruNatomy® Conform Fit™ Gutta-Percha Points matches the correct apical canal diameter. This is achieved by passively placing the last TruNatomy® instrument in the canal.
 ⇨ If the instrument reaches the working length passively, choose and try-in the matching TruNatomy® Conform Fit™ Gutta-Percha Points.
 ⇨ If the fit is not passive, instrument (as described above) with the last TruNatomy® shaping instrument after reconfirming working length, irrigating and confirming patency.
16) Once appropriateness of the shape is confirmed, proceed with 3-D disinfection protocols.
17) Use dedicated paper points to dry and dedicated gutta-percha points to obturate the root canals.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- TruNatomy® instruments are single use devices and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
- Disinfection, cleaning and sterilization are not applicable.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE 2797

Authorised representatives

EU REP
DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP
DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH