Dentsply Sirona

EasyPost®

EN Instructions For Use 2
FR Instructions d'utilisation 9
IT Istruzioni per l'uso 16
DE Gebrauchsanweisung 23
ES Instrucciones de uso 30
PT Instruções de utilização 37
NL Gebruiksinstructies 44
SV Användarinstruktioner 51
DA Brugsanvisning 58
NO Brukerinstruksjoner 65
LT Naudojimo instrukcija 72
LV Lietošanas norādījumi 79
ET Kasutusjuhised 86
CS Návod k použití 93
SK Návod na použitie 100
HU Használati utasítás 107
RO Instrucţiuni de utilizare 114
SR Uputstva za upotrebu 121
BS Uputstva za upotrebu 128
CNR Uputstva za upotrebu 135
SQ Udhëzimet për përdorim 142
EL Οδηγίες χρήσης 149
TR Kullanım Kılavuzu 156
PT-BR Instruções de uso 163
KO 사용 지침 170

Dentsply Sirona

EN

EasyPost®

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
NON-STERILE – MANDATORY PROCESSING BEFORE FIRST USE
SINGLE USE – REPROCESSING PROHIBITED
REUSABLE – LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

POSTS: SINGLE USE NON-STERILE IMPLANTABLE ENDODONTIC DEVICES
REF. C0600 (EasyPost® - Post)
DRILLS FOR POSTS: REUSABLE NON-STERILE ENDODONTIC INSTRUMENTS
REF. C0601 (EasyPost® - Precision Drill)
POSTS SETS[1]
REF. C0602 (EasyPost® Kit Ø X.X mm)
[1.] Posts and Drills for Posts instruments can be sold in sets with other instruments and components. See sections [4] and [17] for details.

1 PRODUCT DESCRIPTION

- Product description:

 EasyPost®
 EasyPost® are single use and non-sterile radiopaque, white fiber post. They are made of glass fibers embedded in an epoxy resin matrix. Their shape is cylindric-conical. These devices have a color code on their upper extremity according to their dimensions. Fiber Posts are available in a full range of sizes (see section [4]).

 Drills for Posts
 EasyPost® precision drills are Stainless steel reusable instruments made of a working part for precision drilling, a shank to be connected to a contra-angle, and a color coding to identify the corresponding post.

- Clinical benefits:
 The clinical benefit for the patient being treated with EasyPost® is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purposes:

- Posts are devices intended to retain the core material in a tooth with extensive loss of coronal structure.
- Drills for Posts instruments are devices used for the creation of hole for post placement.
- Post sets are devices intended to place a post for post-endodontic restoration.

- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

Insufficient residual dentin: at least 2 mm of tooth structure is required around the preparation.

4 DELIVERY FORM

The EasyPost® brand comprises the following devices:

Posts

SKU	Product	Available sizes					Packaging and content
		N°X: Color	Anchor		Head		
			Diameter (mm)	Length (mm)	Diameter (mm)	Length (mm)	
C060000 000X00	EasyPost® - Post	N°1: Yellow	0.8	6.23	1.35	13.77	Blister with 10 devices
		N°2: Red	0.8	7.48	1.47	12.52	
		N°3: Blue	1.0	7.49	1.67	12.51	
		N°4: Green	1.0	9.29	1.83	10.71	
		N°5: Black	1.3	8.37	2.04	11.63	
		N°6: White	1.3	10.37	2.22	9.63	

Drills for Posts

SKU	Product	Available sizes			Packaging and content
		N°X: Color	Ø	Working part length	
C060100 000X00	EasyPost® - Precision Drill	N°1: Yellow	Ø 0.85 mm	6.75 mm	Box with 3 instruments
		N°2: Red	Ø 0.85 mm	8.00 mm	
		N°3: Blue	Ø 1.05 mm	8.00 mm	
		N°4: Green	Ø 1.05 mm	9.75 mm	
		N°5: Black	Ø 1.35 mm	8.75 mm	
		N°6: White	Ø 1.35 mm	10.75 mm	

Posts Sets

- EasyPost® devices can also be sold in sets with Largo® Peeso Reamer instruments. The different configurations of sets are detailed in the section [17].

5 COMPOSITION

EasyPost®

EasyPost® fiber posts are made of glass fibers (81% ± 3% weight) embedded in an epoxy resin matrix (19% ± 3% weight).

Drills for Posts

EasyPost® precision drills are made of Stainless steel.

6 COMPATIBLE DEVICES

- EasyPost® Fiber Posts are used in combination primer adhesive combined with dual cure resin cement and following devices: Largo® Peeso Reamer and EasyPost® Precision Drill.

7 PERFORMANCE CHARACTERISTICS

- Expected life of Posts: 10 years.
- Posts radio-opacity level: 200% Al equivalent (ISO 4049: Dentistry - Polymer-based restorative materials).

8 SAFETY NOTES

8.1 Warnings

- Strictly follow the Instructions For Use to minimize the following risks to the EasyPost® devices, the patient and/or the user:
 - Breakage of instrument.
 - Contamination and cross-contamination.
 - Swallowing of working part of the instrument.
 - Post debonding.
 - Toxic or allergic reactions caused by processing residues.
- The use of inappropriate rotation speed, pressure and contact duration would speed up the wear of drills for posts, increase the risk of breakage and heat creation.
- Please refer to the Largo® Peeso Reamer brand IFU for the Warnings related to these instruments.

8.2 Precautions

General

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Processing step is mandatory before use according to section [10] for all the products mentioned in this IFU.

Posts

- Check the Posts devices before use, do not use and discard the devices if damaged.
- Shortening of the Posts devices should be done outside the mouth using a diamond-tipped disc.
- Do not touch the Posts devices with your fingers after cleaning / sterilization.

Drills for Posts

- Check the Drills for Posts instruments before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended purpose with the required safety level, instruments should be discarded.
- Check the Drills for Posts instruments and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended purpose with the required safety level, instruments should be discarded.
- Before using any Drills for Posts instrument, make sure it is it is well connected to the contra-angle head.
- Always use minimal apical pressure. Never force the Drills for Posts instruments down the canal.

Posts Sets instruments

- Please refer to the Largo® Peeso Reamer brand IFU for the Precautions related to these instruments.

8.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

8.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

Keep away EasyPost® devices from light and heat.

8.5 MRI (Magnetic Resonance Imaging) safety information

Fiber Posts are MR Safe. A patient with EasyPost® devices can be safely scanned without specific conditions.

9 STEP BY STEP INSTRUCTIONS

1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
2) Place the transparent drawing (with full sizes outlines of the different sized posts) over the radiography to select the appropriate EasyPost® - Post size:
 ⇨ At least 4 mm of root canal filling should remain in the apical region. In case of curved canals, the length must be reduced.
 ⇨ Select the post size so that 2/3 of its length will be inserted in the root canal.
3) Select the appropriate Largo® Peeso Reamer and EasyPost® - Precision Drill sizes according to the canal anatomy and the selected EasyPost® - Post size as indicated in the following table:

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4) Remove the gutta percha and perform penetration drilling with the selected Largo® Peeso Reamer size. Drill as deep as necessary for the post to be inserted at 2/3 of its length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canals, the length must be reduced.
5) Perform precision drilling with the selected EasyPost® - Precision Drill size.
6) Rinse the canal with water or air to remove drilling residues.
7) Check the EasyPost® - Post positioning in mouth.
8) Shorten the selected EasyPost® - Post outside the mouth to its final length using a diamond-tipped disc.
9) Immerse or soak the selected EasyPost® - Post with isopropanol for 30 seconds. Let the post dry for 1 minute.
10) For dentin conditioning prior to post cementation, refer to adhesive or cement manufacturer's instructions for use if applicable.
11) Insert the sealing adhesive or cement in the root canal and on the post itself following adhesive or cement manufacturer´s instructions for use[2].
12) Insert the EasyPost® - Post in the root canal down to 2/3 of its length.
13) Wait for the cement to harden, then remove any excess.

[2.] The adhesion between the Posts devices and the Prime&Bond Active with SCA and Core-x Flow has been tested.

10 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- EasyPost® devices are single use implantable devices, re-use and processing more than 1 time is prohibited. Re-use can increase the risk of cross contamination, poor cementation or breakage.
- Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.

11 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

12 EXPIRATION DATE

None.

13 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

14 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for processing. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

15 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Information will be available at both the European database on medical devices (EUDAMED) and the IFU website, where it will be linked to the Basic UDI-DI.

EUDAMED

https://ec.europa.eu/tools/eudamed

IFU

dentsplysirona.com/ifu

The Basic UDI-DI numbers are:

- Posts Fibers: ++J00310056E8
- Post Sets: ++J00310126E4

16 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

17 POSTS SETS

The different configurations of Posts Sets are detailed in the following table:

SKU	Product	XX: Diameter	Packaging and content
C06020XX00100	EasyPost® Kit Ø X.X mm	08: Ø 0.8 mm	Box with: 1 Largo® Peeso Reamer 1 1 EasyPost® Precision Drill 1 1 EasyPost® Precision Drill 2 12 EasyPost® Post 1 12 EasyPost® Post 2
		10: Ø 1.0 mm	Box with: 1 Largo® Peeso Reamer 2 1 EasyPost® Precision Drill 3 1 EasyPost® Precision Drill 4 12 EasyPost® Post 3 12 EasyPost® Post 4
		13: Ø 1.3 mm	Box with: 1 Largo® Peeso Reamer 3 1 EasyPost® Precision Drill 5 1 EasyPost® Precision Drill 6 12 EasyPost® Post 5 12 EasyPost® Post 6

CE 2797

Authorised representatives

EU REP

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP

DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH