Dentsply Sirona

Glyde File Prep™

EN	Instructions For Use	2
FR	Instructions d'utilisation	8
IT	Istruzioni per l'uso	14
DE	Gebrauchsanweisung	20
ES	Instrucciones de uso	26
PT	Instruções de utilização	32
NL	Gebruiksinstructies	38
SV	Användarinstruktioner	44
DA	Brugsanvisning	50
NO	Brukerinstruksjoner	56
LT	Naudojimo instrukcija	62
LV	Lietošanas norādījumi	68
ET	Kasutusjuhised	74
CS	Návod k použití	80
SK	Návod na použitie	86
HU	Használati utasítás	92
RO	Instrucţiuni de utilizare	98
SR	Uputstva za upotrebu	104
BS	Uputstva za upotrebu	110
CNR	Uputstva za upotrebu	116
SQ	Udhëzimet për përdorim	122
EL	Οδηγίες χρήσης	128
TR	Kullanım Kılavuzu	134
PT-BR	Instruções de uso	140
KO	사용 지침	146

Dentsply Sirona

EN

Glyde File Prep™

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
NON-STERILE
DO NOT USE AFTER EXPIRATION DATE
DO NOT USE IF THE PACKAGE IS DAMAGED

INSTRUCTIONS FOR USE

GLYDE FILE PREP™ (REF. A0901x, A0902x, A0903x, A0905x)

1 PRODUCT DESCRIPTION

- Product description:
 Glyde File Prep™ is a product in the form of a gel that allows for cleansing action that facilitates easy removal of vital pulp tissue and necrotic pulp tissue during endodontic treatment.
 The gel may be used in two configurations a) in single use doses or b) in syringes connected to dedicated tips.
- Clinical benefits:
 The clinical benefit for the patient is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Products in the form of a gel for use in the chemical and mechanical cleansing of the root canal preparation during endodontic therapy.
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 These products are intended to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

The Glyde File Prep™ range comprises the following endodontic devices:

Code	Product	Packaging and content		Single Use / Multiple Use
A09 0100 0000 00	Glyde File Prep™ - Syringe Kit	Pack with 3 syringes of 3 ml		Multiple Use
A09 0200 0000 00	Glyde File Prep™ - Intro Package	Pack with	2 syringes of 3 ml	Multiple Use
			25 Tips	Single Use
A09 0300 0000 00	Glyde File Prep™ - Syringe tips	Pack with 50 Tips		Single Use
A09 0500 0000 00	Glyde File Prep™ - Doses pack	Box with 35 ampoules of 0.5 ml		Single Use

5 COMPOSITION

Glyde File Prep™ Root Canal Conditioner is a combination of:

- Ethylenediaminetetraacetic Acid (EDTA).
- Carbamide peroxide in a water-soluble base.
- Lubricating agents.

The EDTA chelates calcium salts from calcified areas within the root canal to facilitate the use of files.

6 COMPATIBLE DEVICES

- Glyde File Prep™ is not intended to be used with any accessories.
- Glyde File Prep™ is designed to be used in combination with:
 - Sodium Hypochlorite irrigation solution (NaOCl).
 - Shaping Files.
 - Disposa-Shield™ Infection Control Barrier.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use and the General Processing Instructions For Endodontic Products (see section [9]) to minimize the following risks to the device, the patient and/or the user:
- The gel shall not be applied to tissue outside of the root canal or on soft tissues such as mucosa or skin.
- Contact of gel on unwanted tissues (eyes, skin etc.) can cause irritation, rash or damage.
 - Immediately wash skin with soap and water, or flush eyes with water.
- Swallowing of the gel can cause a toxic or allergic reaction.
 - Immediately drink plenty of water and provide fresh air.
- In all cases (ingestion, contact with eyes or skin), please contact a doctor immediately.
- Extrusion of gel beyond the apex can cause necrosis of the tissue and may result in eventual tooth extraction.

- Improper use of the gel (deviation from Instructions For Use) can result in ineffective lubrication and/or ineffective cleansing of root canal, thus increasing risk of infection and eventual tooth extraction.
- Not following good dental practices for prevention of contamination and cross-contamination increases the risk of serious infection to the patient and staff.
- Do not use an irrigation solution other than sodium hypochlorite in presence of Glyde File Prep™. Using another irrigation solution may have unintended consequences and may affect the treatment results.

7.2 Precautions

- Review all corresponding instructions and precautions before beginning treatment.
- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use.
- Check the devices before each use for signs of defects such as deformations, breakage or leakage. With these indications the devices are not able to fulfil the intended use with the required safety level, devices should be discarded.
- Do not use an irrigation solution other than sodium hypochlorite. Using another irrigation solution may have unintended consequences and may affect the treatment results.
- A fresh Disposa-Shield (A88013) must be used together with the syringe prior to each use.
- Use only Glyde File Prep™ Syringe Tip with the Glyde File Prep™ Syringe. Use of other tips may cause incorrect application of the gel.
- Do not dispense the gel into the canal directly with the syringe.
- Do not inject gel under pressure into the canal. The gel should extrude easily. DO NOT USE EXCESSIVE FORCE.
- The use of a dappen dish is the appropriate method to dispense the gel and coat the file with the gel.

7.3 Adverse reactions / Undesirable side effects

- Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

⚠ CAUTION

Causes skin irritation (H315)

⚠ CAUTION

Causes serious eye irritation (H319)

⚠ CAUTION

May cause drowsiness or dizziness (H336)

⚠ CAUTION

Harmful if swallowed (H302)

7.4 Storage conditions / Handling

- Close the syringe tightly after each use with the original syringe cap tightly to avoid leakage and evaporation.
- The gel shall be stored between 2°C to 24°C.
- Shelf-life of the gel is 24 months.
- Glyde File Prep™ must be kept away from sunlight and heat.
- Refrigeration has an influence on the viscosity of the gel as it is affected by storage conditions. Store in a refrigerator, if a more viscous or thicker consistency is desired.
- For expiry date, please consult the label on packaging.

8 STEP BY STEP INSTRUCTIONS

1) Open the pulp chamber in the usual manner and obtain adequate access to the root canals.
2) Remove the pulp and flush canal with sodium hypochlorite solution.
3) Dispense Glyde File Prep™.
 - Either directly from the unit dose applicator (or from the syringe with the Glyde File Prep™ syringe tip in place) in and around the canal orifice,
 - Or from the applicator (or from the syringe with the Glyde File Prep™ syringe tip in place) to a dappen dish and coat a glide path file with the gel.
4) Use one or several glide path files and start to explore the canal to define the working length.
5) For handling instruction of the instrument, please consider the corresponding IFU of the used instruments.
6) Irrigate abundantly and frequently the canal with sodium hypochlorite throughout the procedure and after each instrument used (according to good dental practices).
7) The sodium hypochlorite will react with the Glyde File Prep™ and begin to effervesce or bubble. This action allows for pulp tissue, dentinal shavings and debris to be easily removed.
 - PLEASE TAKE NOTE OF THE BUBBLING REACTION.
8) Enhance this reaction by mixing the Glyde File Prep™ with the sodium hypochlorite solution by using a small hand instrument.
9) Continue to use sodium hypochlorite irrigation solution until all evidence of bubbling has ceased.
10) Finalize the endodontic procedure according to the IFU of the used instrument.

11) For mechanized shaping techniques coat the first file in the rotary file technique with Glyde File Prep™ and use to desired depth. Glyde File Prep™ will now act as a file lubricant while instrumenting the canal with a rotary instrument.
 - NOTE: irrigate and evacuate using sodium hypochlorite to remove Glyde File Prep™ mixture and dentinal debris from canal after every 2 files in the rotary file sequence.
12) For manual shaping techniques repeat step 7 for each step in the crown-down technique sequence until canal shaping and cleaning is completed.
13) Once shaping is completed, thoroughly flush and evacuate canal using water to remove all irrigants. Dry canal and obturate.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- Do not sterilize the product.
- When using a dappen dish to dispense the gel ensure that the dish is properly disinfected prior to each use.
- A single-unit dose applicators and syringe tips of Glyde File Prep™ are designed to be used during one endodontic case (tooth) with multiple canals.
- Do not re-use single-unit dose applicator on multiple teeth or multiple patients.
- Do not re-use syringe tip on multiple patients. Re-use increases the risk of cross-contamination.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE 2797

Authorised representatives

EU REP

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP

DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH