Dentsply Sirona

ProFile®

EN	Instructions For Use	2
FR	Instructions d'utilisation	10
IT	Istruzioni per l'uso	19
DE	Gebrauchsanweisung	28
ES	Instrucciones de uso	37
PT	Instruções de utilização	46
NL	Gebruiksinstructies	55
SV	Användarinstruktioner	64
DA	Brugsanvisning	73
NO	Brukerinstruksjoner	82
LT	Naudojimo instrukcija	91
LV	Lietošanas norādījumi	100
ET	Kasutusjuhised	109
CS	Návod k použití	118
SK	Návod na použitie	127
HU	Használati utasítás	136
RO	Instrucţiuni de utilizare	145
SR	Uputstva za upotrebu	154
BS	Uputstva za upotrebu	163
CNR	Uputstva za upotrebu	172
SQ	Udhëzimet për përdorim	181
EL	Οδηγίες χρήσης	189
TR	Kullanım Kılavuzu	198
PT-BR	Instruções de uso	207
KO	사용 지침	216

Dentsply Sirona

EN

ProFile®

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
SINGLE USE – REPROCESSING PROHIBITED / REUSABLE – LIABLE FOR REPROCESSING
STERILE – STERILIZED BY RADIATION / NON-STERILE – STERILIZE BEFORE FIRST USE

INSTRUCTIONS FOR USE

SINGLE USE STERILE ENDODONTIC INSTRUMENTS
REUSABLE NON-STERILE ENDODONTIC INSTRUMENTS

1 PRODUCT DESCRIPTION

- Product description:
 ProFile® are engine-driven tapered endodontic instruments made of nickel titanium alloy with cutting edges. These instruments are used to mechanically shape and prepare the root canals (shaping and debridement of the root canal) during endodontic therapy. ProFile® endodontic files are available either in single use in sterile blister or in multiple use in non-sterile version.
- Clinical benefits:
 The clinical benefit for the patient being treated with ProFile® is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation (shaping and debridement of the root canal).
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

ProFile® instruments are delivered in both versions: sterile blisters or non-sterile plastic boxes.

The ProFile® range comprises the following endodontic instruments:

Code	Product	ISO XXX YYY XXX-ISO size YYY-taper size	Available lengths	Packaging and content
A 0344	ProFile® STERILE .02	ISO 015 002 ISO 020 002 ISO 025 002 ISO 030 002 ISO 035 002 ISO 040 002	21/25 mm	Blister with 6 instruments (sterile)
	ProFile® STERILE .02 Assortment	ISO 015 002 ISO 020 002 ISO 025 002 ISO 030 002 ISO 035 002 ISO 040 002	21/25 mm	Blister with 1 instrument of each ISO (sterile)
	ProFile® .02	ISO 015 002 ISO 020 002 ISO 025 002 ISO 030 002 ISO 035 002 ISO 040 002	21/25 mm	Plastic box with 6 instruments (non-sterile)
	ProFile® .02 Assortment	ISO 015 002 ISO 020 002 ISO 025 002 ISO 030 002 ISO 035 002 ISO 040 002	21/25 mm	Plastic box with 1 instrument of each ISO (non-sterile)
A 011N	ProFile® .04	ISO 015 004 ISO 020 004 ISO 025 004 ISO 030 004 ISO 035 004 ISO 040 004 ISO 045 004 ISO 060 004 ISO 090 004	21/25/31 mm 18 mm (for ISO 015-040)	Plastic box with 6 instruments (non-sterile)
	ProFile® .04 Assortment	ISO 015 004 ISO 020 004 ISO 025 004 ISO 030 004 ISO 035 004 ISO 040 004	21/25/31 mm 18 mm (for ISO 015-040)	Plastic box with 1 instrument of each ISO (non-sterile)
	ProFile® .04 Assortment	ISO 045 004 ISO 060 004 ISO 090 004	21/25/31 mm	Plastic box with 2 instruments of each ISO (non-sterile)

Code	Product	ISO XXX YYY XXX-ISO size YYY-taper size	Available lengths	Packaging and content
A 0345	ProFile® STERILE .06	ISO 015 006 ISO 020 006 ISO 025 006 ISO 030 006 ISO 035 006 ISO 040 006	18/21/25 mm	Blister with 6 instruments (sterile)
	ProFile® STERILE .06 Assortment	ISO 015 006 ISO 020 006 ISO 025 006 ISO 030 006 ISO 035 006 ISO 040 006	18/21/25 mm	Blister with 1 instrument of each ISO (sterile)
	ProFile® .06	ISO 015 006 ISO 020 006 ISO 025 006 ISO 030 006 ISO 035 006 ISO 040 006	18/21/25 mm	Plastic box with 6 instruments (non-sterile)
	ProFile® .06 Assortment	ISO 015 006 ISO 020 006 ISO 025 006 ISO 030 006 ISO 035 006 ISO 040 006	18/21/25 mm	Plastic box with 1 instrument of each ISO (non-sterile)
A 0359	ProFile® O.S.	O.S.1: ISO 020 005 O.S.2: ISO 030 006 O.S.3: ISO 040 006 O.S.4: ISO 050 007 O.S.5: ISO 060 008 O.S.6: ISO 080 008	19 mm	Plastic box with 6 instruments (non-sterile)
	ProFile® O.S.	O.S.1 O.S.2 O.S.3 O.S.4 O.S.5 O.S.6	19 mm	Plastic box with 1 instrument of each ISO (non-sterile)

5 COMPOSITION

ProFile® instruments are made of three main components: the working part made of a nickel-titanium alloy, a colored silicone stopper and a plated brass shank with colored rings.

6 COMPATIBLE DEVICES

- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).
- Torque control devices ensure optimal usage.

- Use with endodontic motors in a constant rotation at a speed of 150 - 350 rpm.
- Set torque at:

File size	Torque [Ncm]	Speed [rpm]
ProFile® ISO 015 002 ProFile® ISO 020 002 ProFile® ISO 025 002 ProFile® ISO 030 002 ProFile® ISO 035 002 ProFile® ISO 015 004 ProFile® ISO 020 004 ProFile® ISO 025 004 ProFile® ISO 030 004	0.6 - 1	150 - 350
ProFile® ISO 040 002 ProFile® ISO 035 004 ProFile® ISO 040 004	1 - 1.5	
ProFile® ISO 045 004 ProFile® ISO 030 006 ProFile® ISO 035 006 ProFile® ISO 040 006 ProFile® O.S.2 (ISO 030 006) ProFile® O.S.3 (ISO 040 006) ProFile® O.S.4 (ISO 050 007)	2 - 3	
ProFile® ISO 060 004 ProFile® ISO 090 004 ProFile® O.S.5 (ISO 060 008) ProFile® O.S.6 (ISO 080 008)	3 - 4	
ProFile® ISO 015 006 ProFile® ISO 020 006 ProFile® ISO 025 006 ProFile® O.S.1 (ISO 020 005)	1.5 - 2	

- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicone stopper to the correct working length is the appropriate method of working length determination.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use and the General Processing Instructions For Endodontic Products (see section [9]) to minimize the following risks to the device, the patient and/or the user:
 - Breakage of instrument.
 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Toxic or allergic reactions caused by processing residues.

- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Similar to all machine-driven root canal instruments, ProFile® instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible and/or NiTi files (this will help avoid canal transportation).
 - Visually inspect the working part for all the defects listed in the former paragraph during use (i.e. after each wave).
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

1) Before use, check if the instrument is sterile or if processing is needed (refer to Processing Instructions).
2) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
3) Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.
4) Irrigate abundantly and frequently after each instrument used (according to good dental practices).

ProFile® sequence for narrow canals:

Step	Instrument	ISO color	Number of rings on handle
Crown down Sequence	ProFile® O.S.3 (ISO 040 006)	Red	3
	ProFile® O.S.2 (ISO 030 006)	Yellow	3
	ProFile® ISO 025 006	Red	2
	ProFile® ISO 020 006	Yellow	2
	ProFile® ISO 025 004	Red	1
	K-File 015	White	N/A
	ProFile® ISO 020 004	Yellow	1
Apical preparation	ProFile® ISO 020 004	Yellow	1
	ProFile® ISO 025 004	Red	1
Flaring	ProFile® ISO 020 006	Yellow	2

5) Start using every instrument for approx. 5-10 seconds and never force a ProFile® into the canal.
6) When progression becomes difficult, withdraw the instrument and go on to the next smaller file size.
7) To determine the exact working length, use a K-File as a depth gauge after the first ProFile® has reached the minimum estimated working length less 3 mm (ex.: 21-3 = 18 mm).
8) The Crown-Down Sequence will be continued until the exact working length is reached.
9) Continue Crown-Down Sequence up until the foramen, using an apex locator as reference.
10) Continue with the apical preparation to the exact working length.
11) Use ProFile® with larger taper 004 if necessary according to the anatomy of canal.
12) Optional finish with the final flaring step according to the obturation method.
13) Use ProFile® with larger taper 006 if necessary according to the anatomy of canal. Penetration to the exact working length is not systematically sought, except if the root canal easily receives the ProFile® with taper 006.

14) Use dedicated Paper Points to dry and dedicated gutta-percha points (like Thermafil® (plastic core coated with gutta-percha)) to obturate the root canals.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- For instruments that are labelled “sterile”:
 - ProFile® instruments are single use devices and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
 - Disinfection, cleaning and sterilization are not applicable.
- For instruments that are labelled “non-sterile”:
 - Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.
 - Do not re-use silicone stops. Remove and discard silicone stops after each use.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

CE 2797

Authorised representatives

EU REP

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP

DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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Dentsply Sirona Europe GmbH