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## EN Vector® cleaner Special detergent for the hose system

Vector® cleaner is a special detergent for the removal of acid-soluble contaminants in the hose system and handpiece of the Dürr Dental Vector® system. Vector® cleaner is characterised by its intensive cleaning action and by the fact that it is very gentle to materials.

**Composition**: Vector® cleaner contains: inorganic acids (5 - 15 %), tensides (< 1 %) and additives in an aqueous solution.

**Directions for use**: Follow the procedure described in the Vector® operating instructions when cleaning the Vector® system. The system should be cleaned after 50 hours of operation and every 4 weeks. Vector® cleaner is used undiluted. Place approximately 500 ml of Vector® cleaner in the water tank of the Vector® and perform cleaning as directed in the operating instructions for the Vector® system. Leave product for an exposure time of 10 minutes (60 minutes maximum). Thereafter, rinse tubes and handpiece thoroughly with water as directed in the operating instructions for the Vector® system.

**Ch.-B./expiry date**: see imprint.

Irritant. Irritating to eyes and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable gloves.

DÜRR DENTAL (Products) UK Ltd.
14 Linnwell Way – Telford Way Industrial Estate
Kettering, Northants NN16 8PS, ☎ +44 1536 526740

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Distributor: DÜRR DENTAL AG, Ing. K. Bažil, Průběžná 74a,
Distribútor: DÜRR DENTAL AG, Ing. K. Bažil, Průběžná 74a,