FD 360

DE	Kunstlederreinigung und -pflege
EN	Cleaning and care of vinyl upholstery
FR	Nettoyage et entretien du simili-cuir
IT	Pulizia e cura dell'ecopelle
ES	Limpieza y mantenimiento del cuero
PT	Limpeza e cuidado da pele artificial
NL	Reiniging en verzorging van kunstleder
SV	Rengöring och vård av konstläder
FI	Tekonahan puhdistus ja hoito
NO	Rensing og behandling av kunstlær
DA	Kunstlæderrensning og -pleje
EL	Καθαριμός και περιποίηση συνθετικού δέρματος
PL	Czyszczenie i pielęgnacja sztucznej skóry
CS	Čištení a péče o umělou kůži
SK	Čistenie a ošetrovanie umelej kože
SL	Čistilno in negolvano sredstvo za umetno usnje
HR	Čišćenje i njega umjetne kože
HU	Műbőrtisztítás és -ápolás
EE	Kunstnaha puhastamine ja hooldus
LV	Mākslīgās ādas tīrīšana un kopšana
LT	Priemonė dirbtinei odai valyti ir prižiūrėti
RU	Очистка и уход за искусственной кожей
BG	Почистване и поддържане на изкуствена кожа
RO	Curățarea și întreținerea pielii sintetice
TR	Suni Deri Bakımı ve Temizliği

DÜRR DENTAL AG
Höpfigheimer Straße 17
74321 Bietigheim-Bissingen, Germany
Tel.: +49 7142 705-0
Fax: +49 7142 705-500
info@duerr.de · www.duerr.de

Z.360-2013-04

DÜRR DENTAL

DE

FD 360

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EN

FD 360

Cleaning and care of vinyl upholstery

FD 360 cleans and cares for all seat covers made from artificial leather and is highly effective in maintaining the value of dental treatment chairs and other artificial leather covers. You can apply FD 360 with the special sponge supplied in the set to rapidly and thoroughly remove stubborn grime, stains, and discoloration, caused for example by non-fast dyestuffs in clothing. FD 360 cleans and cares intensively without leaving any residues; it provides a silky gloss finish.

Composition: FD 360 is based on a mixture of special tensides, silicone compounds and care components derived from avocado oil.

Application: Use FD 360 undiluted. Spray FD 360 onto the surfaces to be cleaned and then immediately wipe off with a cloth. For intensive cleaning rub FD 360 on with the special sponge, which achieves an even more thorough and deeper cleaning effect. Subsequently, wipe off the excess FD 360 with a dry cloth. Use FD 360 daily if required.

Notes: FD 360 can be used without any problems on all artificial leather surfaces. The material tolerability of other leather surfaces should be checked on a concealed surface. The special sponges should **not** be used on leather.

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Türkiye'de satış: DÜRR DENTAL AG