REPROCESSING INSTRUCTIONS

GUIDED BIOFILM THERAPY®

"I FEEL GOOD"

EMS

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BEFORE USE

Please read these Operating Instructions carefully as they explain all the most important details and procedures. Please pay special attention to the safety precautions.

Always keep these Instructions close at hand.

⚠ These Instructions are only applicable to the equipment they were delivered with.

To prevent injury to persons or damage to property, please observe the corresponding directives and symbols.

This instruction follows the recommendations of RKI Robert Koch Institut, Germany.

Translations in this manual are based on the English version. Thus the English text will be referred to in case of discrepancy.

CLEANING AND DISINFECTION FOR DEVICES AND HANDPIECE CORDS


⚠ The following instructions have been validated for all reusable products[1] of EMS dental devices. It remains the responsibility of the processor to ensure that the processing, as actually performed using the equipment, materials and personnel in the processing facility, achieves the desired result. This requires validation and/or validation and routine monitoring of the process. Likewise any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse effects.

AIRFLOW® and PIEZON® devices

Handy

PIEZON® Kit on your unit

CLEAN AND DISINFECT

Use a cleaning and disinfecting wipe (less than 35% alcohol) that complies with the standards applicable in the country.

EMS recommends the use of either

- CleanWipes, IC-100 or PlastiSept Eco Wipes from ALPRO MEDICAL GmbH

- Micro-Kleen3™ Surface Disinfectant Wipes from Micro-Scientific

⚠ Follow carefully the instructions provided by the wipes manufacturer.

⊘ DO NOT use CaviWipes™ or Advantaclear wipes. They damage EMS products.

INTRODUCTION

EMS recommends compliance with cleaning, disinfection, packaging for sterilization and sterilization procedures according to ISO 17664.

❗ Always report adverse events related to device reprocessing directly to EMS.

⚠ Products must be cleaned, disinfected and, if applicable, sterilized prior to first use. Do not reprocess the products over the authorized number of sterilization cycles. To replace, refer to the “Service life” section of these Instructions.

⚠ Concentrations and contact times specified by the manufacturer of the cleaning and disinfection agent must be observed.

⚠ Remember that sterilization mustn't be achieved unless assembly components are cleaned and disinfected first.

If parts of these instructions are not clear or seem to be inadequate, do not hesitate to contact/inform EMS.

❗ **Users shall also observe any legal requirements applicable in their country as well as the hygiene regulations in effect at the hospital or clinic.**

This applies especially with regard to additional requirements for the inactivation of prions.

SERVICE LIFE

The products have been designed for a large number of sterilization cycles. The materials used in their manufacture were selected accordingly. However, with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products.

❗ Always replace products that present signs of wear or early degradation, regardless of the number of sterilization cycles left unused.

⊘ DO NOT expose the products to temperature exceeding 138°C.

⚠ EMS handpieces must not be lubricated or otherwise exposed to oil whether directly or indirectly.

Lifetime : see product's Instructions for Use.

SYMBOLS

⚠ Caution / Danger - Risk of property damage or injury

❗ Please note - Useful additional information and hints

⊘ Prohibited

Wear personal protective equipment

Refer to instruction manual


135°C Sterilizable at up to 135°C in the autoclave

Suitable for thermal disinfection

Required time - Abbreviations used: min (minute) ; s (second)

1 Hereinafter referred to as "products".

PROCESS FOR RELATED PRODUCTS MENTIONED ON PAGE 3


01 PREPARE

At the point of use: Immediately after patient treatment, rinse the channel(s) of the handpiece/instrument with cold tap water. Coarse soiling must be removed by wiping the products.


Special attention need to be paid to cracks, crevices and hard to reach areas.

Place the products on an appropriate rack to transport them to the reprocessing area to avoid any product damage, contamination to the environment and people involved in the reprocessing process.

Always carry out **AIRFLOW® and PERIOFLOW®** Handpiece powder unclogging and check for both channels (water and powder) clearance before proceeding. Use Easy Clean[2].

Products have to be disassembled as shown below:

PIEZON® Endochuck:

Remove any installed file.
Remove the screw and separate the small O-ring.

PIEZON® Handpieces:

Remove any installed instrument.
Remove the nose cap and separate the light guide and O-ring.

PIEZON® Instruments with CombiTorque®:

Separate the PIEZON® Instrument and the CombiTorque® but keep each original instrument and CombiTorque® together throughout their life.

Cleaning is to be performed within 1 hour from use.
Wear personal protective equipment according to the type of preparation.

Any product must be cleaned manually or automatically by washer or disinfector. Follow the dedicated cleaning part according to your dental equipments.

EMS recommends the use of an ISO 15883 compliant automatic washer-disinfector (WD) for optimal effectiveness and service life.

EMS handpieces must not be lubricated or otherwise exposed to oil whether directly or indirectly.

2 An air water syringe can be used with Easy Clean

EITHER 02 03 04 AUTOMATIC CLEAN, DISINFECT AND DRY WITH WASHER DISINFECTOR

The washer-disinfector must comply with ISO 15883, must have suitable baskets to hold small products and must have rinsing connections with diameter of approximately 16mm for the attachment to product lumen.

Correctly place the product into a suitable rack, connect all lumens to the rinsing connectors and start the automated cleaning.

Instructions of the manufacturer of the washer disinfector must also be observed.

Carefully follow the instructions provided by the disinfection solution manufacturer.

The following automated process[3] can be used to achieve an A0 Level of 600:

2 min	Pre-clean with cold tap water[4]. Drain
5 min	Clean at 55°C with tap water and 0.5% neodisher MediClean Dental (Dr Weigert) solution. Drain
3 min	Rinse and neutralize with cold demineralized water[5]. Drain
2 min	Rinse with cold demineralized water[5]. Drain
1 min	Perform thermal disinfection (final rinse) with deionized water at 90°C minimum. Drain
20 min (at least)	Dry at 100°C.

OR 02 MANUAL CLEAN + 03 DISINFECTION + 04 DRY

The following validated EMS process can be used with and without ultrasound.

DO NOT use any ultrasound bath cleaning procedure with the handpieces: it may destroy the products. It is ONLY intended to be used with EMS PIEZON® Instruments and file holders. In case of ultrasound bath cleaning, place products in a metal sieve and sonicate 10 minutes.

<table>
<tr><th>Validated product</th><th>Process</th><th>EMS products to clean/ disinfect</th><th>Process</th></tr>
<tr><td>0.5% neodisher MediClean (Dr. Weigert)</td><td>without ultrasound bath</td><td>all EMS products</td><td rowspan="2">Brush the devices in a cleaning solution of 0.5% neodisher in deionized water at 40°C with a suitable soft bristled brush until all visible residue has been removed.
15 s — For products with lumens, flush all lumens with a spray gun (water jet gun, with a static water pressure of 2 bar) with cold tap water.
15 min (10 min with ultrasound) — Place the product in 0.5% solution of neodisher with deionized water at 40°C. Make certain that all lumens are filled with cleaning solution (use a syringe if needed).
15 s — Rinse all lumens by flushing with a spray gun (water jet gun, with a static water pressure of 2 bar) with cold deionized water[5].
10 s — Rinse the whole product under cold running deionized water.
At room temperature (15-25°C), use an air pistol (compressed air, max. 3 bars) to fully dry the lumen and the whole product until no more residues of water are present (visible or detectable).</td></tr>
<tr><td>0.5% neodisher MediClean (Dr. Weigert)</td><td>with ultrasound bath[6]</td><td>EMS PIEZON® Instruments and file holders</td></tr>
</table>

Then disinfect:

ASP CIDEX® OPA solution must be used undiluted, within its useful and shelf life, and in compliance with the manufacturer's warnings and instructions for use.

5 min	Immerse the product completely into CIDEX® OPA solution at minimum 20°C. Make certain that all lumens are filled with disinfectant solution (use a syringe if needed).
1 min (each rinse)	The ASP CIDEX® OPA disinfectant requires a total of three rinses : keep the product totally immersed, and use a large volume[7] of fresh water.

Do not reuse the water for rinsing or any other purpose. Residues of disinfectant may cause serious side effects.

And dry

At room temperature (15-25°C), use an air pistol (compressed air, max. 3 bars) to fully dry the lumen and the whole product until no more residues of water are present (visible or detectable).

05 INSPECT

If stains are still visible on the product after cleaning/disinfection, the entire cleaning/disinfection procedure must be repeated. Products with visible damage, chip/flake loss, corrosion or deformation must be disposed of (no further use is allowed). Check also the integrity of O-rings and gaskets and replace if damaged or deformed.

Verify the product to be fully dry. In case of detection of residues of water, remove these using an air pistol (clean compressed air). Fully dry the lumen and the whole product until no more residues of water are present (visible or detectable).


3 For example: Miele Professional G 7836 CD with Miele Rack E429
4 cold tap water = 16°C +/- 2°C
5 cold demineralized water = 20°C +/- 2°C
6 For example: Bandelin, Sonorex 1028K, 35kHz
7 e.g. 2 gallons

06 PACK

⚠ Only previously cleaned and disinfected products must be sterilized.

⚠ Effective sterilization must take place only on fully dry products. Ensure each part (internal lumens and any surface) to be perfectly dry before reassembling and packing.

Prior to sterilization, the products need to be placed in suitable sterilization packaging.

PIEZON® Handpieces:

This product must be sterilized either disassembled or reassembled.

If you prefer to reassembled it: reinstall the O-ring first, then put the light guide into the nose cap and screw it on the handpiece.

PIEZON® Endochuck:

Reinstall the small O-ring (gasket) first, and then gently screw the retention bolt without tightening.

PIEZON® Instruments with CombiTorque®:

Reinstall the original instrument with its CombiTorque®.

Package your products with a single or double pouch:
- suitable for prevacuum moist heat sterilization,
- compliant with ISO 11607-1 or EN 868,
- resistant to 138°C
- and with adequate steam permeability (e.g. Wipak STERIKING flat rolls Type R43 and R44).

07 STERILIZE

⚠ Sterilization must be performed immediately after cleaning-disinfection.

The loading pattern defined by the autoclave manufacturer must be follow.

DO NOT exceed a sterilization temperature of 138°C and a holding time of 20 min.

DO NOT use hot-air sterilization and radio-sterilization procedures: they destroy the products.

EMS recommends not to use class N sterilizers to reprocess PIEZON® handpieces unless they include at least 20 minutes of drying cycle. The PIEZON® Handpiece lifetime can be reduced.

Moist heat sterilization of products shall be performed according to ISO 17665 and in compliance with respective country requirements.

The prevacuum moist heat (steam) process can be used with any EMS product packaged in appropriate single or double pouches:

Parameters for the prevacuum moist heat cycle:
- 3 prevacuum phases
- Pressure of 3 bar[8]
- Humidity of 100%
- Minimum temperature of 132°C
- 3 min (minimum) Hold time (full cycle)
- 20 min (minimum) Dry

⚠ Users must ensure that the reprocessing processes, including resources, materials and personnel, are capable to reach the required results and maintained over time: it is the user's responsibility to ensure that validation of the reprocessing procedures is kept up to date at all times.

08 STORE

Store the sterilized products at a temperature up to 40°C in a:
- dry,
- clean
- and dust-free environment.

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VI

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FROM THE VALLEY

▸ For more than 35 years EMS has developed high-end technologies and protocols in conjunction with the world's most advanced clinicians, enabling the patients to enjoy natural teeth and implants for much longer.

„Primum non nocere" and minimal invasiveness are not only a must, they are our way of life.

EMS RESEARCH AND MANUFACTURE CENTER
NYON, SWITZERLAND

TO THE TOP

ONLY USE ORIGINAL EMS PRODUCTS

▶ EMS powders are much less abrasive than traditional mechanical debridement techniques such as rubber cups and handscaling.
We recommend to use AIRFLOW® and PERIOFLOW® first, followed by PIEZON®. GBT: the minimally invasive way for Professional Mechanical Plaque Removal (PMPR).

▶ As part of the GBT disclosing phase, coloration also minimizes AIRFLOW® abrasion. No more color signifies that all biofilm has been removed: the objective of AIRFLOW® is achieved, to preserve the integrity of dental tissues.

⚠ EMS devices have been designed and tested for exclusive use with EMS prophylaxis powders.

ONLY USE EMS POWDERS

So-called "EMS compatible" powders in the marketplace may cause damage to oral tissues and the device. See below a damaged Airflow nozzle after use with a "compatible" powder. If EMS products are damaged by using "compatible" powders the product Warranty will be lost.

Damaged handpiece

ONLY USE THE SWISS ORIGINAL EMS INSTRUMENTS!

The EMS Instrument, the Handpiece and the NO PAIN electronic module have to vibrate in harmony like in a trilogy.
Using so-called “compatible” and copy tips from third parties can damage your patients' teeth as well as EMS handpieces.
They were developed and made for each other by EMS.
One PS Instrument is designed for up to 1000 treatments, thus reducing the cost of one treatment to only a few cents.
Therefore, there is absolutely no reason to "save" money by using so-called "compatible" tips. Moreover, if the Piezon Handpiece gets damaged (the thread) it will lose its warranty so at the end of the day, “compatible” tips will always cost you more.

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© FB-358 rev.I Edition 2024/01
E.M.S. ELECTRO MEDICAL SYSTEMS TRADING (SHANGHAI) CO., LTD.
EMS ELECTRO MEDICAL SYSTEMS KOREA LTD.