# AIRFLOW PERIO

Before using this product, please carefully read and follow the recommendations of the instruction manuals originally delivered with your device and with your Biofilm supplier. Please pay special attention to the safety precautions.
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**Intended use***

This powder is intended to remove biofilm and stains from oral hard tissues and materials.

**Indications***

- Remove supragingival biofilm from oral hard tissues and materials
- Remove subgingival biofilm from oral hard tissues and materials
  - AIRFLOW®: up to 4mm
  - PERIOFLOW®: deeper than 4mm
- Remove stains

**Intended user population**

This powder is intended to be used by qualified dental professionals.

**Intended environment**

The AIRFLOW® Powders are intended to be used in a dental cabinet/setting complying with each country's regulations.

**Intended patient population**

Patients above the age of 2 years old.

**Contraindications***

This powder is contraindicated on patients allergic to Glycine.

**Side effect***

None

**Residual risks***

Allergic reaction.

**Recommendations**

Always disclose biofilm before removing it with AIRFLOW®.
Do not overtreat. Stop the treatment when/where the disclosing agent is gone.
The technique is not recommended for periapical tissues. Prerupturre nevertheless in the clinical stages of periodontal disease, some sites are contraindicated as described in the AIRFLOW® device instructions for use.
Please refer to the instructions for use for further information.

**Storage and shelf-life**

Storage conditions at 25°C and up to 90% relative humidity.
Transport conditions between -29°C to 38°C and up to 90% relative humidity.
Shelf life: 26 months.

**Disposal**

There is no specific safety risk related to the disposal of the product and its container. Disposal must be made according to country specific regulations.

**Main ingredient:**

Glycine

**Contains:**

Hydrophobic modified amorphous silica

Translations in this manual are based on the English version. Thus the English text will be referred to in case of discrepancy.

*AIRFLOW® devices have their own characteristics (intended use, indications, contraindications, side effects, potential complications, recommendations, ...) please refer to their instruction for use for further information.

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© Copyright EMS SA
**E.M.S. ELECTRO MEDICAL SYSTEMS TRADING (SHANGHAI) CO., LTD.**
**EMS ELECTRO MEDICAL SYSTEMS KOREA LTD.**