AIRFLOW COMFORT

Before using this product, please carefully read and follow the recommendations of the instruction manual originally delivered with your device and system. Please pay special attention to the safety precautions.
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Intended use
This powder is intended to remove supragingival biofilm and stains from oral hard tissues.

Indications
- Remove supragingival biofilm from oral hard tissues
- Remove stains

Intended user population
This powder is intended to be used by qualified dental professionals.

Intended environment
The AIRFLOW® Powders are intended to be used in a dental cabinet/setting complying with each country's regulations.

Intended patient population
Patients over the age of 12 y/o.

Contraindication for use*
This powder is contraindicated on patients allergic to flavor aromas.

Side effect*
None

Residual risks*
Allergic reaction

Recommendations
This powder contain salt and should not be used on patients with a low salt diet or hypertension.
- Do not use the device in the direction of oral soft tissues or any dental tissue other than sound enamel.
- Always close bottle before removing it with AIRFLOW®.
- Do not overfill bottle. Stop the treatment when the dosing agent is gone.

Special considerations where clinicians must use their own discretion or consent from the physician:
- Although there have been no clinical complaints or data showing issue when using COMFORT Powder as recommended, EMS recommends the AIRFLOW PLUS Powder instead to treat patients with any form of low salt diet such as patients with hypertension or patient with renal insufficiency. For patients with acute or uncontrolled hypertension, EMS recommends to get physician's consent before using COMFORT Powder.
- Although there is no specific patient contraindication related to the clinical stages of periodontal disease, the AIRFLOW® is not recommended to be used on subgingival calculus as described in the AIRFLOW® device instructions for use. Please refer to the device instructions for use for more information.

Storage and transport conditions
Storage condition up to 30°C and up to 90% relative humidity.
Transport conditions between -29°C to 38°C and up to 90% relative humidity.

Shelf life*
3 years

Disposal
There is no specific safety risk related to the disposal of the product and its container. Disposal must be made according to country specific regulations.

Main ingredient
Sodium Bicarbonate

Contains:
Sweetener (sucralose for Mint flavor / saccharin for Lemon and Cherry flavors)
Aroma (not present in Neutral flavor)
Hydrophobic modified amorphous silica

Translations in this manual are based on the English version. Thus the English text will be referred to in case of discrepancy.

*AIRFLOW® devices have their own characteristics (intended use, indications, contraindications, side effects, potential complications, recommendations, ...) please refer to their instruction for use for further information.

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