GBT FLOWCONTROL® – HIGH-VACUUM SUCTION CANNULA

LANGUAGES

- (EN) ENGLISH
- (BG) БЪЛГАРСКИ
- (CS) ČESKÝ
- (DA) DANSK
- (DE) DEUTSCH
- (EL) ΕΛΛΗΝΙΚΑ
- (ES) ESPAÑOL
- (ET) EESTI KEEL
- (FI) SUOMI
- (FR) FRANÇAIS
- (HR) HRVATSKI
- (HU) MAGYAR
- (ID) INDONESIA
- (IT) ITALIANO
- (JP) 日本
- (KO) 한국어
- (LT) LIETUVIŠKAI
- (LV) LATVIEŠU
- (MS) BAHASA MELAYU
- (NL) NEDERLANDS
- (NO) NORSK
- (PL) POLSKI
- (PT) PORTUGUÊS
- (RO) ROMÂNĂ
- (RU) РУССКИЙ
- (SK) SLOVENSKY
- (SL) SLOVENŠČINA
- (SV) SVENSKA
- (TH) ไทย
- (TR) TÜRKÇE
- (VI) TIẾNG VIỆT
- (ZH) 中文

HEADQUARTERS AND MANUFACTURER

EMS ELECTRO MEDICAL SYSTEMS

(CH) **NYON, SWITZERLAND**
E.M.S. ELECTRO MEDICAL SYSTEMS S.A.
Chemin de la Vuarpillière 31 - CH - 1260 Nyon
www.ems-dental.com/en/contact
T. +41 22 994 47 00 – F. +41 22 994 47 01

WORLDWIDE EMS AFFILIATES

(DE) **MUNICH, GERMANY**
EMS ELECTRO MEDICAL SYSTEMS GMBH
Stahlgruberring 12 – DE - 81829 München
T. +49 89 42 71 61 0 – F. +49 89 42 71 61 60
Email: info@ems-ch.de

(FR) **FONCINE, FRANCE**
EMS ELECTRO MEDICAL SYSTEMS FRANCE SARL
32, route de Pontarlier – FR - 39460 Foncine-le-Haut
T. +33 3 84 51 90 01 – F. +33 3 84 51 94 00
Email: info@ems-france.fr

(ES) **MADRID, SPAIN**
EMS ELECTRO MEDICAL SYSTEMS ESPAÑA SLU
C/ Tomás Bretón, 50-52 2ª planta – ES - 28045 Madrid
T. +34 91 528 99 89 – Email: info@ems-espana.com

(PT) **LISBON, PORTUGAL**
EMSPOR - ELECTRO MEDICAL SYSTEMS PORTUGAL, UNIPESSOAL LDA
Avenida Miguel Bombarda, 36, 6A – PT - 1050-165 Lisboa
T. +351 912 338 777 – E-mail: info@ems-portugal.com

(IT) **MILANO, ITALY – EMS ITALIA S.R.L**
Via Carlo Pisacane 7B – IT - 20016 Pero (MI)
T. +39 02 3453 8111 – +39 02 3453 2778
Email: dental@ems-italia.it

(UK) **CAMBRIDGESHIRE, ENGLAND – EMS UK LTD SUITE 2**
Xenus House, Sandpiper Court Phoenix Park - Eaton Socon
UK - St Neots, Cambridgeshire – PE19 8EP
T. +44 01480 587260 – Email: info@ems-unitedkingdom.com

(SE) **STOCKHOLM, SWEDEN – EMS NORDIC AB**
Hammarby Fabriksväg 61 – SE - 120 30 Stockholm
T. +46 8 899 102 – Email: info@ems-nordic.se

(PL) **TYCHY, POLAND – EMS POLAND SP. Z O.O.**
Fabryczna 43 - PL - 43-100 Tychy
T. +48 32 4937060 – Email: info@ems-poland.com

(GR) **THESSALONIKI, GREECE EMS EAST MEDITERRANEAN**
10 km. Thessaloniki - Thermi - GR – Thessaloniki 57001
T. +30 231 231 61 30 – Email: info@ems-poland.com

(JO) **AMMAN, JORDAN**
EMS ELECTRO MEDICAL SYSTEMS JORDAN
No.2 - second floor, 195 Arar (Wadi Saqra) str. - JO - Amman

(US) **DALLAS, USA – ELECTRO MEDICAL SYSTEMS L.L.C.**
2150E Continental Boulevard, Southlake, Texas, 76092, US
T. +001 800-367-0367 – E-mail: info@ems-na.com
Aftersales: emsrepairs@ems-na.com

(AU) **SYDNEY, AUSTRALIA – EMS OCEANIA PTY LTD**
The Woolstores - Suite 3, Shed 73 – 4E Huntley Street
AU - Alexandria, NSW 2015
T. +61 405 095 867 – E-mail: info@ems-oceania.com

(JP) **TOKYO, JAPAN – E.M.S. JAPAN BRANCH OFFICE**
2F Takebashi Yasuda Building, 3-13 Kanda Nishiki-cho
JP - Chiyoda-ku, Tokyo, 101-0054
T. +81 (0) 3-6810-9351 – F. +81 (0) 3-6810-9352
Email : emsjapan@ems-ch.com

(CN) **SHANGHAÏ, CHINA**
医迈斯电子医疗系统贸易 (上海) 有限公司 E.M.S. ELECTRO MEDICAL SYSTEMS TRADING (SHANGHAI) CO., LTD.
24A, Jin Sui Mansion, N°379 Pudong Nan Rd. - CN - Shanghai 200120
T. +86 21 3363 2323 – E-mail: emschina@ems-ch.com

(KR) **SEOUL, SOUTH KOREA**
EMS ELECTRO MEDICAL SYSTEMS KOREA LTD.
13, Beobwon-ro 8-gil, Songpa-gu
T. +82 70 7807 2875 – E-mail: emskorea@ems-ch.com


© Copyright EMS SA FB-712 rev.D ed.2024/12

	EN	BG	CS	DA	DE	EL	ES	ET	FI	FR
EMS®	Manufacturer's logo	Лого на производителя	Logo výrobce	Producentens logo	Hersteller-logo	Λογότυπο κατασκευαστή	Logotipo del fabricante	Tootja logo	Valmistajan logo	Logo du fabricant
	Manufacturer	Производител	Výrobce	Producent	Hersteller	Κατασκευαστής	Fabricante	Tootja	Valmistaja	Fabricant
	Fabrication date	Дата на производство	Datum výroby	Fremstillingsdato	Herstellungs-datum	Ημερομηνία παρασκευής	Fecha de fabri-cación	Valmistamis-kuupäev	Valmistuspäivä	Date de fabrication
REF	Catalog reference	Каталожен номер	Katalogové číslo	Katalog reference	Katalog-verweis	Αναφορά καταλόγου	Referencia del catálogo	Kataloogi viitenumber	Luettelo-viite	Référence catalogue
LOT	Lot number	Партиден номер	Číslo šarže	Partinummer	Chargennummer	Αριθμός παρτίδας	Número de lote	Partii number	Eränumero	Numéro de lot
	General Warning	Общо преду-преждение	Obecné varování	Generel advarsel	Allgemeine Warnung	Γενική προειδοποίηση	Advertencia general	Üldine hoiatus	Yleinen varoitus	Avertissement général
MD	Medical Device	Медицинско изделие	Zdravotnický prostředek	Medicinsk udstyr	Medizinisches Gerät	Ιατροτεχνολο-γικό προϊόν	Dispositivo médico	Meditsiiniseade	Lääkintälaite	Dispositif médical
CE	CE marking (Medical Device compliant with EU Regulation 2017/745)	Символът CE се отнася до ЕС регламент 2017/745	Symbol CE ozna-čuje nařízení EU 2017/745	CE symbolet henviser til EU Forordning 2017/745	CE-Symbol bezieht sich auf EU-Verordnung 2017/745	Το σύμβολο CE αναφέρεται στον κανονι-σμό 2017/745 της ΕΕ	El símbolo CE se refiere al Regla-mento de la UE 2017/745	CE-sümbol viitab ELi määrusele 2017/745	CE-symboli viit-taa EU-määräykseen 2017/745	Symbole CE en référence au Règlement (UE) 2017/745
	Storage conditions	Условия за съхранение	Podmínky ucho-vávání	Opbevarings-forhold	Lagerbedin-gungen	Συνθήκες αποθήκευσης	Condiciones de almacena-miento	Hoiustamistingi-mused	Säilytysolo-suhteet	Conditions de stockage
135°C SSS	Sterilizable at up to 135°C in the autoclave	Подлежи на стерилизация при до 135°C в автоклав	Sterilizovatelné v autoklávu za teplot do 135 °C	Kan steril-iseres ved op til 135 °C i autoklave	Bei bis zu 135 °C im Autoklav steri-lisierbar	Αποστειρώσιμο σε θερμοκρασία έως 135 °C σε αυτόκαυστο	Esterilizable en el autoclave a una tempera-tura máxima de 135 °C	Autoklaavis ster-iliseeritav kuni 135 °C juures	Sterilointi enintään 135 °C:n auto-klaavissa	Stérilisable jusqu'à 135 °C en autoclave
	Thermal disinfection	Топлинна дезинфекция	Tepelná dezinfekce	Termisk desinfektion	Thermische Desinfektion	Θερμική απολύμανση	Desinfección térmica	Termiline desin-fektsioon	Terminen desin-fiointi	Thermo-désinfection
	Refer to instruction manual	Направете справка с ръко-водството с инструкции	Přečtěte si návod k použití	Henviser til brugsanvis-ningen	Siehe Bedie-nungsanleitung	Ανατρέξτε στο εγχειρίδιο οδηγιών	Consulte el manual de instrucciones	Vaadake kasu-tusjuhendit	Katso käyttöop-paasta	Se reporter au mode d'emploi
	EC Authorized Representative	Упълномощен представител за ЕО	Autorizovaný zástupce ES	EU-godkendt repræsentant	EG-Bevollmäch-tigter für Medizinpro-dukte	Εξουσιοδοτημέ-νος αντιπρόσω-πος στην ΕΚ	Representante autorizado en la CE	ELi volitatud esindaja	Valtuutettu edustaja EY-al-ueella	Représentant autorisé CE

	HR	HU	ID	IT	JP	KO	LT	LV	MS	NL	NO
EMS	Logo proizvođača	Gyártói embléma	Logo produsen	Logo del produttore	製造元ロゴ	제조업체의 로고	Gamintojo logotipas	Ražotāja logotips	Logo Pengilang	Fabrieks-logo	Produsentens logo
	Proizvođač	Gyártó	Produsen	Produttore	製造元	제조업체	Gamintojas	Ražotājs	Pengilang	Fabrikant	Produsent
	Datum proizvodnje	Gyártás dátuma	Tanggal pembuatan	Data di produzione	製造日	제조 날짜	Pagaminimo data	Ražošanas datums	Tarikh pembikinan	Fabricage-datum	Produksjons-dato
REF	Kataloška oznaka	Katalógus szerinti hivatkozás	Referensi katalog	Riferimento catalogo	カタログ参照	카탈로그 레퍼런스	Katalogo nuoroda	Kataloga atsauce	Rujukan katalog	Catalogus-referentie	Katalog-referanse
LOT	Broj serije	Tételszám (lot number)	Nomor lot	Numero di lotto	ロット番号	제품 번호	Partijos numeris	Partijas numurs	Nombor lot	Lotnummer	Partinummer
	Opće upozorenje	Általános figyelmeztetés	Peringatan Umum	Avviso generale	一般警告	일반적인 경고 사항	Bendro pobūdžio įspėjimas	Vispārīgs brīdinājums	Amaran Umum	Algemene waarschuwing	Generell advarsel
MD	Medicinski uređaj	Orvostechnikai eszköz	Perangkat Medis	Dispositivo medico	医療機器	의료 기기	Medicinos priemonė	Medicīniskā ierīce	Alat Perubatan	Medisch appa-raat	Medisinsk enhet
CE	Oznaka CE odnosi se na Uredbu 2017/745 EU-a	A CE jel értelmezéséhez lásd: EU-rendelet (2017/745)	Simbol CE merujuk pada Regulasi UE 2017/745	Il simbolo CE si riferisce al Regolamento UE 2017/745	CE マーキング (EU 規則 2017/745 準拠の医療機器)	CE 마킹(EU 규정 2017/745을 준수하는 의료 기기)	CE simbolis reiškia nuorodą į ES reglamentą 2017/745	CE simbols norāda uz ES Regulu 2017/745	Simbol CE merujuk kepada Peraturan EU 2017/745	CE-marke-ring verwijst naar EU-verordening 2017/745	CE-symbolet viser til EU-forordning 2017/745
	Uvjeti skladištenja	Tárolási körülmények	Kondisi penyim-panan	Condizioni di conserva-zione	保管条件	보관 환경	Laikymo sąlygos	Uzglabāšanas apstākļi	Keadaan penyimpanan	Opslag-condities	Lagringsbe-tingelser
135°C	Moguće sterilizirati pri temperaturi do 135 °C u autoklavu	Autoklávban, legfeljebb 135 °C-on sterilizálható	Dapat disteril-kan hingga 135°C dalam autoklaf	Sterilizzabile in autoclave fino a 135 °C	オートクレーブで最高 135°C で消毒	오토클레이브에서 최대 135°C에서 살균 가능	Galima ste-rilizuoti auto-klave iki 135 °C temperatūroje	Sterilizējams autoklāvā līdz 135 °C temper-atūrā	Boleh disteril pada hingga 135°C di dalam autoklaf	Steriliseerbaar tot 135°C in de autoclaaf	Steriliserbar opptil 135 °C i autoklav
	Termička dezin-fekcija	Termikus fertő-tlenítés	Disinfeksi termal	Disinfezione termica	熱消毒	열 소독	Terminé dezinfekcija	Termiskā dezin-fekcija	Pensterilan terma	Thermische ontsmetting	Termisk desinfeksjon
	Pogledajte korisnički priručnik	Tájékozódjon a használati útmutatóból!	Lihat buku petunjuk	Consultare il manuale di istruzioni	取扱説明書を参照	지침 설명서를 참고하십시오.	Žr. naudojimo instrukciją	Skatiet lietoša-nas instrukciju	Rujuk manual arahan	Lees de gebruiks-aanwijzing	Se bruksanvis-ningen
	Ovlašteni predstavnik u Europskoj zajednici	Meghatalma-zott képviselő az EK-ban	Perwakilan Resmi EC	Rappresen-tante Autoriz-zato EC	EC 認定代理店	EC Authorized Representative	EB įgaliotasis atstovas	Pilnvarotais pārstāvis EK	Wakil Sah Suruhanjaya Eropah (EC)	Gemach-tigd vertegen-woordiger van de EG	Autorisert representant i EF


	PL	PT	RO	RU	SK	SL	SV	TH	TR	VI	ZH
EMS	Logo producenta	Logótipo do fabricante	Logo producător	Логотип производителя	Logo výrobcu	Logotip proizvajalca	Tillverkarens logotyp	โลโก้ของผู้ผลิต	Üretici logosu	Logo của Nhà sản xuất	制造商标识
	Producent	Fabricante	Producător	Производи-тель	Výrobca	Proizvajalec	Tillverkare	ผู้ผลิต	Üretici	Nhà sản xuất	制造商
	Data produkcji	Data de fabrico	Data fabricării	Дата изготовления	Dátum výroby	Datum proizvodnje	Tillverknings-datum	วันที่ ประกอบขึ้น	Üretim tarihi	Ngày sản xuất	生产日期
REF	Numer katalogowy	Referência de catálogo	Referință catalog	Номер по каталогу	Katalógová referencia	Katalog referenc	Produkt-katalog referens	แคตตาล็อก การอ้างอิง	Katalog referans	Danh mục mã tham chiếu	产品编号
LOT	Numer partii	Número de lote	Nr. de lot	Номер лота	Číslo šarže	Številka serije	Partinummer	หมายเลขล็อต	Lot numarası	Số lô	批号
	Ostrzeżenie ogólne	Aviso geral	Avertisment general	Предупрежде-ние общего характера	Všeobecné varovanie	Splošna opozorila	Allmän varning	คำเตือนทั่วไป	Genel Uyarı	Cảnh báo chung	一般警告
MD	Urządzenie medyczne	Dispositivo médico	Dispozitiv medical	Медицинское изделие	Zdravotnícka pomôcka	Medicinska naprava	Medicinsk utrustning	อุปกรณ์ทางการแพทย์	Tıbbi cihaz	Thiết bị Y tế	医疗器械
CE	Symbol CE odnosi się do UE Rozporządze-nie 2017/745	O símbolo CE é relativo ao Regulamento UE 2017/745	Simbolul CE face trimitere la Regulamentul UE 2017/745	Символ CE обозначает соответствие директиве ЕС 2017/745	Symbol CE označuje nariadenie EÚ 2017/745	Simbol CE se nanaša na EU Uredbo 2017/745	CE-symbol hänvisar till EU Förordning 2017/745	สัญลักษณ์ CE อ้างถึง กฎระเบียบ 2017/745 ของ EU	CE sembolü AB'ne işaret eder 2017/745 Düzenlemesi	Biểu tượng CE đề cập đến EU Quy định 2017/745	CE标识 (符合欧盟法规2017/745的医疗器械)
	Warunki prze-chowywania	Condições de armazena-mento	Condiții pentru depozitare	Условия хранения	Podmienky skladovania	Pogoji shranjavanja	Förvaringsmiljö	เงื่อนไข การจัดเก็บ	Saklama koşulları	Điều kiện bảo quản	储存条件
135°C	Nadaje się do sterylizacji w autoklawie do 135°C	Esterilizável até 135 °C na autoclave	Sterilizabil la o temperatură de până la 135°C în autoclavă	Стерилизовать в автоклаве при температуре ниже 135 °C	Sterilizovateľné v autokláve pri teplote max. 135 °C	Primerno za sterilizacijo v avtoklavu do 135 °C	Steriliserbar vid upp till 135 °C i autoklav	ให้ดำเนินกระบวนการทำให้ปลอดเชื้อภายในเครื่องนึ่งฆ่าเชื้อไม่เกิน 135°C	Otoklavda 135°C'ye kadar sterilize edilebilir	Có thể tiệt trùng ở nhiệt độ tối đa 135°C trong nồi hấp	在高压灭菌器中进行灭菌处理，温度不高于135°C
	Dezynfekcja termiczna	Desinfeção térmica	Dezinfecție termică	Термическая дезинфекция	Tepelná dezinfekcia	Toplotno razkuževanje	Värmedesin-fektion	การฆ่าเชื้อด้วยความร้อน	Termal dezenfeksiyon	Khử khuẩn bằng nhiệt	可用于高温消毒
	Patrz instrukcja obsługi	Consulte o manual de instruções	Consultați manualul de instrucțiuni	См. инструк-ции по эксплуа-тации	Prečítajte si pokyny na používanie	Oglejte si navodila za uporabo	Se bruksanvis-ningen	ดูที่ คู่มือการใช้งาน	Kullanma kılavuzuna bakınız	Tham khảo sách hướng dẫn	请参考使用说明
	Autoryzowany przedstawiciel w Europie	Representante Autorizado CE	Reprezentant autorizat CE	Уполномочен-ный предста-витель ЕС	Splnomocnený zástupca ES	Pooblaščeni predstavnik ES	Auktoriserad representant inom EG	ตัวแทนที่ได้รับอนุญาตของทาง EC	AT Yetkili Temsilcisi	Đại diện ủy quyền tại cộng đồng châu Âu	欧共体授权代表

Composition of the GBT FLOWCONTROL® suction cannula


*Adapter (FV-128 - GBT FLOWCONTROL adapter) is sold separately and plugged on the cannula to allow the connexion with a 9 or 11 mm diameter suction tube instead of 16 mm.


How to insert the skirt on the GBT FLOWCONTROL® suction cannula?

How to remove the skirt from the GBT FLOWCONTROL® suction cannula?

Treatment recommendations - the GBT FLOWCONTROL® helps for lip or tongue retraction

Before using this product, please carefully read and follow these operating instructions. Please pay special attention to the safety precautions.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the country in which the user and/or patient is established.

Intended use

The GBT FLOWCONTROL® suction cannula is intended to be used with the GBT FLOWCONTROL skirt. It provides high-volume oral evacuation and aerosol reduction during dental procedures.

Indications

This product is indicated during dental procedure and specially indicated for prevention procedure such as GBT Guided Biofilm Therapy.

Intended user population

This product is intended to be used by dental professionals.

Intended patient population

This product is intended to be used on patients regardless of age, including children.

Undesirable side effects and contraindications

None.

Recommandations

⚠ Product must be cleaned, disinfected and sterilized prior to first use and between patients. Follow the reprocessing instructions of this manual.

Compatibility

The GBT FLOWCONTROL® suction cannula is intended to be used in combination with the GBT chair (Planmeca Compact i Touch Unit manufactured by Planmeca Oy).

Disposal

There is no specific safety risk related to the disposal of the product and its container. Disposal must be made according to country-specific regulations.


Product is made of PP (polypropylene) base and can therefore be recycled in this specific class (5) according to country-specific regulation.

Reprocessing instructions for GBT FLOWCONTROL® and Adapter

1 - At the point of use: Immediately after patient treatment, rinse the devices with cold tap water. Coarse soiling must be removed by wiping the products.

⚠ Special attention needs to be paid to cracks, crevices and hard-to-reach areas.

⚠ Dismantle the skirt from the body: **both parts must be sterilized separately**.

2 - The product must be cleaned manually or automatically by washer-disinfector.

- For automatic process: The washer-disinfector must comply with ISO 15883. Correctly place the product into a suitable rack and start the automated cleaning.

⚠ Instructions of the manufacturer of the washer-disinfector must also be observed.

⚠ Carefully follow the instructions provided by the disinfection solution manufacturer.

The following automated process[1] can be used to achieve an A0 Level of 3000:

2 min	Pre-clean with cold tap water[2]. Drain
5 min	Clean at 55°C with tap water and 0.5% neodisher MediClean Dental (Dr Weigert) solution. Drain
3 min	Rinse and neutralize with cold demineralized water[3]. Drain
2 min	Rinse with cold demineralized water[3]. Drain
3 min	Perform thermal disinfection (final rinse) with deionized water at 93°C minimum. Drain
20 min (at least)	Dry at 100°C.

- For manual process: The following validated EMS process can be used.

Validated product	Process	Process	
0,5% neodisher MediClean (Dr. Weigert)	ultraheli-vanniga[4]	Brush the devices in a cleaning solution of 0.5% neodisher in deionized water at 40°C with a suitable soft bristled brush until all visible residue has been removed.	
		10 min	Place the product in 0.5 % solution of neodisher with deonized water at 40 °C.
		10 s	Rinse the whole product under cold running deonized water
		At room temperature (15-25°C), use an air pistol (compressed air, max. 3 bars) to fully dry the whole product until no more residues of water are present (visible or detectable).	

Then disinfect:

⚠ ASP CIDEX® OPA solution must be used undiluted, within its useful and shelf life, and in compliance with the manufacturer's warnings and instructions for use.

5 min Immerse the product completely into CIDEX® OPA solution at minimum 20°C.

1 min (each rinse) The ASP CIDEX® OPA disinfectant requires a total of three rinses: keep the product totally immersed, and use a large volume[5] of fresh water.

Do not reuse the water for rinsing or any other purpose. Residues of disinfectant may cause serious side effects.

Then dry:

At room temperature (15-25°C), use an air pistol (compressed air, max. 3 bars) to fully dry the whole product until no more residues of water are present (visible or detectable).

3 - If stains are still visible on the product after cleaning/disinfection, the entire cleaning/disinfection procedure must be repeated. Product with visible damage, chip/flake loss, corrosion or deformation must be disposed of (no further use is allowed).

4 - Effective sterilization must take place only on fully dry products. Package your products disassembled with a single or double moist heat sterilization pouch. Avoid any heavy object interfering with the canulla during sterilization to reduce risk of deformation.

5 - Moist heat sterilization of products shall be performed according to ISO 17665 and in compliance with respective country requirements.

Parameters for the prevacuum moist heat cycle:

- 3 prevacuum phases
- Minimum temperature of 132°C - 3 min (minimum)
- 20 min (minimum) Dry

6 - Store the sterilized products at a temperature up to 40°C in a dry, clean and dust-free environment.

Technical description

Manufacturer	E.M.S. Electro Medical Systems S.A., Chemin de la Vuarpillière 31, 1260 Nyon, Switzerland
Essential Performance	This medical device, in the meaning of the EU Regulation 2017/745 has no Essential Performance
UDI code	(01) = GTIN reference (14 digits) (11) = manufacture date (YYMMDD) (19) = Lot number
Weight	8.2 g (cannula) 3 g (adapter)
Dimensions	120.5 x 28.5 x 26.5 mm (cannula) 54 x 20 mm (adapter)
Validated lifetime	250 sterilization cycles
Storage conditions	Temperature: max 40°C Humidity: max 85% relative humidity
Transport conditions	Not applicable
Operating condition	Temperature: 10°C to 35°C Humidity: 30% to 75%

Translations in this manual are based on the English version. Thus the English text will be referred to in case of discrepancy.

1 For example: Miele Professional G 7836 CD with Miele Rack E429
2 cold tap water = 16°C +/- 2°C
3 cold demineralized water = 20°C +/- 2°C
4 For example: Bandelin, Sonorex 1028K, 35kHz
5 e.g. 2 gallons

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