SybronEndo™

elements™ obturation unit

Obturation System

Instructions For Use

TABLE OF CONTENTS


English 2
Français 10
Español 17
Deutsch 24
Italiano 31
Dansk 38
Svenska 45
Ελληνικα 52
Nederlands 59
Norsk 66
Suomi 73
Português 80
Polski 87
Românâ 94
Magyar 101
Slovenský 108
Český 115
Eesti 122
Latviešu 129
Lietuvių 136
Slovensko 143
Български 150
Hrvatski 157

System B

Extruder

USER INTERFACE

Packaged Contents:	Quantity	Available for purchase
elements obturation unit Console	1	Not sold separately
elements Extruder Handpiece	1	973-0326
elements System B Handpiece	1	973-0327
elements Extruder Handpiece Shield	2	973-0331
elements System B Handpiece Shield	2	973-0332
elements obturation unit Power Supply	1	973-0333
AC Power Cord, US/Canada	1	973-0334
AC Power Cord, Europe	1	973-0335
AC Power Cord, UK	1	973-0336
AC Power Cord, Australia	1	973-0361
Cartridge Bender Tool	1	973-0337
elements obturation unit Instruction Manual	1	Not sold separately
elements Silicone Boot - Extruder Handpiece	2	973-0329
elements Silicone Boot - System B Handpiece	2	973-0330
Buchanan Heat Plugger (Fine .06 Taper)	1	952-0004
Quick Start Guide	1	Not sold separately

Accessories: Required to be used with the elements™ obturation unit.	Quantity	Available for purchase
elements Gutta Purcha Cartridge - 23 Gauge - Medium Body (silver)	10	972-1002
elements Gutta Purcha Cartridge - 25 Gauge - Light Body (silver)	10	972-1003
elements Gutta Purcha Cartridge - 23 Gauge - Heavy Body (silver)	10	972-1005
elements Gutta Purcha Cartridge - 23 Gauge - Medium Body (gold)	10	972-2500
elements Gutta Purcha Cartridge - 25 Gauge - Light Body (gold)	10	972-2501
elements Gutta Purcha Cartridge - 23 Gauge - Heavy Body (gold)	10	972-2502
Buchanan Heat Plugger - Extra Fine .04 Taper	1	952-0031
Buchanan Heat Plugger - Fine .06 Taper	1	952-0004
Buchanan Heat Plugger - Fine-Medium, .08 Taper	1	952-0005
Buchanan Heat Plugger - Medium .10 Taper	1	952-0006
Buchanan Heat Plugger - Medium-Large, .12 Taper	1	952-0007
Buchanan Heat Pluggers (EF, F, FM, M, ML)	5	952-0003
Hot Pulp Test Tip	1	973-0219

NOTE: Only accessories marketed by SybronEndo are to be used with the elements obturation unit.

FRONT PANEL

1.		On/Off Power Button
2.		Setup Button
3.		Extruder Plunger Return Home Indicator
4.		Extruder Plunger Return Home Button
5.		Extruder Speed Indicator
6.		Extruder Speed Button
7.		Extruder Heating Indicator
8.		Extruder No Heat Mode Button
9.		No Heat Indicator
10.		Extruder Preheat Button
11.		Extruder Temperature Indicator
12.		System B Temperature Setpoint Indicator
13.		Check Tip Indicator
14.		Temperature Adjustment Buttons
15.		Hot Pulp Test Mode Button
16.		Cautery Mode Button
17.		Backfill Mode Button
18.		Downpack Mode Button
19.		Hot Pulp Test Indicator
20.		Cautery Mode Indicator
21.		Backfill Mode Indicator
22.		Downpack Mode Indicator
23.		Volume Indicator
24.		Audio Mode Indicator
25.		Setup Mode Indicator
26.		Reset to Defaults Indicator
27.		Cartridge Empty Indicator

DEVICE DESCRIPTION

The elements™ obturation unit combines the System B™ Heat Source with the first cartridge-loaded, backfill extruder. The System B Handpiece provides heating of the Heat Plugger with precision controlled temperature and timing. The Backfill Extruder features single-use standard elements™ Gutta Percha Cartridges. The Obturation Unit's motorized system extrudes gutta percha under precise temperature and speed control.

INDICATIONS FOR USE

The elements™ obturation unit is intended to be used in dentistry to provide continuous heat at the tip of a dental instrument to backfill and downpack gutta percha during endodontic root canal treatment.

Additionally, when the unit is used in the Extruder Mode, the elements™ obturation unit is used only to backfill gutta percha during root canal obturation.

CONTRAINDICATIONS

- Not for use in the presence of flammable anesthetics such as a mixture of oxygen, air and nitrous oxide.
- The obturation unit should NOT be used on a patient with a pacemaker.
- Do not use on patients with a known sensitivity to copper and rubber latex.
- When using natural rubber latex gutta percha with this device, please check if patient is sensitive/allergic to natural rubber latex. Note: SybronEndo gutta percha is not made with natural rubber latex.

⚠ WARNINGS

- The System B tips and elements™ Gutta Percha cartridges are intended to reach high temperatures. Inadvertent contact with patient and operator must be avoided to prevent accidental burns..
- When replacing the extruder cartridge during a procedure, the cartridge may be hot to the touch. Handle empty extruder cartridge by the plastic locknut only. Other portions could be hot and could result in serious burns.
- Do not use the device near water. If unit is accidentally dropped in water, it might cause electric shock, destroy the device and can lead to serious injury or death.
- Only use the specified SybronEndo AC power supply and power cords.
- Use of other accessories that are not authorized for use with this device may cause malfunction and compromise patient safety.
- No modification of this equipment is allowed.
- For servicing refer to warranty information section.

⚠ **WARNING:** This product can expose you to chemicals including Lead and Lead compounds, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

PRECAUTIONS

- A protective dental dam is highly recommended when using this device.
- elements™ Gutta Percha Cartridges are for single-patient use.
- For use by qualified and trained personnel only.
- U.S. Federal law restricts this device to sale by or on the order of a healthcare professional.
- When not in use, remove main plug from wall socket to isolate unit from mains power.
- Do not use in teeth with immature and/or over instrumented root canal apices unless the apex has been sealed.
- Avoid multiple activations of the System B handpiece while operating inside the root canal since prolonged exposure to high temperatures may cause damage to the tissues surrounding the tooth (PDL or bone).
- While using the device in Hot pulp test mode, place a small amount of a gutta percha pellet on the hot pulp test tip prior to application to patient's tooth. The hot pulp test tip should not be placed directly on the tooth surface.
- Do not sterilize the Heat Plugger in a chemical sterilization unit as it might corrode its components.

⚠ **CAUTION:** This device has been tested and found to comply with EMC limits for the Medical Device Directive 93/42/EEC. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. The device generates radio frequency energy and, if not installed and used in accordance with these instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference with other devices, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between devices.
- Connect the device into an outlet on a circuit different from that to which the other device(s) is connected.
- For help with installation and configuration of the device contact Customer Service.

⚠ **CAUTION:** This device has been investigated with regard to safety from electrical shock and fire hazard as well as electromagnetic compatibility (EMC). The device has not been investigated for other physiological effects. Please contact SybronEndo if you have any further questions regarding electrical safety or EMC.

ADVERSE REACTIONS

None known.


STEP-BY-STEP INSTRUCTIONS

INITIAL SETUP

1. Plug the Extruder handpiece cord into the appropriate electrical connector on the right front panel of the unit by lining up the red dotted markings.
2. Plug the System B handpiece cord into the appropriate electrical connector on the left front panel of the unit by lining up the red dotted markings.
3. Plug the correct cord for your country's power into the elements™ Obturation Power Supply.
4. Plug the power cord into wall power.
5. Plug the power supply into the back of the unit.

POWERING UNIT ON

Power the unit ON by pushing the ⏻ button on the unit or the button on either of the handpieces for 2 seconds.

INSTALLING/REMOVING SYSTEM B AND EXTRUDER HANDPIECE SHIELDS


To install the System B or Extruder Handpiece Shield line up the index markings and carefully slide the shield onto the handpiece until a "click" is felt and heard. To remove the System B or Extruder Handpiece Shield, hold the bottom part of the handpiece with one hand, use other hand to grab the shield and pull away the handpiece shield without twisting.

SYSTEM B HANDPIECE

INSTALLING BUCHANAN HEAT PLUGGERS/HOT PULP TEST TIP

Install the desired Heat Plugger/ Hot Pulp Test Tip in the System B Handpiece as shown in the image. The Heat Plugger simply slides into one of 6 positions available by indexing the hex nut on the tip. Heat Pluggers without this hex nut should not be used with the elements obturation unit.

SELECTING SYSTEM B MODE

Press the desired mode button (for Downpack, Backfill, Hot Pulp Test, or Cautery). The temperature, temperature range, and maximum duration for each of these procedures is listed in the chart below. Heat Pluggers should be used for Downpack, Backfill, and Cautery Modes. Hot Pulp Test Tip should be used for Hot Pulp Test Mode.

Setting	Symbol	Preset Temperature (°C)	Temperature Range (°C)	Duration (seconds)
Hot Pulp Test		200	140-400	60
Cautery		600	600	5
Backfill		100	100	15
Downpack		200	140-400	4

After selecting the mode, the unit will display the mode selected and its preset temperature.

CHANGING TEMPERATURE

NOTE: Temperature can only be changed in Downpack and Hot Pulp Test Mode

If a temperature setting other than the preset is desired, the (+) or (-) button can be used to change the temperature of the heat plugger tip in 5 degree increments. Holding either (+) or (-) button for two seconds will scroll the temperature until the desired temperature is obtained as indicated on the display. Pressing and holding the current mode buttons for 4 seconds will save that temperature as a preset for that specific mode. This new temperature preset is retained until manually changed or defaults are renewed. If a button for another mode is selected before saving the new temperature, the unit will switch to that mode and the new temperature will not be saved.


ACTIVATING SYSTEM B HANDPIECE

Activate the System B Handpiece extra orally prior to use by depressing and holding its button. If the Heat plugger is installed correctly, the LED indicator on the Handpiece will illuminate, and the appropriate mode indicator will flash rapidly. The Heat Plugger will only remain heated as long as the button is depressed. To avoid overheating, each of the modes has a preset duration, as shown in the table above. The Handpiece will need to be reactivated to resume heating beyond the preset duration. The automatic sensing circuit will display if an incorrect tip is installed (relative to the selected mode) and will not function until correct tip is installed. If a worn-out tip is encountered, will also be displayed.


NOTE: Heat pluggers must be autoclaved between every use. See the instructions in the "Cleaning and Sterilization" section.

EXTRUDER HANDPIECE

EXTRUDER HANDPIECE - INSERTING AND REMOVING CARTRIDGE

- Prior to inserting the cartridge press the button to calibrate the handpiece.
- Insert a new obturation cartridge into the handpiece by pushing the larger end of the cartridge into the handpiece and turning it a quarter turn clockwise so that it engages into the handpiece. Remove the cartridge by turning the locknut counterclockwise in order to dislodge it. Then pull the cartridge out of the handpiece. The cartridge can then be disposed.
- To replace a cartridge before it is completely empty, press the button, which will retract the plunger. During plunger retraction, both indicators will flash on the screen and the indicator will remain solid after complete plunger retraction.

NOTE: When the extruder handpiece is disconnected from the console, the temperature Indicator will display this .

WARNING: Handle empty extruder cartridge by plastic locknut only. Other portions could be hot and result in a serious burn.

ACTIVATING EXTRUDER HEATING


Press the button or the activation button on the extruder handpiece. Segments of will flash sequentially on the screen during preheat, which should take approximately 45 seconds. The thermometer segments on the screen will fill as temperature rises. When preheat is complete, the and indicators on the screen will be solid. The temperature will be maintained 15 minutes for gutta percha. Heating will then shut off, the indicator will disappear and the handpiece will not activate. The heater can also be shut off at any time by pressing the button. Pressing the button displays on the screen. To exit the no heat mode press Extruder No Heat Mode Button again. The thermometer segments will track actual temperature. The heater can be reactivated at any time by pressing the button or Extruder Handpiece activation button again. Pressing the button while heater is on will disengage heat.

SETTING EXTRUSION SPEED

SLOW – Set extrusion speed using the button. Ensure appears on the screen. When preheat is complete, pressing top or bottom handpiece activation button will cause slow extrusion speed.
MEDIUM - Set extrusion speed using the button. Ensure appears on the screen. When preheat is complete, pressing bottom handpiece activation button will cause medium extrusion speed.
FAST- Set extrusion speed using the button. Ensure appears on the screen. When preheat is complete, pressing top handpiece activation button will cause fast extrusion speed.

EXTRUDING GUTTA PERCHA

Activate the extruder by depressing any of the two buttons on the handpiece. The extruder motor will initiate, which is indicated by a lighted LED on the handpiece and a flashing of . Before placing the extruder needle in the root canal, activate the extruder with the button on the handpiece until material extrudes through the needle. Upon release of the button, gutta percha will stop extruding.The extruder can then be placed into the root canal for obturation.

If the handpiece is deactivated before the cartridge is empty, the plunger will retract slightly to prevent excess material from discharging from the needle. When the cartridge is empty, will display and the plunger will automatically retract. The Handpiece has a mechanical "gauge" which indicates the amount of gutta percha remaining in the cartridge as shown in the image.

CHANGING THE DEFAULT SOUND SETUP

The elements™ obturation unit offers three different sound mode settings. Active sounds associated with various events for each of these modes are displayed in the chart below. To cycle through volume indicator , audio mode indicator 3 and Setup Mode Indicator • options press the button repeatedly until the desired option starts flashing.

1. Activate the unit by pushing the button.
2. Press the button, the will flash on the screen. Press one of the (+) or (-) buttons until desired volume level is reached.

Press the button to enter Sound Mode Setup, 1, 2 or 3 will flash. Press one of the (+) or (-) buttons until desired sound mode is reached. The sound mode chart outlines which sounds are active in each mode. If no sound is desired, adjust the volume accordingly.

NOTE: Press and hold button for 2 seconds and the desired volume level or sound mode will be saved. Unit will also leave the sound set up mode and return to functional mode after 5 seconds of inactivity.

SOUND MODE CHART

EVENT	SOUND NAME — ACTIVE IN MODES:	1	2	3
1	**Any Keystroke** Audio feedback for any user interaction with the console and handpiece.	✓		
2	**Extruder Set to No Heat** A chime when no heat button on the button pad membrane is pressed.	✓		
3	**Initiate Extruder Preheat** A chime when the preheat button on the button pad membrane is pressed.	✓		
4	**Extruder Preheat Timeout** A chime to alert user when the extuder keep warm timer has expired. (Standard Gutta Percha mode 15 minute timeout)	✓		
5	**Invalid Keystroke** Audio prompt for invalid buttons are presses in setup mode.	✓	✓	
6	**System B Downpack 5- and 10-second Timers** In Downpack mode after the auto shutoff or release of the handpiece button, the system has two timers. The system beeps once at 5 seconds and beeps twice at 10 seconds upon heater shutoff (by releasing button or auto shutoff).	✓		
7	**Extruder Reaches Temperature** A «zing» sound at the end of preheating cycle notifying the users that the exturder has reached the target tempearture.	✓	✓	
8	**Extruder Mode Active** A humming sound when the extruder motor is active.	✓	✓	
9	**Empty Extruder Cartridge** A «ding» sound when the extruder motor is at the most forward postion notifying the user that the cartridge is empty now.	✓		
10	**Unit Power On** Sound when the «On» button is pushed.	✓	✓	✓
11	**Reset Defaults** Restore factory setting ackowledgement sound.	✓	✓	✓
12	**Incorrect Heat Plugger** System B alert for defective tips in Downpack, Backfill and Cautery mode, or Hot Pulp Test Tip missing.	✓	✓	✓
13	**New Temperature Saved** When the System B side default temperature is changed and saved.	✓	✓	✓
14	**Extruder Is Under Temperature** A high pitch chime when the user tries to activate the extruder handpiece before the tip reaches the minimum required temperature.	✓	✓	✓
15	**Active System B Handpiece** A humming sound when the System B heater is active.	✓	✓	✓
16	**System B Handpiece Timeout** A «Clink» sound coinciding with either the release of the Handpiece Activation button or the end of a System B mode duration.	✓	✓	

RESETTING TO DEFAULTS

To reset to defaults, press the button repeatedly until the • indicator flashes at the bottom right of the screen. Immediately press the + on the button once. The • indicator will change to ▸•◂ indicator and flash. Immediately press and hold the button for 2 seconds and the unit will reset to defaults, including sound mode, volume and preset temperatures.

Unit will leave setup mode and return to functional mode after five seconds of inactivity. When the unit leaves setup mode, 3 , , and ▸•◂ will all disappear.

POWER UNIT OFF

To power the unit off, press the button.

CLEANING AND STERILIZATION

Clean the elements™ obturation system, Buchanan Heat Pluggers, Hot Pulp Test Tip, Cartridge Needle Bender and Extruder and System B Shields immediately after each patient use following the cleaning steps to prevent drying of soil and contaminants.

The components mentioned above should be inspected prior to cleaning for the appearance of defects such as cracks, deformations, corrosion, which are indications that the instruments are not in conditions to be re-used with the required level of confidence.

The use of automated cleaning devices is not recommended for the cleaning of the above components.

- Remove the used gutta percha cartridge and the Heat Plugger from the respective handpieces.
- Carefully discard examination gloves, rinse and disinfect hands using an appropriate hand disinfectant solution and use a new pair Lof examination gloves.
- Follow standard precautions for personal protection using cleaning agents/disinfectant solution as recommended by the manufacturer.

NOTE: FDA-cleared sterilization pouches are to be used following the recommended sterilization parameters.

Health Care facilities are responsible for making sure that the sterilization equipment is calibrated according to the manufacturers manuals and specifications. In addition, health care facilities are responsible for training their staff on infection control and proper sterilization procedures.

elements™ obturation CONSOLE

The elements™ obturation Console has sensitive electronic circuitry and should not be autoclaved as it will damage the device.

CLEANING:

- The Base Unit should be wiped down with a cloth moistened with any of the common nonabrasive, mild, cleaning liquids containing alcohol. SybronEndo recommends Cavicide.
- Repeat with additional wipes until no soil is visible on the cloth.
- Visually inspect the Base Unit to ensure cleanliness and for any damage or contamination. If any contamination is visible, repeat the cleaning steps.

⚠ Do not submerge the device or let liquids enter the Base Unit.

EXTRUDER AND SYSTEM B HANDPIECES

The Extruder and System B Handpieces have sensitive electronic circuitry and should not be autoclaved as it will damage the device.

CLEANING:

- The handpiece shields must be removed from the handpieces prior to cleaning.
- Clean handpieces and cords with a cloth moistened with any of the common nonabrasive, mild, cleaning liquids containing alcohol. SybroEndo recommends CaviCide.
- Repeat with additional wipes until no visible soil is detected on the cloth.
- Visually inspect the Extruder and System B handpieces to ensure cleanliness and for any damage or contamination. If any contamination is visible, repeat the cleaning steps.

⚠ Do not submerge the device or let liquids enter the enclosure.

BUCHANAN HEAT PLUGGERS AND HOT PULP TEST TIP

Heat Pluggers and Hot Pulp Test Tips must be cleaned and sterilized after each use. For the reuse of the Heat Pluggers, please refer to the Instructions For Use enclosed in the individual Heat Pluggers.

CLEANING:

- Remove the Heat Plugger/Hot Pulp Test Tip from the System B Handpiece.
- Clean with water and mild, nonabrasive detergent such as dishwashing liquid using a soft bristled brush until all visible soil is removed.
- Rinse with clean running water for a minimum of 30 seconds.
- Visually inspect for cleanliness and for any damage or contamination.
- If any contamination is visible, repeat the cleaning steps until there is no visible contamination.
- Dry with a clean, lint free cloth until there is no visible moisture.

STERILIZATION:

- Place the Heat Plugger/Hot Pulp Test Tip in a sterilization pouch suitable for steam sterilization.
- For Heat Plugger/Hot Pulp Test Tip operate Gravity Displacement Sterilization Cycle at 121 °C (250 °F) for 30 minutes followed by 30 minutes of drying time or Pre-Vacuum sterilization cycle 132 °C (270 °F) for 4 minutes, followed by 20 minutes of drying time as noted in the table titled: Recommended Sterilization Parameters.

EXTRUDER AND SYSTEM B SHIELDS

The Extruder and System B shields must be cleaned and sterilized after each use.

CLEANING:

- The shields must be removed from the handpieces for cleaning and sterilization.
- Note: Do not remove the Silicone Boot from either handpiece shield.
- The Extruder and System B shields should be wiped down with a cloth moistened with any of the common nonabrasive, mild, cleaning liquids containing alcohol. SybronEndo recommends CaviCide.
- Use a cotton swab or a small, soft brush moistened with CaviCide to remove any soil that may have accumulated in crevices (e.g., between the silicone boot and the body of the shield)
- Repeat with additional wipes until no soil is visible on the cloth. After cleaning soil, the shields should be rinsed with clean running water for a minimum of 30 seconds to remove residues of cleaning agents.
- Visually inspect the Extruder and System B Shields to ensure cleanliness and for any damage or contamination. If any contamination is visible, repeat the cleaning steps.
- Dry with a clean, lint free cloth until there is no visible moisture.

STERILIZATION:

- Place the handpiece shields in a sterilization pouch suitable for steam sterilization. For handpiece shields, operate Gravity Displacement Sterilization Cycle at 121 °C (250 °F) for 30 minutes followed by 30 minutes of drying time or Pre-Vacuum sterilization cycle 132 °C (270 °F) for 4 minutes, followed by 20 minutes of drying time as noted in the table titled: Recommended Sterilization Parameters.

CARTRIDGE NEEDLE BENDER

The Cartridge Needle Bender must be cleaned and sterilized after each use.

CLEANING:

- Prepare an ultrasonic bath with an enzymatic cleaning solution at the concentration and temperature specified by the manufacturer of the enzymatic solution.
- Place Cartridge Needle B Hender in the ultrasonic bath for 10 minutes.
- Remove the Cartridge Needle Bender from the bath and rinse with clean running water for a minimum of 30 seconds
- Visually inspect for cleanliness and for any damage or contamination.
- If any contamination is visible, repeat the cleaning steps until there is no visible contamination.
- Dry with a clean, lint free cloth until there is no visible moisture.

STERILIZATION:

- Place the Cartridge Needle Bender in a sterilization pouch suitable for steam sterilization. For Cartridge Needle Bender operate Gravity Displacement Sterilization Cycle at 121 °C (250 °F) for 30 minutes followed by 30 minutes of drying time or Pre-Vacuum sterilization cycle 132 °C (270 °F) for 4 minutes, followed by 20 minutes of drying time as noted in the table titled: Recommended Sterilization Parameters.

Recommended Sterilization Parameters

Moist Heat Sterilization		
Cycle	**Gravity Displacement***	**Pre-Vacuum**
Temperature (°C)	121 °C (250 °F)	132 °C (270 °F)
Exposure Time (minutes)	30 minutes	4 minutes
Drying Time (minutes)	30 minutes	20 minutes

*The sterilization process of the Shields and Needle Benders were validated within fully loaded chamber using steam gravity displacement sterilization cycle.

STORAGE:

- After sterilization, place the FDA-cleared sterilization pouches containing the devices in a dry and dark place such as a closed cupboard or a drawer.
- Follow the instructions provided by the FDA-cleared sterilization pouches' manufacturer regarding storage conditions and maximum allowed time in storage.

SYMBOLS

SN	Serial Number
	Manufacturer
	Manufacture Date
	This way up
	Fragile, handle with care
	Keep dry
	Temperature limits for storage and transportation
	Humidity Limitation
	Pressure Limitation
	Do not throw into trash
C US MC 157498	CSA mark with «C/US» indicator for certified products
	Warning or Precautions
	Consult instructions for use
	Type B applied part
	Class II Medical Electrical Equipment
	Do not autoclave
REF	Part Number
121°C (250°F)	Gravity Displacement Sterilizable in a steam sterilizer (autoclave) at 121 °C (250 °F).
132°C (270°F)	Pre-Vacuum Sterilizable in a steam sterilizer (autoclave) at 132 °C (270 °F).
Rx only	CAUTION: USA federal law restricts this device to sale by or on the order of a dentist.

DISPOSAL INSTRUCTION

Dispose of Properly

When you plan to dispose of the elements™ obturation unit please return it to SybronEndo for environmentally sound recycling. Please contact your sales representative for details. The Waste Electrical and Electronic Equipment (WEEE) Directive forbids the disposal of waste electrical and electronic equipment as unsorted municipal waste, requiring that it be collected and recycled or disposed of separately. Hazardous substances are present in electrical and electronic equipment and present potential risks to human health and the environment. Electrical equipment must be disposed according to local regulations and regional laws.

Origin: SybronEndo
200 S. Kraemer Blvd.
Brea, CA 92821 USA
1-800-KERR-123
www.sybronendo.com

TECHNICAL SPECS

DIMENSION/WEIGHT

elements™ obturation unit Console:	Weight:	0.93 kg
	Height:	153 mm
	Width:	179 mm
	Depth:	133 mm
Extruder Handpiece:	Weight:	0.23 kg
	Outside Diameter:	24 mm
	Length:	193 mm
System B Handpiece:	Weight:	0.14 kg
	Outside Diameter:	19 mm
	Length:	173 mm

Classification: IEC 60601-1, Class II, patient contact, type B

Mode of operation: Continuous operation

Degree of safety of application in presence of a flammable anesthetic mixture: Unit unsuitable for use in presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Storage and Transport conditions:	Conditions for use:
40°C (104°F) –10°C (14°F) Temperature	35°C (95°F) 10°C (50°F) Temperature
95% 10% Relative Humidity	75% 30% Relative Humidity
1400 hPa 500 hPa Atmospheric pressure	1013 hPa 697 hPa Atmospheric pressure
	0–10,000 ft. Altitude Limitation

Model: elements™ obturation unit

Power Supply: INPUT: 100-240V ~ 47-63Hz 1.35-0.53A
OUTPUT: 24V ⎓ 2.08A max.
Timing Mechanism: ±1 second

WARRANTY INFORMATION

RETURNS AND WARRANTY

The unit is designed for lasting durability and predictable results. The unit comes with a one-year warranty. If the unit fails to operate correctly, contact your SybronEndo customer service representative at 800-346-3636. Outside the United States, call customer service at +1-714-516-7979 or contact your local dealer or supplier to coordinate the returns/repairs with SybronEndo. Repair services are available after one-year warranty period at the customer's own cost. The same RMA process can be followed.

Warranty repairs can only be performed by SybronEndo or by authorized agencies with original factory parts. Authorized service personnel will have access to any necessary technical specifications to ensure the repair is performed correctly. Any unauthorized repairs will void the warranty.

1. Many problems can be solved over the phone. Do not hesitate to contact us if you experience difficulties when using the Unit.
2. For returns, call us for a Return Merchandise Authorization (RMA) number.
3. Use the original packaging when shipping or storing your Unit. If the original box and/or foam packing are not available, call SybronEndo for packing instructions. Damages incurred during shipping due to improper packaging are not covered by the warranty.
4. Mark all outer packaging with the RMA number, your name, address, and phone number.
5. Shipping costs are not covered by the warranty.

DISCLAIMER

The responsibility of SybronEndo, as a manufacturer of Medical Electrical devices, extends to the technical safety features of the device only if maintenance and repairs are carried out by SybronEndo or by authorized agencies with original factory parts.

For safety reasons, this product should be used with accessories manufactured and sold by SybronEndo. Any use of non-authorized accessories or not following any of the instructions for use is done so at the operator's risk and voids the warranty.

SybronEndo does not assume any responsibility for incorrect diagnosis due to operator error or equipment malfunction.

ELECTROMAGNETIC COMPATIBILITY

Medical Electrical Equipment requires special precautions regarding electromagnetic compatibility. The elements™ obturation unit needs to be installed and put into service according to the electromagnetic compatibility information provided in these instructions. Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

⚠ WARNINGS:

- The use of accessories and options other than those specified by the manufacturer may result in increased emissions or decreased immunity of the elements™ obturation unit.
- In the event of intermittent function of the device due to exposure to electromagnetic interference, power cycle the unit as per IFU instructions provided above.
- In the event of permanent damage or display blackouts of the device due to electromagnetic interference, stop use of the unit. Contact Customer Service.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the elements™ obturation unit, including compatible cable. Otherwise, degradation of the performance of this equipment could result.

Guidance and manufacturer's declaration - electromagnetic emissions

The elements™ obturation unit is intended for use in the electromagnetic environment specified below. The customer or the user of the elements™ obturation unit should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The elements™ obturation unit uses RF energy only for its internal function. The Element Obturation Unit delivers energy to the patient in a form other than RF electromagnetic.
RF Emissions CISPR 11	Class B	The elements™ obturation unit is suitable for use install establishments with public low-voltage network that supplies buildings used for domestic purposes and professional healthcare facility environ-ment (e.g. hospital) with dedicated power supply system
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity

The elements™ obturation unit is intended for use in the electromagnetic environment specified below. The customer or the user of the elements™ obturation unit should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conducted disturbances induced by RF fields IEC 61000-4-6	3* Vrms 0.15 MHz – 80 MHz 80% AM at 1 kHz	3* Vrms 0.15 MHz – 80 MHz 80% AM at 1 kHz	The distance between transmitting antenna and EUT was at least 1 meter for proximity field RF wireless communication equipment testing. The testing was performed with 1.7 m System B handpiece cable, 1.7 m extruder cable and 1.85 m AC/DC adapter cable bundled to 1 meter. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EM Fields IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz, 80%	3 V/m 80 MHz – 2.7 GHz, 80% AM at 1 kHz	
Proximity Fields from RF Wireless communications equipment IEC 61000-4-3	Table 9 of IEC 60601-1-2 ed. 4.0	Table 9 of IEC 60601-1-2 ed. 4.0	

NOTE 1: *6 Vrms in the ISM bands
The ISM band is 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz
NOTE 2: The elements™ obturation unit was tested at 120 VAC/60 Hz and 230 VAC/50 Hz
NOTE 3: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity

The elements™ obturation unit is intended for use in the electromagnetic environment specified below. The customer or the user of the elements™ obturation unit should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ± 15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	The elements™ obturation unit complies antistatic material and synthetic material with relative humidity level as low as 10%.
Electrical fast transient/burst (Power Lines) IEC 61000-4-4	±2 kV, 100 kHz repetition frequency	±2 kV, 100 kHz repetition frequency	The elements™ obturation unit is suitable for use in all establishments with public low-voltage network that supplies buildings used for domestic purposes and professional healthcare facility environment (e.g. hospital) with dedicated power supply system.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line(s) to line(s) ±0.5 kV, ±1 kV, ±2 kV line(s) to ground	±0.5 kV, ±1 kV line(s) to line(s) ±0.5 kV, ±1 kV, ±2 kV line(s) to ground	The elements™ obturation unit is suitable for use in all establishments with public low-voltage network that supplies buildings used for domestic purposes and professional healthcare facility environment (e.g. hospital) with dedicated power supply system.
Voltage dips and Voltage interruptions on power supply input lines IEC 61000-4-11	Voltage dips: 0% UT for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT for 1 cycle and 70% UT for 25/30 cycles at single phase (0°) Voltage interruptions: 0% UT, 250/300 cycle	Voltage dips: >95% UT for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. >95% UT for 1 cycle and 70% UT for 25/30 cycles at single phase (0°) Voltage Interruptions: 100% UT for 250 cycles	The elements™ obturation unit is suitable for use in all establishments with public low-voltage network that supplies buildings used for domestic purposes and professional healthcare facility environment (e.g. hospital) with dedicated power supply system.
Rated power frequency magnetic fields IEC 61000-4-8	30 A/m; 50 Hz or 60 Hz	30 A/m; 50 Hz and 60 Hz	The elements™ obturation unit is suitable for use in all establishments with public low-voltage network that supplies buildings used for domestic purposes and professional healthcare facility envi-ronment (e.g. hospital) with dedicated power supply system.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

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Kerr Australia Pty. Ltd.