**SYMBOL DESCRIPTION**
Full explanation of symbols used on Kerr packaging is located at: **http://www.kerrdental.com/symbols-glossary**

| | |
|---|---|
| | Caution |
| | Manufacturer |
| 2797 | European Conformity |
| EC REP | Authorized representative in the European Community |
| | Consult instructions for use |
| | Do not re-use |
| LOT | Batch code |
| | Catalogue number |
| MD | Medical device |
| 2°C 28°C (36°F 82°F) | Temperature storage limitations: low 2°C/36°F and high 28°C/82°F |
| | Keep away from sunlight |
| | Importer |

**CONTENTS**

1. **Syringe (4 g)**
2. **Unidose (0.2 g)**
3. **Syringe (1.7 g)**
4. **Bottle (5 ml)**
5. **Syringe (3 g)**
6. **Dispensing tip**
7. **Applicators**
8. **Mixing wells**
9. **Polisher**

**ADVERSE EVENTS**
If a serious incident occurs with this medical device, report it to the manufacturer and to the competent authority for that country in which the user and/or patient is established.

**STEP-BY-STEP INSTRUCTIONS**
1. Clean the tooth with pumice and water to remove any debris or stain that may interfere with the shade selection process.
2. Select the desired shade. A composite mock-up or VITA-matched shade guide may be used to check the color of the tooth to be restored and any adjacent teeth. In many cases, a single shade of Herculite Ultra in enamel opacity will be sufficient to restore the tooth. The shade of the restorative material is affected by the thickness of the composite. Shade selection and mock-ups should match the thickness of the tooth to be restored.
3. After shade selection, isolate the tooth with your preferred isolation method.
4. Prepare the tooth using your preferred technique.
5. Complete the bonding procedure. Herculite XRV Ultra is compatible with all methacrylate dental adhesive systems. Herculite XRV Ultra kits include Optibond Solo Plus.
6. If desired, place a liner with a low viscosity flowable composite, such as Herculite Ultra Flow and light cure according to manufacturer's instructions.
7. Herculite XRV Ultra features optimal, non-sticky handling. It permits the use of hand instruments to define most of the restoration anatomy, which will help minimize the time and effort associated with rotary instrumentation.
8. For Unidose delivery: Insert Unidose tip into any composite dispenser. Unit-dose should be firmly seated in a dispenser. Remove the cap and inject composite material from the Unidose tip directly into the cavity, starting with the deepest portion.
For Syringe delivery: After extruding material onto the pad, back off syringe at least 2 full turns to relieve pressure on the material. Recap the syringe immediately after use.
9. Apply composite layers up to 2 mm in thickness for enamel/translucent (incisal) shades or 1.5 mm in thickness for dentin shades, taking care to adapt the material to the cavity surface. Avoid exposure of Herculite Ultra to intense operatory light during shaping of the material.
10. Cure each layer. Refer to the table at the end of these instructions for curing recommendations.
11. After curing, complete the finishing and polishing procedure. Adjust occlusion as needed.

**RECOMMANDED CURING TIMES**
Herculite XRV Ultra is intended to be cured by an LED light with a minimum light intensity of 1,000 mW/cm² and a light wavelength output within the 400–520 nm range. Please refer to the table below for Kerr curing recommendations. For other brands of curing lights, follow their manufacturing recommendations. In any posterior restoration, Kerr recommends additional 10 second cures on both the buccal and lingual surfaces of the tooth after removal of the matrix band.

| Curing Time | Dentin Shades | Enamel/Translucent (Incisal) Shades |
|---|---|---|
| Incremental Depth | 1.5 mm | 2.0 mm |
| Demi™ Plus (or light with output > 1,000 mW/cm²) | 10 seconds | 10 seconds |

**CLEANING AND DISINFECTION**
Unidose Tips: Do not re-use Unidose tips.
Unidose tip guns: Refer the reprocessing instructions for the manufacturer of the tip gun you are using. Composite syringes: To prevent syringes from getting contaminated with saliva or other body fluids, use protective barriers. Always handle the syringe with clean gloves. Do not reuse if it becomes contaminated.

**STORAGE AND DISPOSAL**
Product should be stored between 2°C to 28°C (36°F to 82°F). Avoid excess heat. If product is stored in a refrigerator, allow the product to reach room temperature prior to use. Every effort should be made to protect the composite from visible light prior to use.
The expiration date is marked on packaging. Do not use restorative after expiration date.
For proper disposal, follow federal or local laws applicable to that country.

**LIMITED WARRANTY - LIMITATION OF KERR'S LIABILITY**
Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr'sliability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages. EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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