Instructions for Use

K-Files™, K-Flex™, Triple-Flex™, Hedstrom™, and Reamers

Hand Files

Kerr Italia S.r.l.
CE 2797 Via Passanti, 174
Scafati (SA), 84018 Italy
+39 081 850 8311
kerrdental.com

Distribution Facility:
Kerr Australia Pty. Ltd.
Unit 6, 12 Mars Road
Lane Cove West
NSW 2066, Australia
+61 2 8870 3000

077-1129e_Rev.00
2026-06

UK CA 0086 Authorised Rep in UK:
Kerr UK Ltd.
Orega Stockley Park
4 Longwalk Road
Stockley Park
Uxbridge UB 11 1FE
United Kingdom

©2026 Kerr Corporation. All Rights Reserved.

Kerr™ ENDODONTICS

SYMBOLS DESCRIPTION

Full explanation of symbols used on Kerr packaging is located at:
http://www.kerrdental.com/symbols-glossary

121°C (250°F)	Gravity displacement Sterilizable in a steam sterilizer (autoclave) at 121 °C (250 °F).
132°C (270°F)	Pre-Vacuum Sterilizable in a steam sterilizer (autoclave) at 132 °C (270 °F).
134°C (273°F)	Pre-Vacuum Sterilizable in a steam sterilizer (autoclave) at 134 °C (273 °F).

en - ENGLISH

K-Files™, K-Flex™, Triple-Flex™, Hedstrom™, and Reamers
Hand Files

GENERAL INFORMATION

Endodontic hand files are used for root canal cleaning and shaping during endodontic treatment.

K-Files, K-Flex, and **Triple-Flex** files are used for the negotiation of canals and for the initial enlargement of canals to make way for rotary instrumentation increasing flexibility to negotiate with curved canals and stressful situations.
Hedstrom files are used for rapid enlargement of canals because of the aggressive cutting. They are also used for the removal of root canal filling materials e.g., gutta percha during retreatment of the root canal.
Reamers are used for widening the root canals to clean and shape canals using a rotary motion providing flexibility with integrity.

All hand files are available in lengths of 21, 25 and 30 mm, in packets of 6 files and in sizes ranging from 06 to 80.

Intended Purpose: Endodontic hand files are used for root canal cleaning and shaping during endodontic treatment.
Intended Users: This product is to be used by qualified dental professionals.
Intended Patient Population: The device is for use in a general population of patients in need of a root canal treatment.
Clinical Benefit: Allows the dentist to locate and prepare the root canal for rotary instrumentation.

PRODUCT COMPOSITION

These instruments contain Stainless Steel, Polyetherimide, printing inks, and non-latex Silicone.

INDICATIONS FOR USE

Hand Files are instruments used for shaping, cleaning, and enlarging root canals for injured pulp conditions.

CONTRAINDICATIONS

None known.

WARNINGS

- Hand Files are non-sterile products and must be cleaned and sterilized prior to first use to prevent possible biohazard exposure.
- These products contain stainless steel and should not be used for individuals with known allergic sensitivity to this metal. The dental professional should inform the patient that stainless steel is known to cause allergic reactions.
- Use of a rubber dam is recommended when using these instruments for root canal procedures to avoid accidental aspiration or ingestion and avoid injury to the patient's oral soft tissues.

PRECAUTIONS

- Inspect the packaging and the instrument before use: if the packaging or the instrument are damaged, do not use and discard it.
- Determining the working length with the Endo Stop is imperative to ensure proper instrumentation.
- Ensure straight-line access and do not force a file down the canal as it might separate, get stuck, create transportation an/or ledge formation.
- Ensure that excessive stress is not placed on file during instrumentation as it might separate inside the root canal.
- Ensure that hand files are used in the correct sequence to avoid ledge formation, canal transportation or file separation.
- Ensure the final file used has a right size; it is not too large or too small to avoid excessive or insufficient instrumentation as it might cause clinical failure due to perforation, canal transportation or presence of infected pulp tissue left inside the root canal.
- Ensure the file is used in a root canal with proper irrigation/lubrification.
- Use personal protective equipment such as gloves and safety glasses.

ADVERSE REACTIONS

If a serious incident occurs with this medical device, report it to the legal manufacturer, Kerr Italia, S.r.l., and to the competent authority for that country in which the user and/or patient is established.

STEP-BY-STEP INSTRUCTIONS

Limitations on Use:

- Hand Files can be used a maximum of five (5) times.
- For extremely curved canals, use a new file for each canal.
- When the file shows: untwisted threads, dull cutting surfaces, bent instrument, signs of corrosion, loss of color marking and unable to sterilize the instrument prior to use.

Hand Files Technique

a. Obtain straight-line access.
b. Locate the canal orifices and obtain coronal patency.
c. Begin crown down by taking a coronal shaper to resistance.
d. Establish working length using hand files and apex locator. A radiograph may be taken to assist in working length determination.
e. In more difficult, curved, and narrow canals you must go through several recapitulations of the sequence.

Clinical Tips

- Irrigate the canal using NaOCl, 17% EDTA or a suitable endodontic lubricant agent to facilitate debris removal.
- Use light pressure and avoid force. When the instrument no longer advances apically, proceed to a smaller file.
- Prior to use, check instruments for defects; including deformation, breakage, fractures, corrosion, handle damage, loss of color coding or size marking. Discard any product with defects.
- If an instrument failes to function normally, it must not be used.

- Irrigate the canal frequently and abundantly to facilitate effective debris removal.

CLEANING AND DISINFECTION

NOTE: Hand Files can be used a maximum of five (5) times and must be cleaned and sterilized prior to first use and after each patient treatment.

- The instructions provided in the section "Cleaning and Sterilization Instructions" have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse.
- Health Care facilities are responsible for calibration of sterilization equipment and training of staff on infection control and sterilization according to the manufacturers' instructions.

Cleaning

NOTE: Clean Hand Files immediately after each patient use to prevent drying of contaminants and cross contamination between patients. The use of automated cleaning devices is not recommended.

1. Clean the hand files in an ultrasonic device for ten (10) minutes using a neutral-pH enzymatic ultrasonic detergent (Enzol® Enzymatic Detergent, or equivalent) is recommended. Follow the detergent manufacturer's instruction for the dilution, water quality, temperature and how to prepare the ultrasonic bath. Follow the ultrasonic bath manufacturer's instruction for frequency and power density for cleaning.
2. If any contamination is visible, repeat the cleaning steps until there is no visible contamination.
3. Rinse thoroughly for a minimum of 2 minutes with warm, tap water.
4. Allow the devices to air dry on a clean, disposable lint-free cloth.
5. Inspect the files for any signs of damage. Discard files if there are visible signs of deformation, breakage, fractures, corrosion, handle damage, loss of color coding or size marking. Discard any product with defects.

Sterilization

Note: Do not sterilize the Hand File in a chemical sterilization unit as it might corrode.

1. Place the Hand File(s) in a FDA-cleared or ISO 11607 compliant sterilization pouch.
2. The following moist heat sterilization cycles can be used:
 a. **Gravity Displacement** sterilization cycle at 121 °C (250 °F) for 30 minutes followed by 30 minutes of drying time.
 b. **Pre-Vacuum** sterilization cycle 132 °C (270 °F) for 4 minutes, followed by 20 minutes of drying time.
 c. **Pre-Vacuum** sterilization cycle 134 °C (273 °F) for 3 minutes, followed by 20 minutes of drying time.
3. After sterilization, inspect the pouch for damage and, in case of a perforated pouch, do not use the product as the product's sterilization could be compromised. RE-pack the product and repeate the sterilization procedure.

STORAGE AND DISPOSAL

After sterilization, place the pouches containing the files in a clean dry environment. Follow the pouch manufacturer's instructions for storage conditions and maximum storage time. Follow local and national regulations for disposal of sharps.

LIMITED WARRANTY - LIMITATION OF KERR'S ITALIA LIABILITY

Kerr Italia's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr Italia's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr Italia's product. In the event of a defect in material or workmanship, Kerr Italia's liability is limited, at Kerr Italia's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr Italia. In no event shall Kerr Italia, be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR ITALIA, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

©2026 Kerr Corporation. All trademarks applicable to Kerr are property of Kerr Corporation. All other trademarks are property of their respective owners. All Rights Reserved.

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