UltraSeal XT™ plus/UltraSeal XT™ hydro/PrimaDry™

Pit and Fissure Sealant and Flowable Composite/Pit and Fissure Sealant and Flowable Composite/Drying and Priming Agent

EN

Product Description
UltraSeal XT plus sealant and flowable composite – UltraSeal XT plus sealant is a 50% filled, light-cured, radiopaque, methacrylate-based, thixotropic resin sealant. It is cured with PrimaDry drying and priming agent so aid in avoiding moisture contamination—which can cause microleakage and poor retention. It is a Type 1 Class 2 Group 1 polymer-based dental restorative flowable resin which is suitable for occlusal surfaces. It can be cured where the use of a flowable is indicated. It is compatible with adhesives and bonding agents. UltraSeal XT plus sealant has a radiopacity value of 1.44 mm of aluminum. Aluminum has a radiopacity equivalent to that of dentin. Thus, 1 mm of material (having a radiopacity equivalent to 1 mm of aluminum) has a radiopacity equivalent to that of dentin—and 1 mm of aluminum = 1 mm of dentin equivalent. The dimensions of filler particles are as follows:

Component	Particle size (µm)
Aluminum Oxide	≤ 0.05
Glass	1.5 x 0.25

The percentage by volume of filler particles is 33% (total fillers, 19.2 mL/total volume, 57.90 mL * 100)

UltraSeal XT plus sealant	VITA 3D-MASTER™ Shade Guide
OP White	OM
A2	A2
A3	A3
Clear – Unshaded base	

* Refer to registered usernames for Ultradent Products, Inc.
UltraSeal XT hydro-sealant – UltraSeal XT hydro sealant is a Class 2, 53% filled, light-cured, radiopaque, methacrylate-based, thixotropic resin sealant. It is partly hydrophobic and hydrophilic before it is cured, hydrophobic once cured, and has a self-adhesive quality. UltraSeal XT hydro sealant primarily provides the option of a visual verification for a material rinse—with the use of a black light upon placement and after it is cured etc. UltraSeal XT hydro sealant is offered in two shades: Opaque White and Natural. The UV light is not included in the UltraSeal XT hydro sealant kit.
PrimaDry drying agent – PrimaDry drying agent is a hypersonic drying and priming agent—used by dental professionals (in clinical settings)—to assure maximum drying and priming of etched fissure walls prior to resin placement. Utilizing ethanol alcohol, PrimaDry drying agent rapidly volatilizes the moisture content of fissures.
For all products described, carefully read and understand all instructions and SDS information prior to use.
Delivery Form(s)
- 1.2 ml syringe
- Delivery tip (UltraSeal XT plus tip and UltraSeal XT hydro-dispensing tip)
- Endless tip (PrimaDry drying agent): Black Micro™ Applicator Tip

Indications for Use/Intended Purpose
- **UltraSeal XT plus sealant** – UltraSeal XT plus sealant for prophylactic sealing of pits and fissures and can be used where the use of a flowable composite is indicated.
- **PrimaDry drying agent** – As an element of prophylactic pit and fissure sealant procedures, PrimaDry drying agent is used to virtually eliminate residual water and microleakage that can cause failure. PrimaDry drying agent is designed to be used prior to resin placement and is followed by drying with moisture-free and oil free air
- **UltraSeal XT hydro sealant** – Use UltraSeal XT hydro sealant for prophylactic sealing of pits and fissures.

Contraindications
- For patients or users with allergy concerns, refer to product allergen document available at www.ultradent.com. If allergic reaction is observed, rinse exposed area thoroughly with water and have the patient consult their physician.

Warnings and Precautions
- Resins can be sensitizing. Avoid repeated contact of uncured dental resin with skin. If contact occurs, thoroughly rinse skin. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis, or rash develops have patient consult their physician.
- Redirect overhead light, or use UV protective light shield, to prevent premature polymerization.
- Re-cap immediately following use to avoid polymerization in syringe.
- Tips are single use. Avoid cross-contamination between patients.

Stepwise Instructions
UltraSeal XT plus and UltraSeal Xtreme
Preparation
1) Etch recommended tip to Ultra-Etch™ 35% etchant.
2) Apply Universal Brush Tip to UltraSeal XT plus sealant syringe. If desired, the Inspiral Brush Tip bristle length can be adjusted by gripping and moving inward or outward.
3) Attach recommended tip to PrimaDry drying agent syringe.
4) Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
5) To prepare tooth prior to sealant, use one of the following methods:
- an air abrasion unit with appropriate particles
- a thin diamond finishing/prep bur
- air polishing prophylaxis system
- Prophyjet™ system
- Prophy paste
- Flour of pumice

NOTE: Use high-speed suction when using UltraSeal XT plus sealant for restoratives.
5) Apply PrimaDry drying agent to manufacturer's instructions.
6) Attach Inspiral Brush Tip to UltraSeal XT plus sealant syringe.
7) Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
8) Apply sealant UltraSeal XT plus sealant 0.5–1.5 mm thick and light cure 10 seconds.

UltraSeal XT hydro sealant – Pit and Fissure Sealant
Preparation
1) Attach recommended tip to Ultra-Etch 35% etchant.
2) Apply Universal Brush Tip to UltraSeal XT hydro-dispensing syringe. If desired, the Inspiral Brush Tip bristle length can be adjusted by grasping and moving inward or outward.
3) Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
4) To prepare tooth prior to sealant, use one of the following methods:
- an air abrasion unit with appropriate particles
- a thin diamond finishing/prep bur
- air polishing prophylaxis system
- Prophyjet™ system
- Prophy paste
- Flour of pumice

NOTE: Use high-speed suction when using UltraSeal XT plus sealant for restoratives.
5) Apply 35% phosphoric acid etchant (Ultra-Etch etchant) to fissures for 30 seconds on uncoat enamel or 15 seconds on cut enamel.
NOTE: A fissure tear and some air abrasion particles are the only preparations for removing the aprismatic enamel. All other preparations require the 30-second etch.
6) Thoroughly rinse.
Lightly air dry the tooth to remove visible moisture. Do not desiccate the tooth.
Placering
NOTE: Do not use adhesive bonding agent when using UltraSeal XT hydro sealant.
1) Express a small drop of UltraSeal XT hydro sealant into the grooves. Work the resin down into the grooves with the brush tip.
2) Light-cure each fissure using light-cure standard power for 10 seconds on Xtra Power mode.
3) Use an explorer to ensure pits and fissures are sealed and adjust occlusion if necessary.

UltraSeal XT hydro sealant – Fissure Sealant
NOTE: Do not use PrimaDry drying agent when using UltraSeal XT plus sealant for restoratives.
1) Apply dentin bonding agent according to manufacturer's instructions.
2) Attach Inspiral Brush Tip to UltraSeal XT plus sealant syringe.
3) Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
4) Apply sealant UltraSeal XT plus sealant 0.5–1.5 mm thick and light cure 10 seconds.

UltraSeal XT hydro sealant – Pit and Fissure Sealant
Preparation
1) Attach recommended tip to Ultra-Etch 35% etchant.
2) Apply Universal Brush Tip to UltraSeal XT hydro-dispensing syringe. If desired, the Inspiral Brush Tip bristle length can be adjusted by grasping and moving inward or outward.
3) Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
4) To prepare tooth prior to sealant, use one of the following methods:
- an air abrasion unit with appropriate particles
- a thin diamond finishing/prep bur
- air polishing prophylaxis system
- Prophyjet™ system
- Prophy paste
- Flour of pumice

NOTE: Use high-speed suction when using UltraSeal XT plus sealant for restoratives.
5) Apply 35% phosphoric acid etchant (Ultra-Etch etchant) to fissures for 30 seconds on uncoat enamel or 15 seconds on cut enamel.
NOTE: A fissure tear and some air abrasion particles are the only preparations for removing the aprismatic enamel. All other preparations require the 30-second etch.
6) Thoroughly rinse.
Lightly air dry the tooth to remove visible moisture. Do not desiccate the tooth.
Placering
NOTE: Do not use adhesive bonding agent when using UltraSeal XT hydro sealant.
1) Express a small drop of UltraSeal XT hydro sealant into the grooves. Work the resin down into the grooves with the brush tip.
2) Light-cure each fissure using light-cure standard power for 10 seconds on Xtra Power mode.
3) Use an explorer to ensure pits and fissures are sealed and adjust occlusion if necessary.

UltraSeal XT hydro sealant – Fissure Sealant
NOTE: Do not use PrimaDry drying agent when using UltraSeal XT plus sealant for restoratives.
1) Apply dentin bonding agent according to manufacturer's instructions.
2) Attach Inspiral Brush Tip to UltraSeal XT plus sealant syringe.
3) Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
4) Apply sealant UltraSeal XT plus sealant 0.5–1.5 mm thick and light cure 10 seconds.

UltraSeal XT hydro-sealant – Fissure Sealant
Preparation
1) Etch recommended tip to Ultra-Etch 35% etchant.
2) Apply Universal Brush Tip to UltraSeal XT hydro-dispensing syringe. If desired, the Inspiral Brush Tip bristle length can be adjusted by grasping and moving inward or outward.
3) Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
4) To prepare tooth prior to sealant, use one of the following methods:
- an air abrasion unit with appropriate particles
- a thin diamond finishing/prep bur
- air polishing prophylaxis system
- Prophyjet™ system
- Prophy paste
- Flour of pumice

NOTE: Use high-speed suction when using UltraSeal XT plus sealant for restoratives.
5) Apply 35% phosphoric acid etchant (Ultra-Etch etchant) to fissures for 30 seconds on uncoat enamel or 15 seconds on cut enamel.
NOTE: A fissure tear and some air abrasion particles are the only preparations for removing the aprismatic enamel. All other preparations require the 30-second etch.
6) Thoroughly rinse.
Lightly air dry the tooth to remove visible moisture. Do not desiccate the tooth.
Placering
NOTE: Do not use adhesive bonding agent when using UltraSeal XT hydro sealant.
1) Express a small drop of UltraSeal XT hydro sealant into the grooves. Work the resin down into the grooves with the brush tip.
2) Light-cure each fissure using light-cure standard power for 10 seconds on Xtra Power mode.
3) Use an explorer to ensure pits and fissures are sealed and adjust occlusion if necessary.

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