VDW® | Endo Easy Efficient

RECIPROC® blue

EN	Instructions For Use	2
FR	Instructions d'utilisation	9
IT	Istruzioni per l'uso	16
DE	Gebrauchsanweisung	23
ES	Instrucciones de uso	31
PT	Instruções de Utilização	38
NL	Gebruiksinstructies	45
SV	Användarinstruktioner	52
DA	Brugsanvisning	59
NO	Brukerinstruksjoner	66
FI	Käyttöohjeet	73
PL	Instrukcja obsługi	80
LT	Naudojimo instrukcija	87
LV	Lietošanas norādījumi	94
ET	Kasutusjuhised	101
CS	Návod k použití	108
SK	Návod na použitie	115
HU	Használati utasítás	122
RO	Instrucțiuni de utilizare	129
SR	Uputstva za upotrebu	136
HR	Upute za uporabu	144
BG	Инструкции за употреба	151
BS	Uputstva za upotrebu	159
MK	Упатство за употреба	166
EL	Οδηγίες χρήσης	174
UK	Інструкції із застосування	181
TR	Kullanım Kılavuzu	189
KO	사용 지침	196
ID	Petunjuk Penggunaan	204
PT-BR	Instruções de uso	211

VDW®
Endo
Easy
Efficient

RECIPROC® blue

Instructions
For Use

www.vdw-dental.com/ifu

vdw-dental.com

EN

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
SINGLE USE – REPROCESSING PROHIBITED
STERILE – STERILIZED BY RADIATION

INSTRUCTIONS FOR USE

SINGLE USE STERILE ENDODONTIC INSTRUMENTS

REF. V040252, V040052

0. INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation and retreatment (shaping and debridement of the root canal).
- Intended User:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as Endo specialists (Endodontist) and Dental Assistants.
- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).

1. CONTRAINDICATIONS

None.

2. CONTENT, COMPOSITION AND COMPATIBLE DEVICES

RECIPROC® blue instruments are made of three main components: the working part made of a nickel-titanium alloy with a heat treatment, a colored silicone stopper and a plated brass shank with an ABS plastic ring.

- The files are manufactured with a process that results in a file that has a colored appearance. The files may appear slightly curved due to this process. This is not a manufacturing defect, and it is not necessary to straighten the file prior to use. Once inside the canal, the file will follow the root canal anatomy.
- Use in a reciprocing motion with a designated drive system having the approved RECIPROC® settings.
- Lubricants such as NaOCl, EDTA, ProLube®, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicon stopper to the correct working length is the appropriate method of working length determination.

3. WARNINGS

- Strictly follow this Instruction for Use to minimize following risks to the devices, the patient and/or the user:
 - Breakage of instrument.

 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Toxic or allergic reactions caused by processing residues.
- Similar to all machine-driven root canal instruments, RECIPROC® blue instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

4. PRECAUTIONS

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible or/and NiTi files (this will help avoid canal transportation).
 - Visually inspect the working part for all the defects listed in the former paragraph during use (i.e after each wave).
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

5. ADVERSE REACTIONS / UNDESIRABLE SIDE EFFECTS

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

6. STORAGE CONDITIONS

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

7. STEP BY STEP INSTRUCTIONS

Shaping procedure with RECIPROC® blue

Select the correct RECIPROC® blue instrument:

In most cases, size R25 will be suitable for the root canal treatment.

1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
 - If the canal is partially or completely invisible on the radiograph: The canal is considered narrow = use size R25.
 - If the canal is completely visible on the radiograph: Take an ISO size 30 hand instrument and insert it passively into the canal:
 - If it reaches working length, the canal is considered large = use size R50.
 - If an ISO size 30 hand instrument does not go passively to working length, take an ISO size 20 hand instrument. If file goes passively to length, the canal is considered medium = use size R40.
 - If an ISO size 20 hand instrument does not go passively to length, the canal is considered narrow = use size R25.

Note: Passively means that the instrument goes directly to working length without any watch-winding and can be moved up-and-down without shaping action. The amplitude of the up-and-down movement should be 4-5 mm.

2) Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.
3) Irrigate abundantly and frequently after each instrument used (according to good dental practices).
4) Starting preparation depending on the instrument size:
 - Size R25: Prior to preparation with an R25, estimate the working length with the help of a pre-operative radiograph. Set the stopper at 2/3 of that length.
 - Sizes R40 and R50: Before selecting the appropriate RECIPROC® blue instrument, check the working length electronically by using a C-PILOT® or K-File. Set the stopper at that length.
5) With irrigant in the canal, introduce the RECIPROC® blue instrument into the canal. Press the ON/ OFF button or foot pedal of the motor when the instrument is at the root canal orifice.
6) Move the instrument in a slow in-and-out pecking motion until approx. 2/3 of the working length has been reached (indicated by the silicone stopper). One in-and-out movement = one peck. Clean the instrument after each pecking motion. The amplitude of the in-and-out movements should not exceed 3 mm. Only very light pressure should be applied. The instrument will advance easily in the canal.
7) If necessary, the instrument can be used in a brushing motion to enlarge the root canal entrance.
8) After each pecking motion, make sure the canal is free with an ISO size 010 C-PILOT® file.
9) Reconfirm the working length, before continuing shaping until full working length has been reached.
10) Gauge the foramen with a hand instrument which is one ISO size larger than the RECIPROC® blue instrument that reached full working length.

EN

- If the instrument is snug at length, preparation is finished and the canal is ready to be obturated.
- If the instrument is loose at length, use next size of RECIPROC® blue instrument, and, when necessary proceed with the next file size, gauging after each instrument respectively with a hand files of the same size. Use the same shaping technique up to working length.

11) Use dedicated paper points to dry and dedicated gutta-percha points to obturate the root canals.

Retreatment procedure with RECIPROC® blue

1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
2) Remove the bulk of the gutta-percha in the coronal third of the canal with an appropriate instrument e.g. Gates Enlarger or with the help of an ultrasonic device.
3) If necessary, soften gutta-percha with a solvent (e.g. eucalyptus oil).
4) Use size R25 as described until working length has been reached. If resistance is encountered, do not apply pressure.
5) Remove the instrument from the canal, apply a little solvent and repeat the process. The carrier may be removed in one piece during the use of the RECIPROC® blue instrument; otherwise, it will be removed in small pieces with the gutta-percha.
6) Use a brushing motion against lateral walls to remove residual obturation material.
7) After reaching working length with size R25, use size R40 or R50 for an increased apical enlargement if mandatory.

8. HYGIENE, DISINFECTION, CLEANING AND STERILIZATION

- RECIPROC® blue instruments are single use devices and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
- Sterility cannot be guaranteed if packaging is open, damaged or wet.
- Disinfection, cleaning and sterilization is not applicable.
- Products shall be disposed according to local regulations for the safe disposal of sharp and contaminated devices.

9. ADDITIONAL INFORMATION

- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
- To get a free printed copy of IFU with appropriate language please see section “Request an IFU” on website www.vdw-dental.com/ifu.

10. LIST OF APPLICABLE SYMBOLS

Symbols	Description
	**Catalogue number:** Indicates the manufacturer’s catalogue number (SKU). This symbol is accompanied by the manufacturer’s catalogue number adjacent to the symbol.
	**Use-by-date:** Indicates the date after which the medical device is not to be used. This symbol is accompanied by a date to indicate that the medical device should not be used after the end of the year, month or day shown. The date is expressed in YYYY-MM-DD (or YYYY-MM) format.
	**Manufacturer:** Indicates the medical device legal manufacturer [for a Medical Device]. Indicates the procedure pack producer [for a Procedure Pack].
	**Consult instructions for use or consult electronic instructions for use:** Indicates the need for the user to consult the instructions for use. This symbol could be accompanied by the electronic IFU website address.

Symbols	Description
STERILE R	**Sterilized using radiation:** Indicates a medical device that has been sterilized using irradiation.
LOT	**Batch code:** Indicates the manufacturer's batch code. This symbol is accompanied by the manufacturer's batch code adjacent to the symbol.
	Assortment: Identifies the package includes an assortment of types or sizes.
NiTi	**NiTi:** Identifies a material or substance contained in a product or the material from which the product is made. "XXX" is replaced by material letter code (e.g., "SST" for Stainless Steel) or a word (e.g., "Latex"). NiTi: Nitinol / Nickel titanium
	Do not use if package is damaged: Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.
X	**Packaging unit:** Indicates the number of pieces in the package. "X" is replaced by the number of pieces in the package.
	Do not re-use: Indicates a medical device that is intended for one single use only or for use on a single patient during a single procedure.
	RECIPROC: Indicates reciprocating movement.
CE 0123	**CE marking:** CE Marking of Conformity according to EU MDD/MDR requirements. Note: When a notified body is involved, the CE marking is accompanied by the identification number (four-digit number) of the notified body.
MD	**Medical device:** Indicates the item is a medical device.
	Caution: Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
	Keep away from sunlight: Indicates a medical device that needs protection from light sources. This symbol can also mean "Keep away from heat".
	Single sterile barrier system: Indicates a single sterile barrier system.
	Date of manufacture: Indicates the date when the medical device was manufactured. This symbol is accompanied by a date to indicate the date of manufacture. The date is expressed in YYYY-MM-DD (or YYYY-MM) format.
CC	**Country of manufacture:** Indicates the country of manufacture of the products. "CC" is replaced by the two or three letter country code as defined in ISO 3166-1.

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VDW GmbH