Root Filler

EN	Instructions For Use	2
FR	Instructions d'utilisation	13
IT	Istruzioni per l'uso	24
DE	Gebrauchsanweisung	35
ES	Instrucciones de uso	47
PT	Instruções de utilização	58
NL	Gebruiksinstructies	70
SV	Användarinstruktioner	82
DA	Brugsanvisning	93
NO	Brukerinstruksjoner	104
LT	Naudojimo instrukcija	115
LV	Lietošanas norādījumi	126
ET	Kasutusjuhised	137
CS	Návod k použití	148
SK	Návod na použitie	159
HU	Használati utasítás	170
RO	Instrucțiuni de utilizare	181
SR	Uputstva za upotrebu	192
BS	Uputstva za upotrebu	203
CNR	Uputstva za upotrebu	214
SQ	Udhëzimet për përdorim	225
EL	Οδηγίες χρήσης	236
TR	Kullanım Kılavuzu	248
KO	사용 지침	259

Root Filler

Endodontic instruments

Instructions For Use

A brand of Dentsply Sirona

dentsplysirona.com/ifu

www.dentsplysirona.com

EN

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
STERILE REFERENCES – STERILIZED BY GAMMA IRRADIATION
NON-STERILE REFERENCES – MANDATORY PROCESSING BEFORE FIRST USE
REUSABLE – LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

VDW®STERILE ROOT FILLER -L-	(STERILE REF. V040393)
ANTAEOS® ROOT FILLER -L-	(NON-STERILE REF. V010093)
BEUTELROCK® ROOT FILLER -L-	(NON-STERILE REF. V020093)
ZIPPERER® ROOT FILLER -L-	(NON-STERILE REF. V030093)
KENDO® ROOT FILLER	(NON-STERILE REF. V201709)

1 INDICATIONS FOR USE

- Indication for use:
 Paste carriers are indicated for treatment of endodontic diseases.
- Intended purpose:
 Rotary endodontic paste carrier instruments used to introduce sealer, cement or paste in the root canal.
- Intended patient population:
 Patients who need root canal treatment.
- Intended user:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

2 CONTRAINDICATIONS

None.

3 COMPOSITION

Root Filler instruments are made of a working part made of stainless steel, and a nickel-plated brass shank with a colored ring.

4 COMPATIBLE DEVICES

- Root Filler instruments shall be connected to a contra-angle and used in combination with endodontic motors at constant rotation speed within the 300-600 rotations/min range, clockwise direction.
- Root Filler instruments can be used in combination with sealer, cement or paste.
 Safety and performance of Paste Carriers have been tested with Zinc Oxyphosphate Cement Dentsply Detrey cement REF. 611.13.206/611.13.003, AH Plus™ Jet sealer REF. 60620118 and AH Temp™ Dentsply Detrey paste REF. 60620123.

5 WARNINGS

- Strictly follow the Instructions for Use to minimize the following risks to the device, the patient and/or the user:
 - Serious health consequences (infection, toxic or allergic reactions).
 - Tooth extraction.
 - Swallowing or inhalation of the instrument.
 - Treatment delay.
- Similar to all machine-driven root canal instruments, Root Filler instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation.
- The use of inappropriate rotation speed and pressure would speed up the wear of instruments and increase the risk of breakage.

6 PRECAUTIONS

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, always wear protective gloves, glasses and a mask when handling contaminated instruments.
- For sterile references (V04039): Do not use after expiration date.
- Visually inspect the integrity of the packaging (sterile barrier when applicable) before use. Do not use the instruments if the packaging is damaged.
- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Processing step is mandatory before use according to section [10] for non-sterile devices mentioned in this IFU.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Always use minimal apical pressure. Never force instruments down the canal.
- Limitation and restrictions on processing and on re-use:
 - Instruments can be re-used several times – with due care and if they are not damaged and contaminated. Cleaned and undamaged instruments can be used up to eight times depending on the degree of wear. This maximum number of uses has been validated with the standard methods (automated cleaning and disinfection, fractionated vacuum method for steam sterilization). Yet, the maximum number of processing cycles will not

always be reached. Certain applications may cause the instruments to prematurely reach the end of their useful life. Defective instruments should be immediately discarded (according to section [10.2.4]).

- Each re-use or application of non-validated methods is the sole responsibility of the user. All liability is disclaimed for failure to follow these instructions or use of non-validated methods for the re-use of instruments.
- Please always ensure that sterile packaging/wrapping is undamaged. Do not use the instruments if the packaging is damaged.

7 ADVERSE REACTIONS / UNDESIRABLE SIDE EFFECTS

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

8 STORAGE CONDITIONS

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

9 STEP-BY-STEP INSTRUCTIONS

1. Apply the material onto the tip of a Root Filler.
2. Wet the canal walls with the material through a pumping or simultaneously rotating movement (300-600 rotations/min, clockwise direction).
3. Advance the Root Filler slowly to the apex running at low speed (300-600 rotations/min). Avoid the formation of air bubbles in the material and overfilling of the canal.
4. Withdraw Root Filler very slowly still running at low speed (300-600 rotations/min).

10 HYGIENE, DISINFECTION, CLEANING AND STERILIZATION

10.1 General basics

- All instruments must be cleaned, disinfected and sterilized prior to each use; this applies the first time instruments supplied in a non-sterile condition are used, and to instruments delivered in a sterile condition that are intended for repeated use. Thorough cleaning and disinfection are essential prerequisites for effective sterilization.
- As part of your responsibility for the sterility of instruments, always make sure that only validated methods for cleaning/disinfection and sterilization are used, that devices (washer-disinfector, thermal disinfector or sterilizer) are regularly serviced and inspected, and that the validated parameters are maintained during each cycle.
- In addition, always observe all applicable national legal regulations and regulations on hygiene relating to your practice or the hospital. This applies in particular to the guidelines regarding prion inactivation.
- Recommended classification according to RKI/BfArM/KRINKO directive (Germany only, intended use): Critical A.

10.2 Step-by-step procedure

10.2.1 Pre-treatment at the place of use

The pre-treatment process should be performed on used instruments in every case.

For new non-sterile instruments, refer straight to section [10.2.3].

Contaminants (particularly pulp and dentine remnants) must be removed immediately after the instrument has been used on a patient (within maximum 2 hours). All further steps in the preparation process must be performed on the same day.

The following procedures must be used to ensure that no contamination can dry on the instruments, and to make subsequent preparation more effective:

Step 1

A prepared interim stand (REF. V040495) with a new foam disc must be used for each patient. The interim stand must be filled at least two thirds of the way with disinfecting agent.

Step 2

Place in the interim stand prior to pre- disinfection/cleaning and for transport (minimum storage time according to the disinfecting agent manufacturer's instructions for use: max. two hours).

Please note that the disinfecting agent used during pre-treatment is for personal protection only and is not a substitute for the disinfection stage required after cleaning.

All further steps in the preparation process must be performed on the same day.

Warnings:

- Under no circumstances may instruments that have already come into contact with disinfecting agent be used to treat a patient again.

10.2.2 Preparation before cleaning

Step 3

Clean to remove contamination under flowing water for at least 3 x 1 minute; to remove contamination manually, use a soft, clean brush or soft, clean cloth that is only used for this purpose; never use metal brushes or wire wool.

Step 4

Check that no visible contamination or remnants remain and repeat the pre-cleaning process if necessary.

10.2.3 Cleaning and disinfection

We recommend an automated procedure to clean and disinfect the instruments. A manual method, including the use of an ultrasonic bath, should only be used if it is not possible to use an automated method, as it is less effective and demonstrates a lower reproducibility.

Manual cleaning and disinfection is less effective in direct comparison to the automated method. However, it is effective according to the requirements for a processed instrument. All methods are validated and therefore they are efficient and safe for the processing of VDW instruments.

For processing we do recommend the use of the LavEndo® box (REF. V040479) with mini step module (REF. V040403).

Warnings:

- In the case of chemical disinfection, there is a risk of disinfecting agent residues remaining on the instruments.

Step 5: Automated cleaning/disinfection (washer-disinfector) – Recommended procedure

The following must be taken into account when selecting a washer-disinfector:

- The effectiveness of the washer-disinfector must have been verified (DGHM approval, FDA clearance or CE mark according to EN ISO 15883).
- Where possible, a tested thermal disinfection program must be used (A0 value ≥ 3000 or at least five minutes at 90°C, or for older equipment at least 10 min. at 93°C).
- The program used must be suitable for the instruments and include the prescribed rinsing cycles.
- Only sterile or low-germ (< 10 CFU/ml) and low-endotoxin (< 0.25 EU/ml) water (ideally highly purified water HPW) must be used for subsequent rinsing.
- The washer-disinfector must be regularly maintained and inspected.

1) Sort the instruments into the step modules (REF. V040403).
2) Place the endo module in the blue lower section of the LavEndo® box (REF. V040479) and close it (click into place).
 NOTE: Preparation in the socket module is not permitted.
3) Insert the LavEndo® box horizontally into the washer-disinfector.
4) Start the program.
5) After the program has finished, remove the LavEndo® box from the washer-disinfector.
6) Check and package the instruments as soon as possible after removing them (see sections [10.2.4] and [10.2.5]), after leaving them to dry further in a clean place if necessary.

An independent, accredited, recognized test laboratory has demonstrated the intrinsic suitability of the instruments for effective automated cleaning and disinfection using the G 7836 CD washer-disinfector (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the Neodisher Medizym cleaning agent (Dr. Weigert, Hamburg). The laboratory used program D-V-MEDIZYM (based on the program DES-VAR-TD (Miele) under worst-case conditions) according to the procedure described above to demonstrate this effectiveness. The Cidex OPA disinfecting agent and the Cidezyme cleaning agent (both Johnson & Johnson GmbH, Norderstedt) were used for pre-treatment. Cleaning and disinfection validation were performed under worst-case conditions (low temperature, low concentration of agent, short soaking time and no drying).

Step 5: Or Manual cleaning and disinfection – If automated procedure is not possible

<u>Manual Cleaning:</u>

1) Sort the instruments into the step modules (REF. V040403).
2) Place the endo module in the blue lower section of the LavEndo® box (REF. V040479) and close it (click into place).
 NOTE: Preparation in the socket module is not permitted.
3) Insert the LavEndo® box into the cleaning bath for the prescribed contact time, ensuring that the instruments are sufficiently covered.

4) Then remove the LavEndo® box from the cleaning bath and rinse thoroughly with water for at least 3 x 1 min.

Manual disinfection:

5) Insert the LavEndo® box horizontally into the disinfection bath for the prescribed contact time, ensuring that the instruments are sufficiently covered.
6) Then remove the LavEndo® box from the disinfection bath and rinse thoroughly with water at least 5 x 1 min.
7) Dry the LavEndo® box by blowing them down fully with oil-free, filtered compressed air (or medical compressed air from a can) and then leaving them to dry further in a clean place.
8) Check (see section [10.2.4]) package the instruments as soon as possible (see section [10.2.5]).

An independent, accredited, recognized test laboratory has demonstrated the intrinsic suitability of the instruments for effective manual cleaning and disinfection using the cleaning agent Cidezyme/Enzol and disinfecting agent Cidex OPA (Johnson & Johnson GmbH, Norderstedt (Germany)). The laboratory used the procedure described above to demonstrate this. The Cidex OPA disinfecting agent and the Cidezyme cleaning agent (both Johnson & Johnson GmbH, Norderstedt) were used for pre-treatment.

10.2.4 Inspection and maintenance

Open the LavEndo® boxes and remove the step modules. Check all instruments, modules and LavEndo® boxes after cleaning/disinfection. Defective instruments, boxes and modules should be immediately discarded.

- These defects include:
 - Plastic deformation (e.g. caused by an excessively high temperature during sterilization).
 - Breakage.
 - Loss of color coding or marking.
 - Bent instrument.
 - Corrosion.
- Numerical restrictions on re-use are listed under section [6], “Limitation and restrictions on processing and on re-use”.
- Instruments that are still contaminated must be cleaned and disinfected again.

Warnings:

- Instrument lubricants must not be used.

10.2.5 Packaging

Place the step module (REF. V040403) in the lower section of the black sterilization mini box (REF. V040471) and close it with the matching cover. Then package the sterilization trays into disposable sterilization pouches (disposable packaging). Single package using hospital common sterilization packaging shall be used (e.g. single Steriking paper/film packagings from Wipak Medical was used for the validation).

The sterilization pouches shall meet the following requirements:

- Compliance with DIN EN 11607/ANSI AAMI ISO 11607.

- Suitable for steam sterilization (withstands temperatures of up to 142°C (288°F) or more, sufficient vapor permeability).

Warnings:

- Sterilization in the sterilization trays without additional packaging is not permitted. The autoclave paper in the boxes is for added safety only.

10.2.6 Sterilization

Only use the sterilization methods listed below; other sterilization methods are not permitted.

Steam sterilization

- Fractionated vacuum/pre-vacuum method (at least three vacuum cycles) or gravity displacement method[1] with sufficient product drying.
- Steam sterilizer in accordance with DIN EN 13060 or DIN EN 285, ANSI AAMI ST79.
- Validated in accordance with DIN EN ISO 17665 (valid IQ and OQ plus product-specific performance qualification (PQ)).
- The maximum sterilization temperature of 138°C (280°F) must not be exceeded; the maximum sterilization temperature includes a tolerance according to DIN EN ISO 17665.
- The drying time that is actually required depends directly on parameters that are the sole responsibility of the user (loading configuration, how many items are loaded and how closely together they are loaded, condition of the sterilizer, etc.) and must therefore be established by the user. However, the drying time must never be less than 20 minutes.

Sterilization procedure	Sterilization temperature	Minimum sterilization time Exposure time at sterilization temperature
Fractionated vacuum/pre-vacuum method	134°C (273°F)	3 minutes[2]
	121°C (250°F)	20 minutes
Gravity method[1 & 3]	134°C (273°F)	15 minutes
	121°C (250°F)	60 minutes

[1] The less effective gravity method should only be used if the fractionated vacuum method is not available. The gravity method is less effective in direct comparison to the fractionated vacuum method. However, it is effective according to the requirements for a processed instrument. All methods are validated and therefore they are efficient and safe for the processing of VDW instruments.

[2] Or 18 min. (prion inactivation).

[3] Gravity method is not applicable for processing within the European Union.

An independent, accredited, recognized test laboratory demonstrated the instruments' intrinsic suitability for effective steam sterilization using the HST 6x6x6 steam sterilizer (Zirbus Technology GmbH, Bad Grund) together with the fractionated vacuum method and the gravity method.

10.3 Storage and transport

After sterilization, devices must be stored in the sterilization packaging and kept dry and dust-free. In case of damage to the packaging during storage or transport, the processing shall be repeated. Check the instructions for use given by the pouch manufacturer to determine the shelf life of the sterile packaging.

10.4 Disposal

Instruments shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

10.5 Important Information on material resistance

- When selecting cleaning and disinfecting agents, make sure that they do not contain any of the following substances:
 - Phenol.
 - Strong acids (ph < 6) or strong alkalis (ph > 8); neutral enzymatic cleaning agent recommended.
 - Aldehydes.
 - Anti-corrosive substances (especially di- or triethanolamine).
 - Oxidants (hydrogen peroxide, sodium hypochlorite over 5% strength).
 - Solvents.
 - Oils.

Warnings:

- Never clean the instruments, boxes, modules or the interim stand with metal brushes or wire wool.
- Never subject any instruments, boxes, modules or the interim stand to temperatures above 142°C (288°F). It is particularly important to ensure that the products to be sterilized are not stored too close to the walls or floor of the steam sterilizer (risk of excessive temperature and deformation).
- The blue foam insert for the interim stand must only be used once and used blue foam inserts must not be either cleaned/disinfected or sterilized.

The following must be taken into account when selecting cleaning and disinfecting agents:

- They must be suitable for cleaning and disinfecting instruments made from metal and plastic.
- The disinfecting agent must be aldehyde-free (Cidex OPA is permitted due to its special recipe).
- It must be compatible with the instruments.
- A disinfecting agent with verified effectiveness (VAH/DGHM approval, FDA clearance or CE mark) must be used and this must be compatible with the cleaning agent used.
- If a thermal disinfection process is not used, a suitable disinfecting agent with verified effectiveness (VAH/DGHM approval, FDA clearance or CE mark) must also be used and this must be compatible with the cleaning agent used.
- Neutralization must not be necessary (cleaning agent).
- The cleaning agent, if applicable, must be suitable for ultrasonic cleaning (no foaming).
- Combined cleaning agents/disinfecting agents must not be used.

The concentrations, temperatures and contact times specified by the manufacturer of the cleaning agent and disinfecting agent as well as the minimum specifications for subsequent rinsing must be strictly adhered to. Rinse aids must not be used.

Only use freshly prepared solutions and low-germ (< 10 CFU/ml) water; tap water that is particularly hard (≥ 14°dH) is not suitable for this (risk of lime residue).

11 ADDITIONAL INFORMATION

- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
- Disclaimer:
 Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for processing. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the instructions for use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.
- See webpage dentsplysirona.com/ifu to:
 - Get a free print copy of the IFU.
 - Get an explanation of non-harmonized symbols for IFUs and labels.

CE* 2797

EU REP

Authorised representative

DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

Manufacturer

Maillefer Instruments Holding Sàrl
Chemin du Verger 3
1338 Ballaigues
Switzerland

dentsplysirona.com/ifu

EU Importer

VDW GmbH
Bayerwaldstr. 15
81737 München
Germany

*CE is not applicable to the interim stand and all sterilization boxes

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