VDW® | Endo Easy Efficient

VDW.ROTATE®

EN	Instructions For Use	2
FR	Instructions d'utilisation	15
IT	Istruzioni per l'uso	28
DE	Gebrauchsanweisung	41
ES	Instrucciones de uso	54
PT	Instruções de Utilização	67
NL	Gebruiksinstructies	80
SV	Användarinstruktioner	93
DA	Brugsanvisning	106
NO	Brukerinstruksjoner	119
FI	Käyttöohjeet	132
PL	Instrukcja obsługi	145
LT	Naudojimo instrukcija	158
LV	Lietošanas norādījumi	170
ET	Kasutusjuhised	183
CS	Návod k použití	196
SK	Návod na použitie	209
HU	Használati utasítás	222
RO	Instrucțiuni de utilizare	235
SR	Uputstva za upotrebu	248
HR	Upute za uporabu	261
BG	Инструкции за употреба	274
BS	Uputstva za upotrebu	289
MK	Упатство за употреба	303
EL	Οδηγίες χρήσης	318
UK	Інструкції із застосування	332
TR	Kullanım Kılavuzu	347
KO	사용 지침	360
ID	Petunjuk Penggunaan	373
PT-BR	Instruções de uso	386

Endo
Easy
Efficient

VDW.ROTATE®

Instructions
For Use

www.vdw-dental.com/ifu

vdw-dental.com

EN

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
STERILE – STERILIZED BY RADIATION
LIABLE FOR REPROCESSING

INSTRUCTIONS FOR USE

REUSABLE STERILE ENDODONTIC INSTRUMENTS

REF. MSTVROT, MDDVROT

0. INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation and retreatment (glide path, shaping and debridement of the root canal).
- Intended User:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as Endo specialists (Endodontist) and Dental Assistants.
- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).

1. CONTRAINDICATIONS

None.

2. CONTENT, COMPOSITION AND COMPATIBLE DEVICES

VDW.ROTATE® instruments are made of three main components: the working part made of a nickel-titanium alloy with a heat treatment, a colored silicone stopper and a plated brass shank with colored rings.

The files are manufactured with a process that results in a file that has a colored appearance.

Due to this proprietary processing, the files may appear slightly curved. This is not a manufacturing defect. While the file can be easily straightened using only your fingers, it is not necessary to straighten the file prior to use. Once inside the canal, the file will follow the root canal anatomy.

- Use with endodontic motors in a constant rotation at a speed of:
 - 400 rpm for VDW.ROTATE® Retreatment ISO 025 005.
 - 300 - 400 rpm for other references.
- Torque control devices ensure optimal usage.
- Set torque depending on the size and according to indications in the table below for VDW.ROTATE® files:

File Size	Torque [Ncm]
VDW.ROTATE® ISO 015 004	1.3
VDW.ROTATE® ISO 020 005	2.1

EN

File Size	Torque [Ncm]
VDW.ROTATE® ISO 025 004, ISO 025 006, ISO 035 004, ISO 040 004, ISO 050 004, ISO 060 004, ISO 035 006, ISO 040 006	2.3
VDW.ROTATE® ISO 030 004, ISO 030 006	2.0
VDW.ROTATE® Retreatment ISO 025 005	3.5

- Lubricants such as NaOCl, EDTA, ProLube®, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicon stopper to the correct working length is the appropriate method of working length determination.

3. WARNINGS

- Strictly follow this Instruction for Use and the General Processing Instructions For Endodontic Products (see section [8.]) to minimize following risks to the device, the patient and/or the user:
 - Breakage of instrument.
 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Toxic or allergic reactions caused by processing residues.
- Similar to all machine-driven root canal instruments, VDW.ROTATE® instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

4. PRECAUTIONS

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.

- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible or/and NiTi files (this will help avoid canal transportation).
 - Visually inspect the working part for all the defects listed in the former paragraph during use (i.e after each wave).
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.
- Limitation and restrictions on processing and on re-use:
 - Do not re-use silicone stops. Remove and discard silicone stops after each use (during processing preparation before cleaning, see section [8.2.2.]), and replace with a new VDW stopper either before or after automated cleaning/disinfection.
 - Maximum number of processing cycles: 8 times (validated with the standard methods (automated cleaning and disinfection, fractionated vacuum method for steam sterilization).
 - Cleaned and undamaged instruments can be used up to eight times, defective instruments should be immediately discarded (according to section [8.2.4.]).
 - Instruments (only multiple use instruments) can be re-used several times – with due care and if they are not damaged and contaminated. Each re-use or application of non-validated methods is the sole responsibility of the user.
 - Certain applications may cause the instruments to prematurely reach the end of their useful life. The maximum number of processing cycles will not always be reached.
 - All liability is disclaimed for failure to follow these instructions or use of non-validated methods for the re-use of instruments.
 - Please always ensure that sterile packaging/wrapping is undamaged. Do not use the instruments if the packaging is damaged.

5. ADVERSE REACTIONS / UNDESIRABLE SIDE EFFECTS

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

6. STORAGE CONDITIONS

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

7. STEP BY STEP INSTRUCTIONS

Shaping procedure with VDW.ROTATE®

Although the VDW.ROTATE® instruments are versatile, the 2 following basic techniques are recommended:

- Pecking motion: advance of the instrument in 2-3 gentle amplitudes approximately 1-3 mm in-and-out of the canal.
- Brushing motion using the following technique:
 - Insert the rotating instrument into the root canal without touching the walls of the canal.
 - Exerting gentle pressure, allow the instrument to touch the canal wall.
 - Make small, stroking/brushing movements (over 1-3 mm) in a coronal direction, without taking the instrument out of the canal.

- Allow the instrument to move apically for 1-3 mm and then repeat the movement described under c. Gradually advance step-by-step towards the apex with up and down movements.

1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
2) Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.
3) Irrigate abundantly and frequently after each instrument used (according to good dental practices).
4) Scout coronal 2/3 of canals with an ISO 010 C-PILOT® file.
5) Scout the root canal with an ISO 010 hand file, determine the working length by using an electronic apex locator in combination with radiographs, irrigate and confirm patency.
6) Set the silicone stopper accordingly.
7) Choose the appropriate preprogrammed torque-controlled electric motor setting according to the VDW.ROTATE® instrument size.
8) With irrigant in the canal, create and confirm a reproducible glide path by using the VDW.ROTATE® ISO 015 004 until the full working length is reached.
9) Routinely irrigate the canal and clean the files cutting flutes of debris upon removal with an alcohol moistened gauze.
10) Expand and secure this glide path by using the next instrument in the sequence, VDW.ROTATE® ISO 020 005 until full working length is reached.
11) Determine the final working length using an electronic apex locator.
12) Depending on the canal curvature and size, choose between the VDW.ROTATE® ISO 025 004 or ISO 025 006 size as following:
 - If the canal is narrow and highly curved, use the VDW.ROTATE® ISO 025 004 until full working length is reached.
 - For all other canals, use the VDW.ROTATE® ISO 025 006 until full working length is reached.
13) Gauge the foramen with an ISO 025 hand file (NiTi or stainless steel):
 - If the instrument is snug at length, preparation is finished and the canal is ready to be obturated.
 - If the instrument is loose at length, continue shaping with the next size ISO 030, and, when necessary proceed with the next file sizes up to ISO 060, gauging after each instrument respectively with a hand file of the same size. Use the same shaping technique up to working length.
14) Use dedicated paper points to dry and dedicated VDW.ROTATE® gutta-percha points or GuttaFusion® to obturate the root canals.

Retreatment procedure with VDW.ROTATE®

1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
2) Open the canal orifice with a gates glidden drill or an ultrasound instrument and remove gutta-percha from the coronal third.
3) If necessary, soften gutta-percha with a solvent (e.g. eucalyptus oil).
4) Choose the appropriate preprogrammed torque-controlled electric motor setting according to the VDW.ROTATE® instrument size.
5) Place the tip of the VDW.ROTATE® Retreatment ISO 025 005 instrument on the gutta-percha and allow it to rotate. VDW.ROTATE® Retreatment ISO 025 005 file has a cutting tip so that the instrument should progress easily in the gutta-percha without pressure.
6) Use VDW.ROTATE® Retreatment instrument ISO 025 005 without downward pressure.
7) Gradually remove the root canal filling with circumferential filing movements.
8) Use a hand instrument (e.g.C-PILOT®) size ISO 008 to negotiate the root canal to full working length.

9) Once full working length to the apex has been negotiated with the previous ISO 008 hand instrument, prepare the root canal with hand files up to size ISO 010 until full working length is reached.
10) Continue the preparation with VDW.ROTATE® instruments for treatment until the required size has been reached.

8. HYGIENE, DISINFECTION, CLEANING AND STERILIZATION

8.1. General basics

- All instruments must be cleaned, disinfected and sterilized prior to each use; this applies the first time instruments supplied in a non-sterile condition are used, and to instruments delivered in a sterile condition that are intended for repeated use. Thorough cleaning and disinfection are essential prerequisites for effective sterilization. Sterility is not guaranteed once the package has been opened, damaged or wet.
- As part of your responsibility for the sterility of instruments, always make sure that only validated methods for cleaning/disinfection and sterilization are used, that devices (washer disinfector, thermal disinfector or sterilizer) are regularly serviced and inspected, and that the validated parameters are maintained during each cycle. For your own safety, always wear protective gloves, glasses and a mask when handling contaminated instruments.
- In addition, always observe all applicable national legal regulations (KRINKO/RKI/BfArM Processing recommendations) and regulations on hygiene relating to your practice or the hospital. This applies in particular to the guidelines regarding prion inactivation.
- Disclaimer: The instructions for processing products prior to use/re-use herein have been validated by VDW GmbH. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for processing. VDW GmbH accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the instructions for use set forth below. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

8.2. Step by Step procedure

Recommended classification according to RKI/BfArM/KRINKO directive (Germany only, intended use): Critical B.

8.2.1. Pre-Treatment at the place of use

Contaminants (particularly pulp and dentine remnants) must be removed immediately after the instrument has been used on a patient (within maximum 2 hours). All further steps in the preparation process must be performed on the same day.

The following procedures must be used to ensure that no contamination can dry on the instruments, and to make subsequent preparation more effective:

Step 1

A prepared interim stand with a new foam disc must be used for each patient. The interim stand must be filled at least two thirds of the way with disinfecting agent.

Step 2

Place in the interim stand prior to pre- disinfection/cleaning and for transport (minimum storage time according to the disinfecting agent manufacturer’s instructions for use: max. two hours).

Please note that the disinfecting agent used during pre-treatment is for personal protection only and is not a substitute for the disinfection stage required after cleaning.

All further steps in the preparation process must be performed on the same day.

Warnings:

- In the case of chemical disinfection, there is a risk of disinfecting agent residues remaining on the instruments.

- Under no circumstances may instruments that have already come into contact with disinfecting agent be used to treat a patient again.
- The blue foam insert for the interim stand must only be used once and must neither be cleaned/ disinfected nor sterilized.

8.2.2. Preparation before cleaning

Step 3

Remove the stopper from the instrument and dispose of it.

Step 4

Then clean to remove contamination under flowing water for at least 3 x 1 minute; to remove contamination manually, use a soft, clean brush or soft, clean cloth that is only used for this purpose; never use metal brushes or wire wool.

Step 5

Check that no visible contamination or remnants remain and repeat the pre-cleaning process if necessary.

8.2.3. Cleaning and disinfection

We recommend an automated procedure to clean and disinfect the instruments. A manual method- including the use of an ultrasonic bath, should only be used if it is not possible to use an automated method, as it is less effective and demonstrates a lower reproducibility.

For processing we do recommend the use of the LavEndo® box (REF. V040479) with Mtwo® modules (REF. V040408).

Warnings:

- In the case of chemical disinfection, there is a risk of disinfecting agent residues remaining on the instruments.

Step 6: Automated cleaning/disinfection (washer-disinfector) – Recommended procedure

The following must be taken into account when selecting a washer-disinfector:

- The effectiveness of the washer-disinfector must have been verified (DGHM approval, FDA clearance or CE mark according to EN ISO 15883).
- Where possible, a tested thermal disinfection program must be used (A0 value ≥ 3000 or at least five minutes at 90°C, or for older equipment at least 10 min. at 93°C).
- The program used must be suitable for the instruments and include the prescribed rinsing cycles.
- Only sterile or low-germ (< 10 CFU/ml) and low-endotoxin (< 0.25 EU/ml) water (ideally highly purified water HPW) must be used for subsequent rinsing.
- The washer-disinfector must be regularly maintained and inspected.

1) If present fit new VDW stoppers to the pre- cleaned instruments or keep the VDW stopper on the new non-sterile instrument.
2) Sort the instruments into the Mtwo® modules (REF. V040408).
3) Place the endo module in the blue lower section the LavEndo® box (REF. V040479) and close it (click into place).
 NOTE: Preparation in the socket module is not permitted.
4) Insert the LavEndo® box horizontally into the washer disinfector.
5) Start the program.
6) After the program has finished, the LavEndo® box from the washer-disinfector.
7) Check and package the instruments as soon as possible after removing them (see section [8.2.4.] and [8.2.5.]), after leaving them to dry further in a clean place if necessary.

An independent, accredited, recognized test laboratory has demonstrated the intrinsic suitability of the instruments for effective automated cleaning and disinfection using the G 7836 CD washer-disinfector (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the Neodisher Medizym cleaning agent (Dr. Weigert, Hamburg). The laboratory used program D-V-MEDIZYM (based on the program DES-VAR-TD (Miele) under worst-case conditions) according to the procedure described above to demonstrate this effectiveness. The Cidex OPA disinfecting agent and the Cidezyme cleaning agent (both Johnson & Johnson GmbH, Norderstedt) were used for pre-treatment. Cleaning and disinfection validation were performed under worst-case conditions (low temperature, low concentration of agent, short soaking time and no drying).

Step 6: Or Manual cleaning and disinfection – If automated procedure is not possible

If the manual method is used for new non-sterile instruments, the stopper needs to be removed and processed separately in small parts basket.

Manual Cleaning:

1) Sort the instruments, without stoppers, into the Mtwo® modules (REF. V040408).
2) Place the endo module in the blue lower section of the LavEndo® box (REF. V040479) and close it (click into place).

NOTE: Preparation in the socket module is not permitted.

3) Place new VDW stoppers in a small parts basket with a sufficiently small mesh size.
4) Insert the LavEndo® box horizontally and, if present, the small parts basket with the new VDW stoppers into the cleaning bath for the prescribed contact time, ensuring that the instruments are sufficiently covered.
5) Then remove the LavEndo® box and, if present, the small parts basket with the VDW stoppers from the cleaning bath and rinse thoroughly with water for at least 3 x 1 min.

Manual disinfection:

6) Insert the LavEndo® box horizontally and the small parts basket with the new stoppers into the disinfection bath for the prescribed contact time, ensuring that the instruments are sufficiently covered.
7) Then remove the LavEndo® box and the small parts basket with the stoppers from the disinfection bath and rinse thoroughly with water at least 5 x 1 min.
8) Dry the LavEndo® box and the small parts basket with the stoppers by blowing them down fully with oil-free, filtered compressed air (or medical compressed air from a can) and then leaving them to dry further in a clean place.
9) Check (see section [8.2.4.]), fit VDW stoppers to the instruments, and package the instruments as soon as possible (see section [8.2.5.]).

An independent, accredited, recognized test laboratory has demonstrated the intrinsic suitability of the instruments for effective manual cleaning and disinfection using the cleaning agent Cidezyme/Enzol and disinfecting agent Cidex OPA (Johnson & Johnson GmbH, Norderstedt (Germany)). The laboratory used the procedure described above to demonstrate this. The Cidex OPA disinfecting agent and the Cidezyme cleaning agent (both Johnson & Johnson GmbH, Norderstedt) were used for pre-treatment.

8.2.4. Inspection and maintenance

Open the LavEndo® boxes and remove the step modules. Check all instruments, modules and LavEndo® boxes after cleaning/disinfection. Defective instruments, boxes and modules should be immediately discarded.

- These defects include:
 - Plastic deformation (e.g. caused by an excessively high temperature during sterilization).
 - Breakage.
 - Loss of color coding or marking.
 - Bent instrument.
 - Untwisted threads.

- Damaged cutting surfaces.
- Dull cutting blades.
- Missing size marking.
- Corrosion.

- Numerical restrictions on reuse are listed under section [4.], “Limitation and restrictions on processing and on re-use”.
- Instruments that are still contaminated must be cleaned and disinfected again.

Warnings:

- Instrument lubricants must not be used at this stage (i.e.for protection purpose before packaging, sterilization and storage).

8.2.5. Packaging

Place the Mtwo® module (REF. V040408) in the lower section of the black sterilization mini box (REF. V040471) and close it with the matching cover. Then package the sterilization trays into disposable sterilization pouches (disposable packaging) that meet the following requirements:

- Compliance with DIN EN 11607/ANSI AAMI ISO 11607.
- Suitable for steam sterilization (withstands temperatures of up to 142°C (288°F) or more, sufficient vapor permeability).

Warnings:

- Sterilization in the sterilization trays without additional packaging is not permitted. The autoclave paper in the boxes is for added safety only.

8.2.6. Sterilization

Only use the sterilization methods listed below; other sterilization methods are not permitted.

Steam sterilization

- Fractionated vacuum/pre-vacuum method (at least three vacuum cycles) or gravity displacement method[1] with sufficient product drying.
- Steam sterilizer in accordance with DIN EN 13060 or DIN EN 285, ANSI AAMI ST79.
- Validated in accordance with DIN EN ISO 17665 (valid IQ and OQ plus product-specific performance qualification (PQ)).
- The maximum sterilization temperature of 138°C (280°F) must not be exceeded; the maximum sterilization temperature includes a tolerance according to DIN EN ISO 17665.
- The drying time that is actually required depends directly on parameters that are the sole responsibility of the user (loading configuration, how many items are loaded and how closely together they are loaded, condition of the sterilizer, etc.) and must therefore be established by the user. However, the drying time must never be less than 20 minutes.

Sterilization procedure	Sterilization temperature	Minimum sterilization time Exposure time at sterilization temperature
Fractionated vacuum/pre-vacuum method	134°C (273°F)	3 minutes[2]
	121°C (250°F)	20 minutes
Gravity method[1 & 3]	134°C (273°F)	15 minutes
	121°C (250°F)	60 minutes

[1] The less effective gravity method should only be used if the fractionated vacuum method is not available. The gravity method is less effective in direct comparison to the fractionated vacuum method. However, it is effective according to the requirements for a processed instrument. All methods are validated and therefore they are efficient and safe for the processing of VDW instruments.

[2] Or 18 min. (prion inactivation).

[3] Gravity method is not applicable for processing within the European Union.

An independent, accredited, recognized test laboratory demonstrated the instruments' intrinsic suitability for effective steam sterilization using the HST 6x6x6 steam sterilizer (Zirbus Technology GmbH, Bad Grund) together with the fractionated vacuum method and the gravity method.

8.3. Storage and transport

- After sterilization, devices must be stored in the sterilization packaging and kept dry and dust-free. In case of damage to the packaging during storage or transport, the processing shall be repeated. Check the instructions for use given by the pouch manufacturer to determine the shelf life of the sterile packaging.

8.4. Disposal

- Instruments shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

8.5. Important Information on material resistance

- When selecting cleaning and disinfecting agents, make sure that they do not contain any of the following substances:
 - Phenol.
 - Strong acids (ph < 6) or strong alkalis (ph > 8); neutral enzymatic cleaning agent recommended.
 - Aldehydes.
 - Anti-corrosive substances (especially di- or triethanolamine).
 - Oxidants (hydrogen peroxide, sodium hypochlorite over 5% strength).
 - NiTi instruments may only be placed in oxidants (< 5% strength sodium hypochlorite) for a maximum of 5 minutes.
 - Solvents.
 - Oils.

Warnings:

- Never clean the instruments, boxes, modules or the interim stand with metal brushes or wire wool.
- Never subject any instruments, boxes, modules or the interim stand to temperatures above 142°C (288°F). It is particularly important to ensure that the products to be sterilized are not stored too close to the walls or floor of the steam sterilizer (risk of excessive temperature and deformation).
- The blue foam insert for the interim stand must only be used once and used blue foam inserts must not be either cleaned/disinfected or sterilized.

The following must be taken into account when selecting cleaning and disinfecting agents:

- They must be suitable for cleaning and disinfecting instruments made from metal and plastic.
- The disinfecting agent must be aldehyde-free (Cidex OPA is permitted due to its special recipe).
- It must be compatible with the instruments.
- A disinfecting agent with verified effectiveness (VAH/DGHM approval, FDA clearance or CE mark) must be used and this must be compatible with the cleaning agent used.
- If a thermal disinfection process is not used, a suitable disinfecting agent with verified effectiveness (VAH/DGHM approval, FDA clearance or CE mark) must also be used and this must be compatible with the cleaning agent used.
- Neutralization must not be necessary (cleaning agent).
- The cleaning agent, if applicable, must be suitable for ultrasonic cleaning (no foaming).

EN

- Combined cleaning agents/disinfecting agents must not be used.

The concentrations, temperatures and contact times specified by the manufacturer of the cleaning agent and disinfecting agent as well as the minimum specifications for subsequent rinsing must be strictly adhered to. Rinse aids must not be used.

Only use freshly prepared solutions and low-germ (<10 CFU/ml) water; tap water that is particularly hard (≥14°dH) is not suitable for this (risk of lime residue).

9. ADDITIONAL INFORMATION

- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
- To get a free printed copy of IFU with appropriate language please see section “Request an IFU” on website www.vdw-dental.com/ifu.

10. LIST OF APPLICABLE SYMBOLS

Symbols	Description
	**Catalogue number:** Indicates the manufacturer’s catalogue number (SKU). This symbol is accompanied by the manufacturer’s catalogue number adjacent to the symbol.
	**Use-by-date:** Indicates the date after which the medical device is not to be used. This symbol is accompanied by a date to indicate that the medical device should not be used after the end of the year, month or day shown. The date is expressed in YYYY-MM-DD (or YYYY-MM) format.
	**Manufacturer:** Indicates the medical device legal manufacturer [for a Medical Device]. Indicates the procedure pack producer [for a Procedure Pack].
	**Consult instructions for use or consult electronic instructions for use:** Indicates the need for the user to consult the instructions for use. This symbol could be accompanied by the electronic IFU website address.
	**Sterilized using radiation:** Indicates a medical device that has been sterilized using irradiation.
	**Batch code:** Indicates the manufacturer’s batch code. This symbol is accompanied by the manufacturer’s batch code adjacent to the symbol.
	**Assortment:** Identifies the package includes an assortment of types or sizes.
	**NiTi:** Identifies a material or substance contained in a product or the material from which the product is made. “XXX” is replaced by material letter code (e.g., “SST” for Stainless Steel) or a word (e.g., “Latex”). NiTi: Nitinol / Nickel titanium
	**Do not use if package is damaged:** Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.
	**Packaging unit:** Indicates the number of pieces in the package. “X” is replaced by the number of pieces in the package.
	**Clockwise Rotation:** Identifies the control or the indicator for clockwise rotary motion or for the adjustment of rotation speed.
	**Sterilizable in a steam sterilizer (autoclave) at temperature specified:** indicates that the instrument is sterilizable in a steam sterilizer (autoclave) at the specified temperature. “XXX” is replaced by a temperature value.
 0123	**CE marking:** CE Marking of Conformity according to EU MDD/MDR requirements. Note: When a notified body is involved, the CE marking is accompanied by the identification number (four-digit number) of the notified body.
	**Medical device:** Indicates the item is a medical device.
	**Caution:** Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.

Symbols	Description
	**Keep away from sunlight:** Indicates a medical device that needs protection from light sources. This symbol can also mean “Keep away from heat”.
	**Single sterile barrier system:** Indicates a single sterile barrier system.
	**Date of manufacture:** Indicates the date when the medical device was manufactured. This symbol is accompanied by a date to indicate the date of manufacture. The date is expressed in YYYY-MM-DD (or YYYY-MM) format.
	**Country of manufacture:** Indicates the country of manufacture of the products. “CC” is replaced by the two or three letter country code as defined in ISO 3166-1.

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