VOCO

Bifix® QM

EN **Instructions for use**
MD **EU Medical device**

CE 0482

Product description:
Bifix QM is a radiopaque, dual-curing composite-based luting system for permanent adhesive luting of metal, ceramic and composite restorative pieces.
Bifix QM is available in three shades: universal, transparent and white-opaque. QuickMix syringes deliver bubble-free, perfectly mixed material that can be applied directly from the syringe. The **Bifix QM** system achieves excellent adhesion with a suitable adhesive and **Ceramic Bond** for silanisation of the restoration.

Indications:
Adhesive luting of inlays, onlays, veneers, crowns, partial crowns, bridges, root posts and adhesive bridges (Maryland bridges).

Contraindications:
Bifix QM contains methacrylates, benzoyl peroxide, BHT and amines. **Bifix QM** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
Bifix QM is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
Bifix QM should only be used by a professionally trained dental practitioner.

Use:

For detailed information please also refer to the instructions for use of the used adhesive and the silane **Ceramic Bond**.

1. **Preparation**
 Remove the temporary if necessary. Clean, thoroughly rinse and slightly dry the tooth with air.
 Note: The working area must be kept clean and dry.
2. **Conditioning the restoration**
 a. Full ceramic restorations
 The inner adhesive surfaces of the restoration should be micro-etched with a hydrofluoric acid etchant.
 Place **Ceramic Bond** (silane) in a mixing palette. Apply **Ceramic Bond** to the inside of the restoration with a brush and allow it to take effect for 60 s. Afterwards, dry the surface with air for 5 s. Do not touch or otherwise contaminate the conditioned surfaces.
 b. Metal, composite and zirconium dioxide restorations
 Sandblast the surfaces with a maximum of 1.5 - 2 bar pressure and a suitable blasting agent. Follow the manufacturer's instructions for use for the respective material.
 Carefully remove the blasting abrasive residue and dry the adhesive surfaces. Apply **Ceramic Bond** to the conditioned surfaces with a brush, allow to take effect for 60 s and then dry thoroughly with air.
 Do not touch the adhesive surfaces afterwards and quickly lute the restoration with **Bifix QM** while observing the following points.
3. **Application of the adhesive**
 Use **Bifix QM** in conjunction with a suitable adhesive. Follow the manufacturer's instructions for use.
4. **Luting the restoration with Bifix QM**
 Insulate/block out any removable prosthodontic elements around the connectors before luting. Only mix **Bifix QM** immediately before luting, as described in the following: Due to technical reasons there may be slight differences in the filling levels of the two syringe chambers prior to first use. Therefore remove syringe cap and extract material until both chambers deliver uniformly. Attach a mixing tip type 11. The notches on the syringe and the mixing tip must line up exactly. Turn the mixing tip 90° clockwise to lock the cannula into place. The material is automatically mixed in the mixing tip and can be applied directly onto the prepared contact areas. If required, the optional intraoral tip type 4 can be placed onto the mixing tip. Processing time in reduced workspace lighting is approximately 2.5 – 3.5 min. Setting time is 3 min from the time of inserting the restorative piece.
 Insert the piece and fix it in place with slight pressure.
 Immediately remove excesses of **Bifix QM** with a foam pellet, disposable brush or similar. Use interproximal floss.
 Advice for luting Maryland bridges
 When luting Maryland bridges, depending on the (universal) adhesive used, the enamel must be conditioned by etching with phosphoric acid gel before applying the adhesive.
5. **Luting of fibre-reinforced resin posts**
 1. For the cementation of root posts, the adhesive is applied into the root canal and onto the occlusal surface of the core according to the manufacture's instructions.
 2. It is recommended to apply **Bifix QM** directly into the root canal using fine intraoral tips type 1. Do not use a lentulo! After preparing the root post according to the manufacturer's instructions, it can additionally be moistened with **Bifix QM**. Insert the post making sure to leave a small amount of excess material.
 3. Light-cure for min. 40 s to fix the root post.

 The attached mixing tip is intended for single use only. After use, store the QuickMix syringe tightly sealed. Check that material flows freely through the openings before reusing.
 The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories.
6. **Light-curing**
 The polymerisation time can be controlled by additional light-curing. Light-cure the material at the luting margins for at least 20 s per segment with a polymerisation lamp with a light output of at least 500 mW/cm². Hold the tip of the light guide as close as possible to the surface to be cured. The restorative piece must be fixed in its desired position with a suitable instrument during light-curing.
 Advice on oxygen inhibition
 Like all resin-based materials, **Bifix QM** surfaces in contact with oxygen do not cure completely (approximately 50 µm inhibited layer).
 Therefore, remove the excess material completely and cover the margins of the restorative piece with oxygen-impermeable protective gel.
7. **Finishing/Fluoridation**
 The processing of any excesses is possible directly after the polymerisation. Remove cured excess material with a finishing diamond or flexible polishing disk. Finish approximate areas with finishing and polishing strips.

Warnings, precautionary measures:
- Only **use** suitable adhesive materials.
- Some enamel/dentine adhesives are not compatible with dual-curing composites. The instructions for use of the respective adhesive must be observed.
- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- Irritation to the pulp may occur if dentine thickness is insufficient (less than 1 mm). In such cases, take appropriate measures to protect the pulp.
- Phenolic substances, especially preparations containing eugenol and thymol, lead to curing failures in composites. The use of zinc oxide eugenol cements or other eugenol-containing materials in combination with **Bifix QM** is therefore to be avoided.
- **Bifix QM** is not suitable for fixing brackets due to its high adhesion values.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):
Barium aluminium borosilicate glass, HEDMA, BisGMA, fluorosilicate glass, fumed silica, initiators, stabilisers, pigments

Storage:
Storage at 4°C-23°C. Do not use after the expiry date. Refrigeration is recommended.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Bifix QM** must be reported to VOCO GmbH and the responsible authority.

Note:
The Summary of Safety and Clinical Performance of **Bifix QM** can be found in the European database on medical devices (EUDAMED – https://ec.europa.eu/tools/eudamed). Detailed information can also be found at www.voco.dental.

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VC 60 AA1210 E1 0226 99 © by VOCO