VOCO

VOCO Ionofil® Molar CE 0482

EN Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 9917-1

Product description:
VOCO Ionofil Molar is a radiopaque glass ionomer cement for fillings with chemical adhesion to enamel and dentine and continuously high release of fluoride ions.
VOCO Ionofil Molar is available in the shades A2 and A3.

Indications:
- Restoration of V-shaped defects and cervical enamel or root erosions including class V cavities
- Restoration of deciduous teeth (especially class I)
- Fillings of class II cavities (not occlusion-bearing) by means of tunnel preparation
- Build-up fillings and linings, core build-up
- Temporary long-term treatment of class I and II cavities

Contraindications:
VOCO Ionofil Molar contains polyacrylic acid and parabens. **VOCO Ionofil Molar** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
VOCO Ionofil Molar is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
VOCO Ionofil Molar should only be used by a professionally trained dental practitioner.

Use:

Preparation:
Prepare the cavity using a minimally invasive approach (removing only carious tooth substance). The layer thickness of the restorations must not be less than 1 mm. Avoid thinly tapering margins. Clean the cavity, rinse with water and dab or blow dry. Do not overdry the dentine: it should have a moist, shiny appearance. If the residual dentine thickness is insufficient, the area close to the pulp should be covered with a suitable calcium hydroxide product.
Apply **VOCO Ionofil Molar** at room temperature (19 °C - 23 °C); in case of high temperatures the use of a cooled glass plate is advantageous. Low temperatures prolong the working time.
Furthermore, **VOCO Ionofil Molar** is suitable for use for the purposes of ART-technique (Atraumatic restorative treatment).

Mixing powder/liquid:
Shake powder well before use!
Mix one level measuring scoop of powder with one drop of liquid. This corresponds to a mixing ratio of 4.0 parts powder to 1 part liquid by weight.
Hold dropping bottle in vertical position and place 1 drop of liquid onto a dry glass or plastic plate. Hold bottle high enough so that a single drop can fall down freely. Only use drops without enclosed air bubbles!
Mix the powder into the liquid in two equal portions using a clean metal spatula, or preferably with a clean **plastic spatula**. Only mix in the second portion when the first portion has been mixed in homogenously. The mixed cement has a firm, packable consistency.
Avoid any contact with humidity!
Never add new liquid during mixing!

Mixing time	Working time	Setting time
50 - 60 s	2 - 3 min	4 - 5 min

The mixing time takes approx. 50 - 60 s and should not exceed 1 min.

Application:
During the placing of a filling the working area must be kept dry. The use of a rubber dam is recommended.
Glass ionomer cements are susceptible to humidity (breath, saliva) or, respectively, dehydration during the setting time, especially directly after application. If possible, cover immediately with a strip or matrix or apply a protective varnish.
Finishing and polishing (flexible polishing discs) of the fillings should follow 7 min after application at the earliest. Afterwards, cover again with varnish.

Warnings, precautionary measures:
- Remove any remaining liquid from the outside of the liquid dispenser (by tissue) to keep opening clean for a correct mixing ratio.
- **VOCO Ionofil Molar** should only be mixed immediately before use to ensure that the entire working time is available for placing and modelling. The development of fissures due to delayed use of the mixed cement can thus be prevented.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):
Fluorosilicate glass, polyacrylic acid, water, tartaric acid, colour pigments, parabens

Storage instructions and application method:
Store in dry conditions at 15-23°C. Tightly close powder and liquid container immediately after use. Do not use after the expiry date.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **VOCO Ionofil Molar** must be reported to VOCO GmbH and the responsible authority.

Note:
The Summary of Safety and Clinical Performance of **VOCO Ionofil Molar** can be found in the European database on medical devices (EUDAMED – https://ec.europa.eu/tools/eudamed). Detailed information can also be found at www.voco.dental.

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VC 60 AA1441 E1 0824 99 © by VOCO