zeta 3 FOAM

SURFACES

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EN - Alcohol-free foam spray disinfectant and cleaner for delicate medical device surfaces.
Zeta 3 Foam is an alcohol-free, ready-to-use broad spectrum biocidal disinfectant and cleaner, prepared and tested in compliance with new harmonised European disinfection regulations. Its excellent compatibility with materials makes Zeta 3 Foam particularly suitable for rapid high level disinfection between patients of delicate medical device surfaces (handpieces, contra-angles, dental units, etc.).
Composition:
100 g of Zeta 3 Foam contains 0.315 g dimethyl-didecyl-ammonium chloride, 0.075 g alkylbenzyl-dimethyl ammonium chloride, non-ionic surfactants, additives, auxiliaries and water to 100 g.
Range of action:
Bactericidal, fungicidal (test on Candida albicans), tuberculocidal, virucidal. Suspension and surface activity, tested in accordance with EN13727, EN13624, EN14348, EN14561, EN14562, EN14563. Virucidal activity tested in accordance with EN14476.
Use:
Spray Zeta 3 Foam over the area until the surfaces and areas to be disinfected are sufficiently covered. Leave it to work, then remove the disinfectant emulsion with a wipe and finally leave to dry. Alternatively, it is possible to spray Zeta 3 Foam on a tissue and wipe it over surfaces and areas to be disinfected.
Virucidal activity: 1 minute "dirty" conditions (EN 14476)
Suspension activity: 5 minutes bactericidal (EN13727)
fungicidal (EN1624, test on Candida albicans)
tuberculocidal (EN14348)
Surface activity: 5 minutes bactericidal (EN14561)
fungicidal (EN14562, test on Candida albicans)
tuberculocidal (EN14563)
Warnings, precautions and conservation conditions:
This product has been designed for use by qualifi ed operators with knowledge of the required safety standards. Use protective gloves. Operate in well ventilated areas. Store the closed container in a cool dry place. Prevent Zeta 3 Foam from coming in contact with other disinfectants or cleaners, particularly with anionic or aldehyde surfactants (formaldehyde, glutaraldehyde, glyoxal, etc.). Any slight difference in the perfume is due to the presence of natural fragrances and does not jeopardise quality of the product. Once the container has been opened and the applicator inserted, keep in a vertical position.

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